



Premium Surgiclip™ titanium clips

PRODUCT INFORMATION

Code	Description	Qty Per Case
134031	Premium Surgiclip™ M-11.5, 20 Medium Titanium Clips	6
134044	Premium Surgiclip™ M-9.75, 15 Medium Titanium Clips	6
134048	Premium Surgiclip™ L-13.0, 15 Large Titanium Clips	6

MATERIALS OF MANUFACTURE

Code	Description
134031	ABS GE, 301, 302, 304 Stainless Steel, Aluminum 5052, Tyvek Dupont Adhesive, Titanium
134044	ABS GE, 301, 302, 304 Stainless Steel, Aluminum 5052, Polycarbonate GE, Acetel Dupont Delrin, Titanium composition, Tyvek Dupont Adhesive
134048	ABS GE, 301, 302, 304 Stainless Steel, Aluminum 5052, Polycarbonate GE, Acetel Dupont Delrin, Titanium composition, Tyvek Dupont Adhesive

SHIPPING INFORMATION

	Inner Box Size (mm)
134031	L 356 × W 141 × H 106
134044	L 356 × W 141 × H 106
134048	L 387 × W 154 × H 106

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Ethylene Oxide
Shelf Life	5 years
CE Marked	Class IIb According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-premium-surgiclip-tech-sheet-1516917

Medtronic
Further, Together