



Endobag™

Specimen Retrieval System

PRODUCT INFORMATION

Code	Description	Qty Per Case
250-30	Endobag™ 3" × 6" Specimen Retrieval System – 20.9 cm Shaft Length	20
250-40	Endobag™ 5" × 8" Specimen Retrieval System – 20.9 cm Shaft Length	20

MATERIALS OF MANUFACTURE

ABS, Polyurethane

SHIPPING INFORMATION

	Inner Box Size (mm)
250-30	L 278 × W 146 × H 295
250-40	L 278 × W 146 × H 295

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class IIa According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-endobag-tech-sheet-2338468

Medtronic
Further, Together