

Endo Shears™ 5 mm Single Use Instrument with Monopolar Cautery



PRODUCT INFORMATION

Code	Description	Qty Per Case
174601	Endo Shears™ Long 5 mm, Single Use Instrument with Monopolar Cautery	1
176643	Endo Shears™ 5 mm, Single Use Instrument with Monopolar Cautery	6

MATERIALS OF MANUFACTURE

301, 302, 303, 304, and 416 Stainless Steel, 10% Glass Filled Polycarbonate, Polycarbonate 20% glass, Pvc Grade 6k10, Style Vw-1, Aluminum 6061-T6, Union Carbide

SHIPPING INFORMATION

	Inner Box Size (mm)	Shipper Box Size (mm)
174601	L 162 × W 27 × H 692	L 702 × W 181 × H 171
176643	L 610 × W 191 × H 70	N/A

ADDITIONAL INFORMATION

Latex Free	✓
DEHP Free	✓
Sterilization method	Radiation
Shelf Life	5 years
CE Marked	Class IIb According to Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. EMEA-AI-2000024 - 20-weu-endo-shears-5-mm-tech-sheet-4750109

Medtronic
Further, Together