

OPERATING
ROOM
ESSENTIALS

THE COMPLETE
SOLUTION
YOU CAN TRUST



VersaOne™
Blunt Trocar
VersaOne™ Access System

Medtronic
Further. Together

**EFFICACY
& SIMPLICITY.
EVERY TIME.**

The VersaOne™ Blunt Trocar with threaded anchor is designed for surgeons that use a Hasson or open technique. Available in 12 mm size.

Safe Insertion*


Blunt tip aids in reducing potential injury to internal structures upon insertion

Outstanding Visualization

Universal, clear cannula provides visibility to all tissue layers

Secure Placement

Threaded anchor ensures trocar placement and maintains pneumoperitoneum



Blunt Obturator

- Rounded distal tip designed to allow for a less traumatic entry
- Low profile hub accommodates a variety of patients

Universal Cannula

- Increased inner diameter accommodates standard 12 mm laparoscopic instruments
- Funneled entry, enhanced seal allows smooth instrument exchange¹

Threaded Anchor

- Secure suture tie down with two suturing wings
- Adjust anchor to desired location on cannula using a spring lever

VersaOne™ Blunt Trocar with Threaded Anchor

	12 mm
Standard Length Smooth Cannula (100 mm)	BPT12STS

Sold three units per box

Contact your Medtronic representative today for more information about the VersaOne™ Access System or visit medtronic.com/covidien.

* Compared to bladed trocars

1. Based on internal test report#2143-123. Applies to Medtronic 12 mm trocars, when compared to the Versaport™ Bladeless 12 mm trocar. March 2013.

© 2016 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

16-eu-versaone-blunt-sell-sheet-906623

Medtronic



To contact us, please visit
medtronic.com/covidien/support/emea-customer-service

Use scan app to read