



VersaOne™ optical trocar

PRODUCT INFORMATION

| Code | Description | Ship Case Qty |
|---------------------------|---|---------------|
| ONB5SHF | VersaOne™ optical trocar 5 mm Short (70 mm) Fixation Cannula | 6 per box |
| ONB12SHF | VersaOne™ optical trocar 12 mm Short (70 mm) Fixation Cannula | " |
| ONB5STF | VersaOne™ optical trocar 5 mm Standard (100 mm) Fixation Cannula | " |
| ONB5STF2C | VersaOne™ optical trocar 5 mm Standard (100 mm) with Dual Pack Fixation Cannula | " |
| ONB11STF | VersaOne™ optical trocar 11 mm Standard (100 mm) Fixation Cannula | " |
| ONB12STS | VersaOne™ optical trocar 12 mm Standard (100 mm) Smooth Cannula | " |
| ONB12STF | VersaOne™ optical trocar 12 mm Standard (100 mm) Fixation Cannula | " |
| ONB15STF | VersaOne™ optical trocar 15 mm Standard (100 mm) Fixation Cannula | " |
| ONB5LGF | VersaOne™ optical trocar 5 mm Long (150 mm) Fixation Cannula | " |
| ONB11LGF | VersaOne™ optical trocar 11 mm Long (150 mm) Fixation Cannula | " |
| ONB12LGF | VersaOne™ optical trocar 12 mm Long (150 mm) Fixation Cannula | " |

MATERIALS OF MANUFACTURE

| | |
|-----------|---|
| Cannula | <ul style="list-style-type: none"> ▪ Transparent polyamide ▪ Polysoprene (synthetic rubber) ▪ C-Flex Thermoplastic Elastomer ▪ Marlex 9006 HDPE ▪ Silicone oil |
| Obturator | <ul style="list-style-type: none"> ▪ Transparent polyamide ▪ Acrylonitrile-butadiene-styrene ▪ Versaflex Thermoplastic Elastomer |

DIMENSIONS

| | | |
|---------------|----------|--------|
| Trocar length | Short | 70 mm |
| | Standard | 100 mm |
| | Long | 150 mm |

PACKAGING (APPLIES TO ALL CODES)

| | Material | Qty |
|-----------------|----------------------------|-----|
| Out packaging | Carton box | 1 |
| Inner packaging | Soft breather Tyvek® pouch | 6 |

SHIPPING INFORMATION

| | Shipped case dimensions (cm) |
|--|------------------------------|
| ONB5SHF, ONB5STF | L 30,5 × W 12,7 × D 12,7 |
| ONB5LGF, ONB5STF2C | L 36,8 × W 12,1 × D 14,0 |
| ONB11STF, ONB12SHF, ONB12STS, ONB12STF, ONB15STF | L 22,9 × W 12,3 × D 26,7 |
| ONB11LGF, ONB12LGF | L 28,2 × W 12,3 × D 26,7 |

ADDITIONAL INFORMATION

| | |
|----------------------|---|
| Latex Free | ✓ |
| DEHP free | ✓ |
| Sterilization method | Ethylene Oxide |
| Shelf Life | 5 years |
| CE Marked | Class IIa According to EU directive 93/42/EEC, Annex II.3 |

Medtronic

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
19-weu-versaone-optical-ts-sept19-3588034

[medtronic.com/covidien/en-gb/index.html](https://www.medtronic.com/covidien/en-gb/index.html)