



# VersaOne™

## fascial closure system bladed trocar

### PRODUCT INFORMATION

Code	Description	Ship Case Qty
<b>B12STFCS</b>	VersaOne™ Fascial Closure System Bladed trocar 12 mm Standard (100 mm) Fixation Cannula	6 per box

### MATERIALS OF MANUFACTURE

Cannula	<ul style="list-style-type: none"><li>▪ Nylon</li><li>▪ Acrylonitrile-butadiene-styrene</li><li>▪ Polyisoprene Rubber</li></ul>
Guide	<ul style="list-style-type: none"><li>▪ Acrylonitrile-butadiene-styrene</li></ul>
Obturator	<ul style="list-style-type: none"><li>▪ Acrylonitrile-butadiene-styrene</li><li>▪ Nylon</li><li>▪ Stainless Steel</li></ul>

### DIMENSIONS

Standard	100 mm
----------	--------

### PACKAGING

	Material	Qty
Outer packaging	Carton box	1
Inner packaging	Soft breather Tyvek® pouch	6

Package contains cannula, obturator and suture guide

### SHIPPING INFORMATION

	Shipped case dimensions (cm)
<b>B12STFCS</b>	L 22,9 × W 12,3 × D 26,7

### ADDITIONAL INFORMATION

Latex Free	✓
DEHP free	✓
Sterilization method	Ethylene Oxide
Shelf Life	3 years
CE Marked	Class IIa According to EU directive 93/42/EEC, Annex II

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.  
™\*\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.  
19-weu-versaone-fcs-bladed-trocar-technical-sheet-3742202

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

**Medtronic**  
Further, Together