



VersaOne™

fascial closure system

optical trocar

PRODUCT INFORMATION

Code	Description	Ship Case Qty
ONB12STFCS	VersaOne™ Fascial Closure System Optical trocar 12 mm Standard (100 mm) Fixation Cannula	6 per box

MATERIALS OF MANUFACTURE

Cannula	<ul style="list-style-type: none"> Nylon Acrylonitrile-butadiene-styrene Polyisoprene Rubber
Guide	<ul style="list-style-type: none"> Acrylonitrile-butadiene-styrene
Obturator	<ul style="list-style-type: none"> Acrylonitrile-butadiene-styrene Nylon Stainless Steel

DIMENSIONS

Standard	100 mm
----------	--------

PACKAGING

	Material	Qty
Outer packaging	Carton box	1
Inner packaging	Soft breather Tyvek® pouch	6

Package contains cannula, obturator and suture guide

SHIPPING INFORMATION

	Shipped case dimensions (cm)
ONB12STFCS	L 22,9 × W 12,3 × D 26,7

ADDITIONAL INFORMATION

Latex Free	✓
DEHP free	✓
Sterilization method	Ethylene Oxide
Shelf Life	3 years
CE Marked	Class IIa According to EU directive 93/42/EEC, Annex II

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 19-weu-versaone-fcs-technical-sheet-2964525

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

Medtronic
Further, Together