



Ti-Cron™ suture

Reference Guide

Ti-Cron™ suture has been your trusted brand for more than 20 years — with a proven history of performance

Features and Benefits¹

- Soft and supple
- Silicone (impregnated) coating designed for handling and pliability
- Knot security

Pledgets

Available as firm and soft pledgets in four sizes:

- 3 mm by 3 mm
- 3 mm by 5 mm
- 3 mm by 7 mm
- 4.8 mm by 9.5 mm

Frequent Uses

- Soft tissue approximation
- Ophthalmic
- Cardiovascular
- Neurological surgery

Ti-Cron™ Suture Top Codes

Code	Color	Size USP	Length		Needles
			cm	Inches	
88863323-56	5 Blue, 5 White	2-0*	75	30	CV-316
88863219-56	4 Blue, 4 White	2-0*	75	30	Y-31
88863027-79	Blue	5-0	75	30	HOS-14
88863385-56	5 Blue, 5 White	2-0*	75	30	CV-316
88863218-56	4 Blue, 4 White	2-0*	75	30	Y-31

Product Information

Product	Ti-Cron™ Suture ²	Ethibond™ Excel ³
Structure	Braided	Braided
Suture Type	Non-Absorbable	Non-Absorbable
Composition	Polyester	Polyester
Coating	Silicone	Polybutylate
Color	Blue, White	Green, Undyed (white)
Tensile Strength	Permanent	Permanent
Absorption Profile	Permanent	Permanent
Sizes	7-0 through 5	6-0 through 2, 5
Indications	Ti-Cron™ polyester suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological surgery.	Ethibond™ Excel suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.
Contraindications	None known.	
Box Quantities	1/2, 1, 2 or 3 dozen	

*Double-arm

1. Ti-Cron™ Brochure - P120070.

2. Ti-Cron™ Instructions for Use.

3. Ethibond™ Excel Instructions for Use.

© 2017 Medtronic. All rights reserved. 17-emea-wc-guide-ti-cron-en-1605650

Medtronic
Further, Together