



POWER, CONNECTIVITY AND MOUNTING ACCESSORIES

for the Capnostream™ 20p
Bedside Patient Monitor



Power and Connectivity Accessories	Part Number
 <p>U.S. Power Cord (3 meters)</p>	014256
 <p>Battery Pack</p>	BP01000
 <p>Nurse Call System Cable ▪ Supplied unterminated</p>	011149

Power and Connectivity Accessories	Part Number
 <p>Power Cable Strain Relief</p>	CS08807
 <p>Fuse, 3.14 amp/250 volt (2/package)</p>	010543
 <p>Paper (6 rolls/package)</p>	010516

Mounting Accessories

Part Number

	Wireless Bridge Adapter (Bernoulli-compatible)	011892	Supports all Nellcor™ OxiNet™ technology bridges: <ul style="list-style-type: none"> – Moxa TCC-80 (DB 9) RS-232 to 422 Adapter – Moxa NPort RS 232 to wired ethernet adapter – Nellcor™ brand version of wireless bridge Reuses Capnostream™ 20p patient monitor feet for table-top use. Handle designed for monitor transport and use with GCX 5" Slide-in Plunger Mounting Plate. Hook and loop straps offer storage of the SpO ₂ cabling and power cords.
	Quick Release Pole Clamp Kit	011782	Clamp fits pole mount (0.75"-1.5"/19mm-38mm) or rail mount (0.39"/10mm x 0.98"/25mm). Note: Never use a temporary mount when a falling Capnostream™ 20p monitor could cause patient injury.
	VESA Standard Mounting Adapter Plate	010713	For mounting Capnostream™ 20p monitor on any VESA standard 75 x 75 mounting system as required by various industry mounting solutions.
	Roll Stand, Handle and Basket	10002385	Features tilt/swivel adjustments. Post height is 38"/96.5 cm, overall height with casters, floor to mounting plate is 43"/109.2 cm. Mounting plate sold separately (010713).
	CPU Channel Mount for Nuvon VEGA™ System	GXC Part #RS-0009-40 ¹	When the Capnostream™ 20p patient monitor is connected to the Nuvon VEGA™ System, the IDM-MG 3000 network appliance can be mounted on the post of a GCX Corporation roll stand.
	Wall Channel, Vertical	10003482	Standard channel. Channel permits monitor to be set anywhere in the channel.
	Wall Mount Pivot Arm	10003484	12" pivot arm. Pivots at the slide for lateral adjustability. Wall Channel sold separately (10003482).

¹Please visit gcx.com for product details and availability. Capnostream™ 20p bedside monitor accessories are designed to accommodate several mounting options. Additional mounting options are available through GCX Corporation.

© 2018, 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 03/2019-18-PM-0048(2)-[WF#2559308]

6135 Gunbarrel Avenue
Boulder, CO 80301

800.635.5267

medtronic.com/covidien

Medtronic