



CLIP APPLIER SPECIFICATIONS

Comparing Medtronic and Ethicon clip appliers

Product Name	Product Code	Clip Size	Clip Dimension (Open) (W)	Clip Dimension (Closed) (L)	Number of Clips	Clip Material	Disposable
Medtronic							
Premium Surgiclip™ III Clip Applier S-9.0"	133650	Small	2.25 mm	3.7 mm	20	Titanium	Disposable
Premium Surgiclip™ II Clip Applier M-9.75"	134051	Medium	4.6 mm	6.0 mm	20	Titanium	Disposable
Premium Surgiclip™ II Clip Applier M-11.5"	134053	Medium	4.6 mm	6.0 mm	30	Titanium	Disposable
Premium Surgiclip™ Clip Applier L-13.0"	134048	Large	6.5 mm	11.0 mm	15	Titanium	Disposable
Ethicon†							
Ligaclip™* 9.38"	MCS20	Small	2.1 mm	3.63 mm	20	Titanium	Disposable
Ligaclip™* 9.38"	MSM20	Medium	4.3 mm	6.0 mm	20	Titanium	Disposable
Ligaclip™* 11"	MCM20	Medium	4.3 mm	6.0 mm	20	Titanium	Disposable
Ligaclip™* 11"	MCM30	Medium	4.3 mm	6.0 mm	30	Titanium	Disposable
Ligaclip™* 13"	MCL20	Large	6.3 mm	10.8 mm	20	Titanium	Disposable

† Based on Ethicon US, LLC Website. Accessed May 14, 2020.

For more information, call **800-772-8772**
or visit us at [medtronic.com/covidien](https://www.medtronic.com/covidien)

©2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 12/2020 – US-AI-2000019 – [WF# 2231759]