

SIZING CHART

POCKET GUIDE

Shiley™ Flexible

Medtronic
Further, Together

Shiley™ Flexible Adult Tracheostomy Tubes with Disposable Inner Cannula

Size	SKU#		Specifications					
			Outer Cannula Inner Diameter (mm I.D.)	Inner Cannula Inner Diameter (mm I.D.)	Outer Cannula Outer Diameter (mm O.D.)	Length (mm)	Cuff Resting Diameter (mm)	
Jackson	Cuffed	Cuffless						
4	4CN65H	4UN65H	6.5	5.5	9.4	62	20.6	
5	5CN70H	5UN70H	7.0	6.0	10.1	68	23.0	
6	6CN75H	6UN75H	7.5	6.5	10.8	74	25.4	
7	7CN80H	7UN80H	8.0	7.0	11.4	77	25.4	
8	8CN85H	8UN85H	8.5	7.5	12.2	79	26.6	
9	9CN90H	9UN90H	9.0	8.0	12.7	79	27.6	
10	10CN10H	10UN10H	10.0	9.0	13.8	79	28.6	

Previous Shiley™ Tracheostomy Tubes with Disposable Inner Cannula

SKU#		Specifications			
Cuffed	Cuffless	Outer Cannula Inner Diameter (mm I.D.)	Outer Cannula Outer Diameter (mm I.D.)	Length (mm)	Cuff Resting Diameter (mm)
4DCT	4DCFS	5.0	9.4	62	20.0
6DCT	6DCFS	6.4	10.8	74	24.0
8DCT	8DCFS	7.6	12.2	79	27.0
10DCT	10DCFS	8.9	13.8	79	29.0

The Shiley sizing guide is not intended to replace the product directions for use. It is intended as a guide only and should not replace institutional policies or physicians' orders. Please always refer to other product instructions for use for complete use instructions.

© 2015 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 15-AW-0045 (US)