

DUAL CHAMBER TEMPORARY EXTERNAL PACEMAKER

Model 5392

Specifications

Modes	DDD, DOO, DDI, AAI, AOO, VVI, VOO
Basic Pacing Rates	30-200 ppm
Upper Rate	80-230 ppm
Rapid Atrial Pacing Rates	80-800 ppm
Output Amplitude	Atrial: 0.1 – 20 mA Ventricular: 0.1 – 25 mA
Pulse Width	Atrial: 1.0 ms Ventricular: 1.5 ms
Sensitivity	Atrial: 0.4 – 10 mV Ventricular: 0.8 – 20 mV
A-V interval	Paced A-V (PAV): 50-250 Auto 200 – 300 Manual Sensed A-V (SAV): 50 – 250
Refractory Period	Atrial: 150 – 500 ms (PVARP) Atrial Refractory after an Atrial event is equal to the AV interval Ventricular: NA
Ventricular Blanking	Pace: 200 ms Sense: 120 ms
Height	20.3 cm (8.0 in)
Width	8.6 cm (3.375 in)
Depth	4.45 cm (1.75 in)
Weight	680 g (24 ounces)
Battery Type	Two IEC type LR6-sized (AA-sized) 1.5 V alkaline batteries (Duracell MN1500, Eveready E91 or equivalent)
Battery Life	9 days typical, 7 days minimum



The Next Generation Temporary Pacemaker

- Enhanced user interface with improved visibility and increased responsiveness
- Improved low battery indicator
- Single and dual chamber pacing capability
- 7-year service life

Medtronic

Brief Statement

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

medtronic.ca

99 Hereford Street Tel: (905) 460-3800
Brampton, Ontario Fax: (905) 460-3998
L6Y 0R3
Canada

Medtronic and the Medtronic logo are trademarks of Medtronic.

*Third party brands are trademarks of their respective owners.

All other brands are trademarks of a Medtronic company.

UC201205468b EC © 2017 Medtronic. All Rights Reserved.

Medtronic