

DEVICE AND SOFTWARE COMPATIBILITY: MINIMED™ 770G SYSTEM APPS

Two smart device apps are available with the MiniMed™ 770G system.

The following devices are currently compatible with our apps. Please note, compatibility may vary between the two apps. We will continue to update this page as we gain new compatibility for an operating system (OS).



MiniMed™ Mobile App allows those wearing the system to view their data on a smart device.

MiniMed Mobile app 1.1.1 APPLE DEVICES AND OPERATING SYSTEMS (IOS)	
Operating system	13.0, 13.3.1, 13.4.1, 13.5.1, 13.6, 13.6.1, 13.7, 14.0, 14.0.1, 14.1, 14.2, 14.2.1, 14.3, 14.4
Compatible Phones	iPhone SE (2020) iPod Touch (6th generation) iPod Touch (7th generation) iPhone 6 iPhone 6 Plus iPhone 6S iPhone 6S Plus iPhone 7 iPhone 7 Plus iPhone 8 iPhone 8 Plus iPhone X iPhone XS iPhone XS Max iPhone XR iPhone 11 – (only OS 13 +) iPhone 11 Pro Max – (only OS 13 +) iPhone 12 – (only OS 14 +) iPhone 12 Mini – (only OS 14 +) iPhone 12 Pro – (only OS 14 +) iPhone 12 Pro Max – (only OS 14 +)

MiniMed™ Mobile app 1.1.1
ANDROID OS & DEVICES

Manufacturer	Model	Android OS
Samsung	<ul style="list-style-type: none"> ■ S8 (SC-02J, SCV36, SM-G950N, SM-G950F, SM-G950W, SM-G950U, SM-G9508, SM-G950U1, SM-G9500) ■ S8+ (SC-03J, SCV35, SM-G955U1, SM-G9550, SM-G955N, SM-G955F, SM-G955W, SM-G955U) 	9
Samsung	<ul style="list-style-type: none"> ■ S9 (SC-02K, SCV38, SM-G960F, SM-G960N, SM-G9600, SM-G9608, SM-G960W, SM-G960U, SM-G960U1) ■ S9+ (SC-03K, SCV39, SM-G965F, SM-G965N, SM-G9650, SM-G965W, SM-G965U, SM-G965U1) 	9, 10
Samsung	<ul style="list-style-type: none"> ■ Galaxy A8 (2018) (SM-A530F, SM-A530K, SM-A530L, SM-A530S, SM-A530N, SM-A530W, SM-A530X) ■ Galaxy A10e (SM-A102*, SM-S102DL) ■ Galaxy A20 (SM-A205*, SM-S205DL, SC-02M, SCV46-J, SCV46-U, SCV46) 	9
	<ul style="list-style-type: none"> ■ Galaxy A7 (2018) (SM-A750F, SM-A750FN, SM-A750G, SM-A750GN, SM-A750C, SM-A750X, SM-A750N) ■ Galaxy A20e (SM-A202F, SM-A202K) ■ Galaxy A21s (SM-A217F, SM-A217F/DS, SM-A217F/DSN, SM-A217M, SM-A217M/DS, SM-A217N) 	10
	<ul style="list-style-type: none"> ■ S10 (SC-03L, SCV41, SM-G973F, SM-G973N, SM-G9730, SM-G9738, SM-G973C, SM-G973U, SM-G973U1, SM-G973W) ■ S10 lite (SM-G770*) – only supports OS 10 and OS 11 ■ S10 5G (SM-G977B, SM-G977N, SM-G977P, SM-G977T, SM-G977U) ■ S10e (SM-G970F, SM-G970N, SM-G9700, SM-G9708, SM-G970U, SM-G970U1, SM-G970W) ■ S10+ (SC-04L, SCV42, SM-G975F, SM-G975N, SM-G9750, SM-G9758, SM-G975U, SM-G975U1, SM-G975W) 	9, 10, 11
Samsung	<ul style="list-style-type: none"> ■ S20 (SM-G980F) ■ S20+ (SM-G985F) ■ S20 5G (SC-51A, SCG01, SM-G9810, SM-G981N, SM-G981U, SM-G981U1, SM-G981V, SM-G981W, SM-G981B) ■ S20+ 5G (SC-52A, SCG02, SM-G9860, SM-G986N, SM-G986U, SM-G986U1, SM-G986W, SM-G986B) ■ S20 Ultra 5G (SCG03, SM-G9880, SM-G988N, SM-G988U, SM-G988U1, SM-G988W, SM-G988B) ■ S20 FE (SM-G780*) ■ S20 FE 5G (SM-G7810, SM-G781B, SM-G781N, SM-G781U, SM-G781U1, SM-G781V, SM-G781W) 	10, 11
Samsung	<ul style="list-style-type: none"> ■ Samsung Galaxy S21 (SM-G990*) ■ Samsung Galaxy S21 5G (SM-G991*) ■ Samsung Galaxy S21+ (5G SM-G996*) ■ Samsung Galaxy S21 Ultra (5G SM-G998*) 	11
Samsung	<ul style="list-style-type: none"> ■ Galaxy Note 8 (SC-01K, SCV37, SM-N9500, SM-N9508, SM-N950F, SM-N950N, SM-N950U, SM-N950U1, SM-N950W, SM-N950XN) 	9
	<ul style="list-style-type: none"> ■ Galaxy Note 9 (SC-01L, SCV40, SM-N9600, SM-N960F, SM-N960N, SM-N960U, SM-N960U1, SM-N960W) 	10
	<ul style="list-style-type: none"> ■ Galaxy Note 10 (SCV45, SM-N9700, SM-N970F, SM-N970U, SM-N970U1, SM-N970W) 	10, 11
	<ul style="list-style-type: none"> ■ Galaxy Note 10 5G (SM-N971N, SM-N971U) 	10, 11
	<ul style="list-style-type: none"> ■ Galaxy Note 10 + (SC-01M, SCV45, SM-N9750, SM-N975C, SM-N975F, SM-N975U, SM-N975U1, SM-N975W) 	10, 11
	<ul style="list-style-type: none"> ■ Galaxy Note 10 + 5G (SM-N9760, SM-N976B, SM-N976Q, SM-N976U, SM-N976V) 	
	<ul style="list-style-type: none"> ■ Galaxy Note 10 Lite (SM-N770F, SM-N770F/DS, SM-N770F/DSM) 	10, 11
	<ul style="list-style-type: none"> ■ Galaxy Note 20 (SM-N980F) 	10, 11
	<ul style="list-style-type: none"> ■ Galaxy Note 20 5G (SM-N9810, SM-N981B, SM-N981N, SM-N981U, SM-N981U1, SM-N981W) 	10, 11
	<ul style="list-style-type: none"> ■ Galaxy Note 20 Ultra (SM-N985F) 	10, 11
<ul style="list-style-type: none"> ■ Galaxy Note 20 Ultra 5G (SC-53A, SCG06, SM-N9860, SM-N986B, SM-N986N, SM-N986U, SM-N986U1, SM-N986W) 	10, 11	
Motorola	<ul style="list-style-type: none"> ■ Moto G6 (XT1925, XT1925-10, XT1925-13, XT1925-6, XT1925-5, XT1925-3, XT1925-7, moto g(6)) 	9
	<ul style="list-style-type: none"> ■ Moto G7 (XT1962, XT1962-4, moto g(7)) – only supports OS 9 ■ Moto G7 Power (XT1955, moto g(7) power) ■ Moto G7 Play (XT1952-2, XT1952-3, XT1952-4, moto g(7) play) 	9, 10
	<ul style="list-style-type: none"> ■ OnePlus 6T (ONEPLUS A6010, ONEPLUS A6013) 	9, 10
	<ul style="list-style-type: none"> ■ OnePlus 6T (ONEPLUS A6010, ONEPLUS A6013) 	9, 10
LG	<ul style="list-style-type: none"> ■ G7 ThinQ, V40 Thin Q, G8 ThinQ, V50 ThinQ, V50S ThinQ, V60 ThinQ 5G 	10
Google	<ul style="list-style-type: none"> ■ Pixel 3, Pixel 3A, Pixel 3 XL, Pixel 3A X 	9, 10, 11
	<ul style="list-style-type: none"> ■ Pixel 4, Pixel 4 XL, Pixel 4A, Pixel 4A 5G 	10, 11
	<ul style="list-style-type: none"> ■ Pixel 5 	11



CareLink™ Connect App allows care partners to follow someone wearing the system.

Carelink Connect app 2.1.1 APPLE DEVICES AND OPERATING SYSTEM (IOS)

Operating system	13.0, 13.3.1, 13.4.1, 13.5.1, 13.6, 13.6.1, 13.7, 14.0, 14.0.1, 14.1, 14.2, 14.2.1, 14.4
Compatible Phones	iPhone SE (2020) iPhone 6 iPhone 6S Plus iPhone 7 iPhone 7 Plus iPhone 8 iPhone 8 Plus iPhone X iPhone XS iPhone XS Max iPhone XR iPhone 11 – (only OS 13 +) iPhone 11 Pro – (only OS 13 +) iPhone 11 Pro Max – (only OS 13 +) iPhone 12 – (only OS 14 +) iPhone 12 Mini – (only OS 14 +) iPhone 12 Pro – (only OS 14 +) iPhone 12 Pro Max – (only OS 14 +)

CareLink™ Connect App Compatibility ANDROID OS & DEVICES

Manufacturer	Model	Android OS
Samsung	<ul style="list-style-type: none"> S8 (SC-02J, SCV36, SM-G950F, SM-G950N, SM-G950W, SM-G9500, SM-G9508, SM-G950U, SM-G950U1) S8+ (SC-03J, SCV35, SM-G955F, SM-G955N, SM-G955W, SM-G9550, SM-G955U, SM-G955U1) 	9
Samsung	<ul style="list-style-type: none"> S9 (SC-02K, SCV38, SM-G960F, SM-G960N, SM-G9600, SM-G9608, SM-G960W, SM-G960U, SM-G960U1) S9+ (SC-03K, SCV39, SM-G965F, SM-G965N, SM-G9650, SM-G965W, SM-G965U, SM-G965U1) 	9, 10
Samsung	<ul style="list-style-type: none"> S10 (SC-03L, SCV41, SM-G973F, SM-G973N, SM-G9730, SM-G9738, SM-G973C, SM-G973U, SM-G973U1, SM-G973W) S10 5G (SM-G977B, SM-G977N, SM-G977P, SM-G977T, SM-G977U) S10+ (SSC-04L, SCV42, SM-G975F, SM-G975N, SM-G9750, SM-G9758, SM-G975U, SM-G975U1, SM-G975W) S10e (SM-G970F, SM-G970N, SM-G9700, SM-G9708, SM-G970U, SM-G970U1, SM-G970W) S10 Lite (SM-G770F, SM-G770U1) – only supports OS 10 and OS 11 	9, 10, 11
Samsung	<ul style="list-style-type: none"> S20 (SM-G980F) S20+ (SM-G985F) S20 5G (SC-51A, SCG01, SM-G9810, SM-G981N, SM-G981U, SM-G981U1, SM-G981V, SM-G981W, SM-G981B) S20+ 5G (SC-52A, SCG02, SM-G9860, SM-G986N, SM-G986U, SM-G986U1, SM-G986W, SM-G986B) S20 Ultra 5G (SCG03, SM-G9880, SM-G988N, SM-G988U, SM-G988U1, SM-G988W, SM-G988B) S20 FE (SM-G780F) 	10, 11
Samsung	<ul style="list-style-type: none"> S20 FE 5G (SM-G7810, SM-G781B, SM-G781N, SM-G781U, SM-G781U1, SM-G781V, SM-G781W) S21 5G (SM-G9910, SM-G991U1, SM-G991W, SM-G991B, SM-G991N) S21 Ultra 5G (SM-G9980, SM-G998U, SM-G998U1, SM-G998W, SM-G998B, SM-G998N) S21+ 5G (SM-G9960, SM-G996U1, SM-G996W, SM-G996B, SM-G996N) 	11

Samsung	▪ Galaxy Note 8 (SC-01K, SCV37, SM-N9500, SM-N9508, SM-N950F, SM-N950N, SM-N950U, SM-N950U1, SM-N950W, SM-N950XN)	9
	▪ Galaxy Note 9 (SC-01L, SCV40, SM-N9600, SM-N960F, SM-N960N, SM-N960U, SM-N960U1, SM-N960W)	10
	▪ Galaxy Note 10 (SCV45, SM-N9700, SM-N970F, SM-N970U, SM-N970U1, SM-N970W)	10, 11
	▪ Galaxy Note 10 5G (SM-N971N, SM-N971U)	10, 11
	▪ Galaxy Note 10 + (SC-01M, SCV45, SM-N9750, SM-N975C, SM-N975F, SM-N975U, SM-N975U1, SM-N975W)	10, 11
	▪ Galaxy Note 10 + 5G (SM-N9760, SM-N976B, SM-N976Q, SM-N976U, SM-N976V)	
	▪ Galaxy Note 20 (SM-N980F)	10, 11
	▪ Galaxy Note 20 5G (SM-N9810, SM-N981B, SM-N981N, SM-N981U, SM-N981U1, SM-N981W)	10, 11
	▪ Galaxy Note 20 Ultra (SM-N985F)	10, 11
	▪ Galaxy Note 20 Ultra 5G (SC-53A, SCG06, SM-N9860, SM-N986B, SM-N986N, SM-N986U, SM-N986U1, SM-N986W)	10, 11
LG™	▪ K20 (LM-X120EMW), G6	9
	▪ G8 ThinQ, G7 ThinQ, V40 Thin Q, G8 ThinQ, V50 ThinQ, V50S ThinQ, V60 ThinQ 5G	10
Sony™	▪ Xperia 1 (802SO, J8110, J8170, J9110, J9180, SO-03L, SOV40)	10
	▪ Xperia 5 (901SO, J8210, J8270, J9210, J9260, SO-01M, SOV41)	
Google	▪ Pixel 3, Pixel 3 XL, Pixel 3A, Pixel 3A XL	10, 11
	▪ Pixel 4, Pixel 4 XL, Pixel 4A, Pixel 4A 5G	10, 11
	▪ Pixel 5	11
Huawei	▪ Mate Pro 20	9, 10
Xiaomi	▪ Mi 9	9, 10

Finding your model number on your phone:

ANDROID: Settings > About Phone

APPLE MOBILE DEVICE: Settings > General > About