



DESIGNED WITH PURPOSE

With a long history of innovation and a reputation for quality, we design products with your critical needs in mind.

Concerto™ detachable coil systems are designed to help you easily and precisely get to the embolization area and instantly **detach the coil on target.**¹

99.7%
DETACHMENT
SUCCESS²

PROVEN
RELIABILITY





Concerto™
Detachables Coil System

Medtronic
Further, Together

THE UNIQUE CONCERTO™ COIL DESIGN

- Ball-and-socket design allows the coil implant to move independently of the pusher in a 360° fashion, resulting in **less energy buildup and microcatheter kickout.**¹
- Bevel retainer design allows for **articulation zone flexibility.**¹
- Soft distal pusher **minimizes catheter kickback.**¹



¹ Medtronic data on file.

² Medtronic data on file. Bench testing conducted at 95% confidence, results may not be indicative of clinical performance.

Important Information: Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device.

Indications for Use: The Concerto™ detachable coil system is indicated for arterial and venous embolizations in the peripheral vasculature. Results may vary. Not all patients achieve the same results.

CAUTION: Federal (USA) law restricts these devices to sale, distribution, and use by or on the order of a physician.

[medtronic.com/embolization](https://www.medtronic.com/embolization)

UC202001818 EN ©2019 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. For distribution in the U.S. only. 10/2019

Medtronic