

DESIGNED WITH PURPOSE

With a long history of innovation and a reputation for quality, we design products with your critical needs in mind.


Concerto™ detachable coil systems are designed to access small or large vessels, position precisely, and detach when you decide. You can achieve mechanical occlusion and support biological response with our fibered coils.

Achieve occlusion with our fibered coils.




Your goal is to occlude — quickly and completely.

All Concerto™ coils are designed with nylon or PGLA fibers **to increase thrombogenicity and stop flow.**¹



Concerto™
Detachable Coil System



Concerto™ fibered coils:

ACCELERATE
THROMBUS FORMATION

241%

HIGHER RATE OF
THROMBIN
PRODUCTION^{1*}

NYLON/PGLA
FIBERED COILS

BARE
COILS

MAXIMIZE
THROMBUS FORMATION

80%

INCREASE IN PEAK
THROMBIN
PRODUCTION^{1*}

NYLON/PGLA
FIBERED COILS


BARE
COILS

Medtronic
Further, Together

Concerto™ coils
are available in both helix
and complex shapes,
in a broad range of sizes
from 2 to 20 mm and
lengths from 2 to 50 cm.



Concerto™
Helix Detachable Coil System



Concerto™
3D Detachable Coil System

Feel confident that
Concerto™ coils will
perform for you
and your patients.

REFERENCES

*Percentages are calculated from the average of nylon and PGLA fibered coils.

¹ Girdhar G, Read M, Sohn J, Shah C, Shrivastava S. In-vitro thrombogenicity assessment of polymer filament modified and native platinum embolic coils. *J Neurol Sci.* April 15, 2014;339(1-2):97-101.

The Concerto™ detachable coil system is indicated for arterial and venous embolizations in the peripheral vasculature. Results may vary. Not all patients achieve the same results.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

CAUTION: Federal Law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

Medtronic

Medtronic
9775 Toledo Way
Irvine, CA 92618
USA
Tel: 877.526.7890
Fax: 763.526.7888

medtronic.com/dontgobare

UC201901719b EN ©2018 Medtronic. All rights reserved.

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. For distribution in the United States only. 10/18