

CARDIOINSIGHT™ WORKSTATION

Product specifications

Ordering Information

Order Number

CIT200	CardioInsight™ Workstation
--------	----------------------------

Workstation Specifications

Size

Height	MIN: 1,300 mm (51.18 in) MAX: 1,461 mm (57.52 in)
--------	--

Depth	MIN: 664 mm (26.14 in) MAX: 862 mm (33.94 in)
-------	--

Width	MIN: 615 mm (24.21 in) MAX: 696 mm (27.40 in)
-------	--

Weight	102 kg (225 lbs)
--------	------------------

Other Specs

Power Input Parameters-System	100/120/220/240 V AC
----------------------------------	----------------------

Power Input Parameters-Amplifier	100 - 240 V AC
-------------------------------------	----------------

Temperature and Humidity Limitations	15 °C to 30 °C (59 °F to 86 °F) Humidity 20%-80%
---	---

Workstation Non-Sterile Components

Description

Mapping Amplifier

Workstation Computer

Workstation Monitor

Workstation Cart

Sensor Array Signal Cables (Right, Left, Back right, Back left)

Data Cable (3', 10', 25', 50', 75', 100')

Patient Ground and Reference Cables

Isolation Transformer

Display Cables

Power Cables



The CardioInsight™ noninvasive 3D mapping system is intended for gathering, analysis, display, and storage of cardiac electrophysiological data and maps for analysis by a physician.

The CardioInsight™ workstation is used with the CardioInsight™ mapping vest. Together, they function as a noninvasive mapping system that collects body surface ECG signals and combines these signals with CT scan data to produce and display simultaneous bi-atrial and bi-ventricular, three-dimensional cardiac maps.

The CardioInsight™ software provides various cardiac signal analyses and displays interactive three-dimensional color maps.

CardioInsight™ Noninvasive 3D Mapping System



Essential Product Information

Indications and relevant warnings, and instructions for use can be found in the product labeling supplied with each device.

Brief Statement

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative or consult the Medtronic website at www.medtronic.com.

Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
www.medtronic.eu
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

United Kingdom/Ireland
Medtronic Limited
Building 9
Croxley Park
Hatters Lane
Watford
Herts WD18 8WW
www.medtronic.co.uk
Tel: +44 (0)1923 212213
Fax: +44 (0)1923 241004

UC201604044aEE © Medtronic 2017.
All rights reserved. Printed in Europe.