

Medtronic

IPC™

Integrated Power Console



Integration
and versatility

The IPC™ System

Integration that brings versatility to your OR

The Integrated Power Console's design was inspired by feedback from you, the surgeons and nurses in the field. Our IPC™ system is a powered ENT console with an expansive range of application-specific products with the adaptability to run our future product pipeline.

Versatility in a single console

A broad selection of powered surgical blades, burs and saws, for a wide range of surgical applications/procedures

- Rhinology
- Facial plastics
- Nasal obstruction
- Transnasal procedures including anterior skull base and DCR
- Otology
- Neurotology
- Nasopharyngeal
- Laryngeal, tracheal and bronchial



Integrated lens cleaning system

Endo-Scrub™. Maintain your visualisation at the surgical site without removing the endoscope

- Six cycle settings to suit surgeon preference
- High-flow wash cycle
- Auto-prime designed to reduce drips or air bubbles from forming over the lens
- Additional 10 second reverse cycle to prevent dripping over lens
- A range of sheaths compatible with a variety of endoscopes from multiple manufacturers

Complete

Powering a variety of handpieces and Endo-Scrub™ 2

- Straightshot™ M5 Microdebrider
- Skeeter™ Otologic Drill
- Endo-Scrub™ 2 Lens Cleaning Sheaths
- ENTMR8 Midas Rex™ Otological Drill System



Customisable to surgeon preference

Ability to change default settings

- Speed and irrigation default settings can be customised to surgeon preference
- Continuous irrigation 0.5-100 cc/min
- Pulse flow irrigation (20 cc/minute and below)
- Rotatable blade can be set to:
stop and start in the same position
- 60° and 180°



Integrated facial nerve monitoring

NIM Vital Nerve Monitoring System

- Works complementary to ENTMR8 Otological Drill system in Otology and Neurotology procedures



Ordering Information

IPC™ System and Endo-Scrub™ 2 Accessories

Product	Description	Qty
1898001	Integrated Power Console	1
1898430	Integrated Power Console Multifunction Footpedal	1
1898600	IPC Cart	1
1852000	Footswitch Endo-Scrub 2	1

Straightshot™ M5 Microdebrider Components

Product	Description	Qty
1899200	Straightshot™ M5 Microdebrider	1
1899076	M5 Instrument Tray	1
3318520	Cleaning Cap Handpiece Universal	1

Skeeter™ Components

Product	Description	Qty
3155615	Skeeter Instrument Tray	1
3052460	Skeeter Otologic Metric Scale	1
3155641	Skeeter Oto-Flex Bur	1

Midas™ Rex Otological Drill and Components

Product	Description	Qty
1847DRLMTR	ENTMR8 Motor assembly	1
1847ATTAS	ENTMR8 Straight Attachment	1
1847ATTSS	ENTMR8 Angled Attachment	1
1847TRAY	ENTMR8 Instrument Tray	1

INSTRUCTIONS FOR USE.

Please refer to the product insert or Instructions for Use for a complete list of indications, contraindications, warnings, precautions, and other important medical information. Always follow the Instructions for Use.

Medtronic

Medtronic Australasia Pty Ltd

2 Alma Road, Macquarie Park, NSW 2113 Australia
Toll Free: 1800 668 670 Fax: +61 2 9889 5167

Medtronic New Zealand Ltd

Level 3 - Building 5, Central Park Corporate Centre
666 Great South Road, Penrose, Auckland 1051 New Zealand
Fax: +64 9 918 3742 Toll Free: 0800 377 807

© Medtronic 2022 All Rights Reserved
ENT 006-08-22 ANZ. #11840-112022

medtronic.com.au
medtronic.co.nz