

# SCREENING QUESTIONS

MyCareLink Smart Patient Monitor is the first app-based remote monitoring system for patients with Medtronic implanted pacemakers, including CRT-P.

To see if you qualify for the MyCareLink Smart Monitor, please answer the following two questions.

1. Do you or your care supporter have access to an Apple® or Android™ smartphone or tablet, or plan to in the immediate future?

Yes  No

2. Was the smartphone or tablet you plan to use purchased new within the last two years?

Yes  No

3. If you answered “yes” to questions 1 and 2, go to [mycarelinksmart.com](http://mycarelinksmart.com) to download the app. If you're able to download the app, it means your smartphone or tablet is compatible and could be used with the MyCareLink Smart Monitor.

MyCareLink Smart™  
Patient Monitor



The MyCareLink Smart Monitor is compatible with the following:

Android: 4.1.2 and above

Apple iPhone®: 4S and above

Apple iPad®: 3 and above

Apple iPad mini™: all versions

Apple iPad Air®: all versions

BOTH Apple iPhone and iPad: iOS 7.1 and above

The MyCareLink Smart Monitor combines a hand-held Reader that communicates to your smartphone or tablet using a free app, called the MyCareLink Smart™ Application (App). Once you have the Reader—which is prescribed by your doctor—you simply download the app onto your smartphone or tablet and pair the Reader with the app.

Medtronic

## **Brief Statement**

For all products referenced please see the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential complications/adverse events. Please contact your local representative for further information.

**Medtronic of Canada, Ltd.**  
99 Hereford Street  
Brampton, Ontario L6Y 0R3  
Canada  
Tel: (905) 460-3800  
Fax: (905) 826-6620  
Toll-free: 1 (800) 268-5346

**Patient Services:**  
1 (888) 660-4616  
Monday-Friday  
9:00 am-5:00 pm (EST)