


Medtronic

Evolut™ FX+ TAVI System

Product specifications

Bioprosthesis




Evolut FX+ 23 mm Valve

Model number:
EVFXPLUS-23

Size: 23 mm

Aortic annulus
diameter:
17†/18–20 mm




Evolut FX+ 26 mm Valve

Model number:
EVFXPLUS-26

Size: 26 mm

Aortic annulus
diameter:
20–23 mm



Evolut FX+ 29 mm Valve

Model number:
EVFXPLUS-29

Size: 29 mm

Aortic annulus
diameter:
23–26 mm




Evolut FX+ 34 mm Valve

Model number:
EVFXPLUS-34

Size: 34 mm

Aortic annulus
diameter:
26–30 mm


Loading system



Evolut FX Loading System

Model number:
L-EVOLUTFX-2329

23, 26, and 29 mm TAV



Evolut FX Loading System

Model number:
L-EVOLUTFX-34

34 mm TAV


†Measurement for TAV-in-SAV only.

All items above are single-use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

Delivery catheter system


Evolut FX delivery catheter system with InLine sheath 14 Fr equivalent (6.0 mm capsule outer diameter)

Model number:
D-EVOLUTFX-2329
23, 26, and 29 mm TAV



Evolut FX delivery catheter system with InLine sheath 18 Fr equivalent (7.33 mm capsule outer diameter)

Model number:
D-EVOLUTFX-34
34 mm TAV



Medtronic

99 Hereford Street
Brampton, Ontario, L6Y 0R3
Toll-free: 800.268.5346
Tel: 905.460.3800

medtronic.ca

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

©2025 Medtronic. All rights reserved. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic. TM*Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
UC202402378aEC CA-CTL-0503-E Rev. 03/2025