

Medtronic infusion sets and reservoirs

A set for every body. A set for every lifestyle.

	Medtronic Extended infusion set	MiniMed™ Mio™ Advance	MiniMed™ Quick-Set™	MiniMed™ Silhouette™	MiniMed™ Sure-T™
	Extended reservoir	Standard reservoir			
					
Infusion set design	<ul style="list-style-type: none"> Fewer set changes Dexterity issues Needle anxiety 	<ul style="list-style-type: none"> Dexterity issues Needle anxiety 	<ul style="list-style-type: none"> Alternative site usage Dexterity issues 	<ul style="list-style-type: none"> Lean/active 	<ul style="list-style-type: none"> History of bent cannula History of scar tissue
Recommended body type	Medium/High BMI*	Medium/High BMI	Medium/High BMI	Low BMI	Medium/High BMI
Insertion device	All in one	All in one	Quick-serter™	MiniMed™ Sil-serter™ or Manual insertion	Manual insertion
Cannula/needle type	Soft	Soft	Soft	Soft	Steel
Insertion angle	90°	90°	90°	30°- 45°	90°
Cannula length	6 mm, 9 mm	6 mm, 9 mm	6 mm, 9 mm	13 mm, 17 mm	6 mm, 8 mm
Cannula fill	0.600 units	0.600 units	0.300 units, 0.500 units	0.700 units	No cannula fill
Tubing length	23 in, 32 in	23 in, 43 in	18 in, 23 in, 32 in, 43 in	18 in, 23 in, 32 in, 43 in	18 in, 23 in, 32 in
Change frequency	6-7 days	2-3 days	2-3 days	2-3 days	1-2 days

*Body Mass Index

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the respective infusion set and reservoir Instructions for Use. Refer to your Medtronic insulin pump user guide for a list of compatible insulins and infusion sets.