

# Device and Software Compatibility: MiniMed™ 700-Series System Apps

**Medtronic**

Two smart device apps are available with the MiniMed™ 770G system and MiniMed™ 780G system.

The following devices are currently compatible with our apps. Please note, compatibility may vary between the two apps. We will continue to update this page as we gain new compatibility for an operating system (OS).

We are constantly testing new devices and operating systems (OS) to ensure they will work properly with our apps. As we gain compatibility with more devices and OS, we will continue to update this page. Devices from selected manufacturers may be supported, but they have not been fully tested. Please note, compatibility may vary among the apps.



**MiniMed™ Mobile App allows those wearing the system to view their data on a smart device.**

| MiniMed™ Mobile app 2.1<br>Apple Devices and Operating Systems (iOS)        |   |   |   |   |
|---|---|---|---|---|
| Operating system  | 14, 15 and 16   |   |   |   |
| Compatible Devices  | iPhone SE (2020)<br>iPod Touch (7th generation)<br>iPhone 6S - only iOS 14<br>iPhone 6S Plus - only iOS 14<br>iPhone 7<br>iPhone 7 Plus<br>iPhone 8<br>iPhone 8 Plus<br>iPhone X<br>iPhone XS<br>iPhone XS Max<br>iPhone XR<br>iPhone 11<br>iPhone 11 Pro<br>iPhone 11 Pro Max<br>iPhone 12<br>iPhone 12 Mini |   | iPhone 12 Pro<br>iPhone 12 Pro Max<br>iPhone 13<br>iPhone 13 Mini<br>iPhone 13 Pro<br>iPhone 13 Pro Max<br>iPhone 14<br>iPhone 14 Plus<br>iPhone 14 Pro<br>iPhone 14 Pro Max<br>Apple Watch Series 4<br>Apple Watch Series 5<br>Apple Watch Series 6<br>Apple Watch Series 7<br>Apple Watch Series 8<br>Apple Watch Ultra<br>Apple Watch Series SE (2020)<br>Apple Watch Series SE 2nd Gen (2022) |   |
|   | Compatible Apple Watch Faces  | Activity Analog<br>Activity Digital<br>Astronomy<br>Breathe<br>Explorer | Fire and Water<br>Kaleidoscope<br>Liquid Metal<br>Mickey Mouse<br>Motion  | Nike Analog<br>Photos<br>Portraits<br>Solar Graph |
| Medtronic will continue to work to bring support to additional watch faces. |   |   |   |   |

**MiniMed™ Mobile app 2.1**  
Android Devices and Operating Systems (OS)

| Manufacturer | Model   | Android OS     |
|--------------|---|----------------|
| Samsung      | <ul style="list-style-type: none"> <li>• S9 (SC-02K, SCV38, SM-G960F, SM-G960N, SM-G9600, SM-G9608, SM-G960W, SM-G960U, SM-G960U1)</li> <li>• S9+ (SC-03K, SCV39, SM-G965F, SM-G965N, SM-G9650, SM-G965W, SM-G965U, SM-G965U1)</li> </ul>   | 9, 10          |
| Samsung      | <ul style="list-style-type: none"> <li>• Galaxy A8 (2018) (SM-A530F, SM-A530K, SM-A530L, SM-A530S, SM-A530N, SM-A530W, SM-A530X)</li> </ul>   | 9              |
|              | <ul style="list-style-type: none"> <li>• Galaxy A10e (SM-A102*, SM-S102DL)</li> <li>• Galaxy A20 (SM-A205*, SM-S205DL, SC-02M, SCV46-J, SCV46-U, SCV46)</li> <li>• Galaxy A7 (2018) (SM-A750F, SM-A750FN, SM-A750G, SM-A750GN, SM-A750C, SM-A750X, SM-A750N)</li> <li>• Galaxy A20e (SM-A202F, SM-A202K)</li> <li>• Galaxy A21s (SM-A217F, SM-A217F/DS, SM-A217F/DSN, SM-A217M, SM-A217M/DS, SM-A217N)</li> </ul>   | 10             |
|              | <ul style="list-style-type: none"> <li>• Galaxy A52 5G (SM-A526*)</li> <li>• Galaxy A52s 5G (SM-A528*)</li> </ul>   | 11, 12         |
| Samsung      | <ul style="list-style-type: none"> <li>• S10 (SC-03L, SCV41, SM-G973F, SM-G973N, SM-G9730, SM-G9738, SM-G973C, SM-G973U, SM-G973U1, SM-G973W)</li> <li>• S10 lite (SM-G770*) - only supports OS 10 and OS 11</li> <li>• S10 5G (SM-G977B, SM-G977N, SM-G977P, SM-G977T, SM-G977U)</li> <li>• S10e (SM-G970F, SM-G970N, SM-G9700, SM-G9708, SM-G970U, SM-G970U1, SM-G970W)</li> <li>• S10+ (SC-04L, SCV42, SM-G975F, SM-G975N, SM-G9750, SM-G9758, SM-G975U, SM-G975U1, SM-G975W)</li> </ul>                                   | 9, 10, 11, 12  |
| Samsung      | <ul style="list-style-type: none"> <li>• S20 (SM-G980F)</li> <li>• S20+ (SM-G985F)</li> <li>• S20 5G (SC-51A, SCG01, SM-G9810, SM-G981N, SM-G981U, SM-G981U1, SM-G981V, SM-G981W, SM-G981B)</li> <li>• S20+ 5G (SC-52A, SCG02, SM-G9860, SM-G986N, SM-G986U, SM-G986U1, SM-G986W, SM-G986B)</li> <li>• S20 Ultra 5G (SCG03, SM-G9880, SM-G988N, SM-G988U, SM-G988U1, SM-G988W, SM-G988B)</li> <li>• S20 FE (SM-G780*)</li> <li>• S20 FE 5G (SM-G7810, SM-G781B, SM-G781N, SM-G781U, SM-G781U1, SM-G781V, SM-G781W)</li> </ul> | 10, 11, 12, 13 |
| Samsung      | <ul style="list-style-type: none"> <li>• Samsung Galaxy S21 (SM-G990*)</li> <li>• Samsung Galaxy S21 5G (SM-G991*)</li> <li>• Samsung Galaxy S21+ (5G SM-G996*)</li> <li>• Samsung Galaxy S21 Ultra (5G SM-G998*)</li> </ul>  | 11, 12         |
|              | <ul style="list-style-type: none"> <li>• Samsung Galaxy S22 5G (SM-S901*)</li> <li>• Samsung Galaxy S22+ 5G (SM-S906*)</li> <li>• Samsung Galaxy S22 Ultra 5G (SM-S908*)</li> </ul>   | 12, 13         |

|  |  |                |
|--|--|----------------|
| <b>Samsung</b>   | • Galaxy Note 9 (SC-01L, SCV40, SM-N9600, SM-N960F, SM-N960N, SM-N960U, SM-N960U1, SM-N960W)                                   | 10             |
|  | • Galaxy Note 10 (SCV45, SM-N9700, SM-N970F, SM-N970U, SM-N970U1, SM-N970W)  | 10, 11, 12     |
|  | • Galaxy Note 10 5G (SM-N971N, SM-N971U)   |                |
|  | • Galaxy Note 10 + (SC-01M, SCV45, SM-N9750, SM-N975C, SM-N975F, SM-N975U, SM-N975U1, SM-N975W)                                |                |
|  | • Galaxy Note 10 + 5G (SM-N9760, SM-N976B, SM-N976Q, SM-N976U, SM-N976V)   |                |
|  | • Galaxy Note 10 Lite (SM-N770F, SM-N770F/DS, SM-N770F/DSM)  |                |
|  | • Galaxy Note 20 (SM-N980F)  |                |
|  | • Galaxy Note 20 5G (SM-N9810, SM-N981B, SM-N981N, SM-N981U, SM-N981U1, SM-N981W)  |                |
|  | • Galaxy Note 20 Ultra (SM-N985F)  |                |
| • Galaxy Note 20 Ultra 5G (SC-53A, SCG06, SM-N9860, SM-N986B, SM-N986N, SM-N986U, SM-N986U1, SM-N986W) |  |                |
| <b>Motorola</b>  | • Moto G6 (XT1925, XT1925-10, XT1925-13, XT1925-6, XT1925-5, XT1925-3, XT1925-7, moto g(6))                                    | 9              |
|  | • Moto G7 (XT1962, XT1962-4, moto g(7)) - only supports OS 9   | 9, 10          |
|  | • Moto G7 Power (XT1955, moto g(7) power)  |                |
|  | • Moto G7 Play (XT1952-2, XT1952-3, XT1952-4, moto g(7) play)  |                |
| <b>One Plus</b>  | • OnePlus 6T (ONEPLUS A6010, ONEPLUS A6013)  | 9, 10          |
|  | • OnePlus 9 (LE2110)   | 12             |
|  | • One Plus 10 Pro (NE2210)   |                |
| <b>LG</b>  | • G7 ThinQ, V40 Thin Q, G8 ThinQ, V50 ThinQ, V50S ThinQ, V60 ThinQ 5G  | 10             |
| <b>Google</b>  | • Pixel 3, Pixel 3A, Pixel 3 XL, Pixel 3A XL   | 10, 11, 12     |
|  | • Pixel 4, Pixel 4 XL, Pixel 4A, Pixel 4A 5G   | 10, 11, 12, 13 |
|  | • Pixel 5, Pixel 5A, Pixel 5A 5G   | 11, 12, 13     |
| <b>Huawei</b>  | • Mate Pro 20  | 9, 10          |
| <b>Sony™</b>   | • Xperia 1 (80250, J8110, J8170, J9110, J9180, SO-03L, SOV40)<br>• Xperia 5 (90150, J8210, J8270, J9210, J9260, SO-01M, SOV41) | 10             |
| <b>Xiaomi</b>  | • Mi 6   | 9, 10          |
|  | • Xiaomi 12 (2201123)  | 12             |
|  | • Xiaomi 12 Pro (2201122)  |                |



## CareLink™ Connect App allows care partners to follow someone wearing the system.

| CareLink™ Connect app 3.1<br>Apple Devices and Operating Systems (iOS) |  |
|--|--|
| Operating system   | 14, 15 and 16  |
| Compatible Phones  | <ul style="list-style-type: none"> <li>iPhone SE (2020)</li> <li>iPhone 6S</li> <li>iPhone 6S Plus</li> <li>iPhone 7</li> <li>iPhone 7 Plus</li> <li>iPhone 8</li> <li>iPhone 8 Plus</li> <li>iPhone X</li> <li>iPhone XS</li> <li>iPhone XS Max</li> <li>iPhone XR</li> <li>iPhone 11- only OS 13 +</li> <li>iPhone 11 Pro - only OS 13 +</li> <li>iPhone 11 Pro Max - only OS 13 +</li> <li>iPhone 12 - only OS 14 +</li> <li>iPhone 12 Mini - only OS 14 +</li> <li>iPhone 12 Pro - only OS 14 +</li> <li>iPhone 12 Pro Max - only OS 14 +</li> <li>iPhone 13</li> <li>iPhone 13 Mini</li> <li>iPhone 13 Pro</li> <li>iPhone 13 Pro Max</li> <li>iPhone 14</li> <li>iPhone 14 Plus</li> <li>iPhone 14 Pro</li> <li>iPhone 14 Pro Max</li> </ul> |

| CareLink™ Connect app 3.1<br>Android Devices and Operating Systems (OS) |   |                |
|---|---|----------------|
| Manufacturer  | Model   | Android OS     |
| Samsung   | <ul style="list-style-type: none"> <li>• S9 (SC-02K, SCV38, SM-G960F, SM-G960N, SM-G9600, SM-G9608, SM-G960W, SM-G960U, SM-G960U1)</li> <li>• S9+ (SC-03K, SCV39, SM-G965F, SM-G965N, SM-G9650, SM-G965W, SM-G965U, SM-G965U1)</li> </ul>   | 9, 10          |
| Samsung   | <ul style="list-style-type: none"> <li>• S10 (SC-03L, SCV41, SM-G973F, SM-G973N, SM-G9730, SM-G9738, SM-G973C, SM-G973U, SM-G973U1, SM-G973W)</li> <li>• S10 5G (SM-G977B, SM-G977N, SM-G977P, SM-G977T, SM-G977U)</li> <li>• S10+ (SSC-04L, SCV42, SM-G975F, SM-G975N, SM-G9750, SM-G9758, SM-G975U, SM-G975U1, SM-G975W)</li> <li>• S10e (SM-G970F, SM-G970N, SM-G9700, SM-G9708, SM-G970U, SM-G970U1, SM-G970W)</li> <li>• S10 Lite (SM-G770F, SM-G770U1) - only supports OS 10 and OS 11</li> </ul>                       | 9, 10, 11      |
| Samsung   | <ul style="list-style-type: none"> <li>• S20 (SM-G980F)</li> <li>• S20+ (SM-G985F)</li> <li>• S20 5G (SC-51A, SCG01, SM-G9810, SM-G981N, SM-G981U, SM-G981U1, SM-G981V, SM-G981W, SM-G981B)</li> <li>• S20+ 5G (SC-52A, SCG02, SM-G9860, SM-G986N, SM-G986U, SM-G986U1, SM-G986W, SM-G986B)</li> <li>• S20 Ultra 5G (SCG03, SM-G9880, SM-G988N, SM-G988U, SM-G988U1, SM-G988W, SM-G988B)</li> <li>• S20 FE (SM-G780F)</li> <li>• S20 FE 5G (SM-G7810, SM-G781B, SM-G781N, SM-G781U, SM-G781U1, SM-G781V, SM-G781W)</li> </ul> | 10, 11, 12, 13 |

|   |   |                |
|---|---|----------------|
| Samsung   | <ul style="list-style-type: none"> <li>• S21 5G (SM-G9910, SM-G991U1, SM-G991W, SM-G991B, SM-G991N)</li> <li>• S21 Ultra 5G (SM-G9980, SM-G998U, SM-G998U1, SM-G998W, SM-G998B, SM-G998N)</li> <li>• S21+ 5G (SM-G9960, SM-G996U1, SM-G996W, SM-G996B, SM-G996N)</li> </ul> | 11, 12, 13     |
|   | <ul style="list-style-type: none"> <li>• Samsung Galaxy S22 5G (SM-S901*)</li> <li>• Samsung Galaxy S22+ 5G (SM-S906*)</li> <li>• Samsung Galaxy S22 Ultra 5G (SM-S908*)</li> </ul>   | 12, 13         |
| Samsung   | <ul style="list-style-type: none"> <li>• Galaxy Note 9 (SC-01L, SCV40, SM-N9600, SM-N960F, SM-N960N, SM-N960U, SM-N960U1, SM-N960W)</li> </ul>  | 10             |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 10 (SCV45, SM-N9700, SM-N970F, SM-N970U, SM-N970U1, SM-N970W)</li> </ul>   | 10, 11         |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 10 5G (SM-N971N)</li> </ul>  |                |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 10 + (SC-01M, SCV45, SM-N9750, SM-N975C, SM-N975F, SM-N975U, SM-N975U1, SM-N975W)</li> </ul>   |                |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 10 + 5G (SM-N9760, SM-N976B, SM-N976Q, SM-N976U, SM-N976V)</li> </ul>  |                |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 20 (SM-N980F)</li> </ul>   |                |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 20 5G (SM-N9810, SM-N981B, SM-N981N, SM-N981U, SM-N981U1, SM-N981W)</li> </ul>   |                |
|   | <ul style="list-style-type: none"> <li>• Galaxy Note 20 Ultra (SM-N985F)</li> </ul>   |                |
| <ul style="list-style-type: none"> <li>• Galaxy Note 20 Ultra 5G (SM-N9860, SM-N986B, SM-N986N, SM-N986U, SM-N986U1, SM-N986W)</li> </ul> |   |                |
| LG™   | <ul style="list-style-type: none"> <li>• K20 (LM-X120EMW), G6</li> </ul>  | 9              |
|   | <ul style="list-style-type: none"> <li>• G8 ThinQ, G7 ThinQ, V40 Thin Q, G8 ThinQ, V50 ThinQ, V50S ThinQ, V60 ThinQ 5G</li> </ul>   | 10             |
| Sony™   | <ul style="list-style-type: none"> <li>• Xperia 1 (802SO, J8110, J8170, J9110, J9180, SO-03L, SOV40)</li> <li>• Xperia 5 (901SO, J8210, J8270, J9210, J9260, SO-01M, SOV41)</li> </ul>  | 10             |
| Google  | <ul style="list-style-type: none"> <li>• Pixel 3, Pixel 3 XL, Pixel 3A, Pixel 3A XL</li> </ul>  | 10, 11, 12, 13 |
|   | <ul style="list-style-type: none"> <li>• Pixel 4, Pixel 4 XL, Pixel 4A, Pixel 4A 5G</li> </ul>  | 11, 12, 13     |
|   | <ul style="list-style-type: none"> <li>• Pixel 5, Pixel 5A</li> </ul>   |                |
| Huawei  | <ul style="list-style-type: none"> <li>• Mate Pro 20</li> </ul>   | 9, 10          |
| Xiaomi  | <ul style="list-style-type: none"> <li>• Mi 9</li> </ul>  | 9, 10          |
|   | <ul style="list-style-type: none"> <li>• Xiaomi 12 (2201123)</li> <li>• Xiaomi 12 Pro (2201122)</li> </ul>  | 12             |
|   | <ul style="list-style-type: none"> <li>• OnePlus 6T (A6010, A6013)</li> </ul>   | 10, 11         |
| One Plus  | <ul style="list-style-type: none"> <li>• OnePlus 9 (LE2110)</li> <li>• OnePlus 10 Pro (NE2210)</li> </ul>   | 12             |

## Finding your model number on your phone:

Android: Settings > About Phone

Apple Mobile Device: Settings > General > About