

# Medtronic

Engineering the extraordinary

Sonicision™ 7 curved jaw cordless ultrasonic dissector

## Cordless freedom. Less waste.<sup>1,†</sup>



The reusable, autoclavable generator with integrated transducer works for

# 150 procedures

before you need to replace it.<sup>2</sup>



The Harmonic Ace™\* +7 Shears with Advanced Hemostasis device's transducer only works

## 95 times

before needing replacement.<sup>3</sup>



The Harmonic™\* HD1000i Shears and 1100 Shears instruments and their integrated transducers are

## single-use –

requiring replacement after every procedure.<sup>4,5</sup>



The Thunderbeat™\* Front-Actuated Grip Type S device's transducer only works

## 100 times

before needing replacement.<sup>6</sup>

†Mass comparison made with dissector and assembly tool.

# The freedom to do more<sup>7,†,‡</sup> — while throwing away less.<sup>1,§</sup>



The reusable battery works for

200  
procedures

before requiring replacement.<sup>8,9</sup>

Recycling programs available in some regions



**Contact your Medtronic sales rep today.**

©2024 Medtronic. All rights reserved.  
Medtronic, Medtronic logo and  
Engineering the extraordinary are  
trademarks of Medtronic. ™\*Third  
party brands are trademarks of their  
respective owners. All other brands are  
trademarks of a Medtronic company.  
CA-OR-0023-E Rev. 01/2024

99 Hereford Street  
Brampton, Ontario, L6Y 0R3  
Toll-free: 800.268.5346  
Tel: 905.460.3800  
**medtronic.ca**

**Medtronic**

†25 out of 27 surgeons surveyed after use agreed.

‡26 out of 27 surgeons surveyed after use agreed. Compared to the Sonicision™ curved jaw device.

§Mass comparison made with dissector and assembly tool. Sonicision™ 7 generator and battery not included.

#### References:

1. Based on internal test report #RE00439207 rev A, Soni comparison report. December 8, 2022.
2. Based on internal test report #RE00108733 rev A, Reusable life test verification report, Project Sabre. August 2017.
3. Harmonic™\* Gray Hand Piece and Gray Test Tip [instructions for use]. Guaynabo, PR: Ethicon Endo-Surgery; 2020.
4. Harmonic™\* 1100 Shears [value analysis summary]. Ethicon; 2022.
5. Harmonic™\* HD1000i Shears [sales brochure]. Cincinnati, OH: Ethicon Endo-Surgery, Inc.; 2018.
6. Thunderbeat™\* System Reference Guide [summary of instructions for use]. Southborough, MA: Olympus America Inc.; 2014.
7. Based on internal test report #RE00329878 rev A, Marketing evaluation of surgeon experience using the Sonicision™ 7 curved jaw cordless ultrasonic dissector. Apr. 14-15 and 20-22, 2021.
8. Based on internal test report #RE00113195 rev A, Reusable life test verification report for the SCBA reusable battery pack – Oxivir. September 13, 2017.
9. Based on internal test report #R0042707 rev A, Reusable life test verification report for the SCBA reusable battery pack. April 26, 2017.