

Medtronic

ClearView™ Technology

MR8™ Midas Rex ClearView™



**Advancing surgical
possibilities**

MR8™


Variety of shapes, lengths,
and cutting geometries
for various applications

Visibility and
protected shaft

Extended, self-irrigating sleeve protects surrounding
tissue from minimally exposed rotating shaft

Spine solutions

- Improved surgical visibility and access for minimally invasive surgeries
- Navigate select ClearView™ tools with the StealthStation™ S8 surgical navigation system



Midas Rex ClearView™

Low-profile integrated irrigation that minimizes obstructions while delivering irrigation directly near the cutting site



Direct connect base for one-step assembly

Cranial and ENT solutions



- Lateral Skull Base tools now available ~ 5 cm and 7 cm
- Improved surgical visibility and access for skull base and endoscopic surgery
- Navigate select transnasal and lateral skull base tools in cranial applications with the StealthStation™ S8 system



Ordering information: Spine

12cm length, proximal bend	12cm length, distal bend	14cm length, proximal bend
MR8-SP12BA30	MR8-SD12BA30	MR8-SP14BA30
MR8-SP12BA40	MR8-SD12BA40	MR8-SP14BA40
MR8-SP12MH25	MR8-SD12MH25	MR8-SP14MH25
MR8-SP12MH25D	MR8-SD12MH25D	MR8-SP14MH25D
MR8-SP12MH25T	MR8-SD12MH25T	MR8-SP14MH25T
MR8-SP12MH30	MR8-SD12MH30	MR8-SP14MH30
MR8-SP12MH30D	MR8-SD12MH30D	MR8-SP14MH30D
MR8-SP12MH30T*	MR8-SD12MH30T	MR8-SP14MH30T*

Transnasal

13cm length
MR8-ND13MH30D
MR8-ND13MH30
MR8-ND13BA30DC*
MR8-ND13BA30
MR8-ND13BA40DC*
MR8-ND13BA40
MR8-ND13BA45DC

* Indicates Stealth-Midas™ MR8™ compatibility

Lateral skull base

5cm length, proximal bend	7cm length, proximal bend
MR8-LP05BA05D	MR8-LP07BA10
MR8-LP05BA08D	MR8-LP07BA10D*
MR8-LP05BA10	MR8-LP07BA15
MR8-LP05BA10D*	MR8-LP07BA20
MR8-LP05BA12D	MR8-LP07BA20D
MR8-LP05BA15	MR8-LP07BA20DC
MR8-LP05BA15D	MR8-LP07BA30
MR8-LP05BA20	MR8-LP07BA30D
MR8-LP05BA20D	MR8-LP07BA30DC
MR8-LP05BA20DC	MR8-LP07BA40
MR8-LP05BA30	MR8-LP07BA40D
MR8-LP05BA30D	MR8-LP07BA40DC
MR8-LP05BA30DC	
MR8-LP05BA40	
MR8-LP05BA40D	
MR8-LP05BA40DC	

Irrigation

IRD870	TUBING IRD870 MR8 IRIGTN CV LOW-PRFL EA
IRD875	TUBING IRD875 MR8 IRIGTN CV LOW-PRFL 5PK

Medtronic

This material is for Healthcare Professionals in countries with applicable health authority product registrations.

Important: Always refer to the Instructions For Use (IFU) packaged with the product/e-IFU for complete instructions, indications, contraindications, warnings, and precautions.

© 2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

UC202003024cEE - cst-mr8-clearview-brochure.pdf - 9578123 - 06/2023

cema.medtronic.com