


**CONFIDENTLY**  
CROSS



**Intuition™**  
Advanced Workhorse Guidewire



# Intuition™

Advanced Workhorse Guidewire

The Intuition™ guidewire's torque performance and tip weight give you precision and confidence as you control and direct the wire.

- Improved torque control
- Designed to provide steerability and control in tortuous anatomy

## ACCU-CORE TECHNOLOGY OFFERS YOU:


### TORQUE PERFORMANCE

Our core wire process is designed to provide torque performance to control and direct the wire.

### TACTILE FEEL AND CROSSABILITY

The rounded tip and body provide control.

#### ACCU-CORE TECHNOLOGY




Centered core wire

#### PREVIOUS DESIGN



Uncentered core wire

## Intuition™ Advanced Workhorse Guidewire



### ADDITIONAL FEATURES

- Femoral (proximal) markers on wire, 90 cm and 100 cm from distal tip
- DOC™\* compatible proximal ends on 180-cm wires

A **smooth**, responsive feel you can **depend** on

# Intuition™ Advanced Workhorse Guidewire

**A choice of two coatings to meet your needs** — Available in Pro/Pel™ silicone coating or Hydro-Track™ hydrophilic coating, for a more lubricious wire without increasing tip stiffness.

**Ordering Information** — Contact your Medtronic sales representative.

Product	Product Code	Length (cm)	Coating	Tip
Intuition™ Guidewire	INTU180HS	180	Hydro-Track™	Straight tip
	INTU180SS	180	Pro/Pel™	Straight tip
	INTU300HS	300	Hydro-Track™	Straight tip
	INTU300SS	300	Pro/Pel™	Straight tip

**Note:** All guidewires are 0.014-in wires and 180-cm wires are DOC™\* compatible.

## Medtronic

**Medtronic**  
3576 Unocal Place  
Santa Rosa, CA 95403  
USA  
Tel: 707.525.0111

**LifeLine Customer Support**  
Tel: 877.526.7890  
Tel: 763.526.7890

**Product Services**  
Tel: 888.283.7868  
Fax: 800.838.3103

**Europe**  
**Medtronic Intl. Trading SARL**  
Route du Molliat 31  
Case Postale  
CH-1131 Tolochenaz  
Switzerland  
Tel: 41.21.802.7000  
Fax: 41.21.802.7900

**Canada**  
**Medtronic of Canada Ltd.**  
99 Hereford Street  
Brampton, Ontario L6Y 0R3  
Canada  
Tel: 905.460.3800  
Fax: 905.460.3998  
Toll-free: 800.268.5346

**Asia Pacific**  
**Medtronic Intl. Ltd.**  
49 Changi South Avenue 2  
Nasaco Tech Centre  
Singapore 486056  
Singapore  
Tel: 65.6436.5000  
Fax: 65.6776.6335

**Latin America**  
**Medtronic**  
9850 Doral Blvd., 4th Floor  
Miami, FL 33178  
USA  
Tel: 786.709.4200  
Fax: 786.709.4244

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. ™\* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. For distribution only in markets where the Intuition™ crossing guidewire has been approved. Not for distribution in Japan or France. UC201707399IE 12/18