

Set your sights on gold.

Beacon™ fine needle fiducial (FNF) system


The Beacon™ FNF needle's solid gold markers are designed to help you diagnose with confidence. Its unique, functional design:

- Produces relatively small artifact on MR imaging to assist in visualization and patient treatment planning¹
- Provides a clear echogenic response²
- Visualizes distinctly across imaging modalities like CT, x-ray, and kV³


The solid-gold
fiducial markers
deliver clear
visualization
across CT, x-ray,
and kV³



Precise. Predictable. Preloaded.



Simplify your endoscopic ultrasound (EUS) procedures with a single solution. The Beacon™ FNF system's exchangeable design and two pre-loaded fiducial markers allow you to facilitate tissue sampling and fiducial placement in the same procedure,⁴ eliminating manual fiducial loading,² reducing overall procedure time,⁵ and lessening potential needle stick exposure for staff.⁶



Sharpen your diagnostic vision

Order code	Product description	Unit of measure
DSF-19-01	Beacon™ FNF pre-loaded - 19ga nitinol	Each
DSF-22-01	Beacon™ FNF pre-loaded - 22ga	Each
F-19-05	Beacon™ FNF individually packaged needles - 19ga nitinol	Box of 5
F-22-05	Beacon™ FNF individually packaged needles - 22ga	Box of 5

Explore the power of the Beacon™ FNF system



Contact your sales rep
to learn more or scan the QR code
to see the device in action on YouTube.



Visit Medtronic Academy
for training and educational resources

CAUTION

Federal law restricts this device to sale by or use under the order of a physician.

1. Shellock R & D Services, Inc. Evaluation of Magnetic Field Interactions, Heating, and Artifacts for the Covidien Gold Seed Fiducial. 2015. 2. Covidien. TR-20015 (01) Design Verification Report - FNF. 2015. 3. Covidien. TR-20019 (01) Design Verification Report - FNF. 2015. 4. Choi JH, Seo DW, Park DH, Lee SK, Kim MH. Fiducial placement for stereotactic body radiation therapy under only endoscopic ultrasonography guidance in pancreatic and hepatic malignancy: practical feasibility and safety. *Gut Liver*. 2014;8(1):88-93. doi:10.5009/gnl.2014.8.1.88. 5. Kothari S, Kothari T, Wang G, Gonda TA, Gress FG, Kaul V. EUS-guided fiducial placement for stereotactic body radiotherapy (SBRT): First global experience with a new 22-gauge 'pre-loaded' EUS fiducial needle. *Am J Gastroenterol*. 2015;110(Suppl 1):S102. 6. Covidien. TR-20018 (01) Design Verification Report - FNF. 2015.

©2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 07/2025 - MC-DPG-2500007 - [WF#14075294]

[medtronic.com/EUS](https://www.medtronic.com/EUS)

Medtronic