

**Medtronic**

Engineering the extraordinary

# ENT product catalog

2025 - 2026



# Medtronic

Engineering the extraordinary

ENT product catalog

# Directory

# Terms and conditions†

Toll free phone 800-874-5797

Toll free fax 800-678-3995

[medtronicent.com](http://medtronicent.com)

## Payment terms and prices

**Payment terms are net 30 days from date of invoice.** Prices do not include any applicable federal, state, or other taxes and will be added to the invoice unless an exemption certificate is supplied. Prices are subject to change at Medtronic ENT's discretion without notice.

## Shipments

FOB, Jacksonville, FL, USA. Freight is prepaid and added to invoice. Date of shipment is based on best estimates using commercially reasonable efforts to ship on schedule. Medtronic ENT shall not be liable for damage or loss due to any delay in delivering. Partial shipments will be made as required. Claims for short shipment are to be made within 10 days of shipment. Claims for goods damaged in shipment are to be made to the delivery carrier. Do not accept damaged shipments unless noted on the Bill of Lading at time of receipt. Medtronic ENT shall not be liable for any failure or delay in delivery occasioned in whole or in part by strikes, lockouts, fire, inability to obtain materials or shipping space, wars, civil commotion, destruction of plants, delay of carriers, any governmental act or regulation, or any other cause beyond its control.

## Warranty

All merchandise shall be covered by the Medtronic ENT standard limited warranty applicable to that product.

**Medtronic ENT specifically disclaims all other warranties, including, without limitation, those of merchantability and fitness for a particular purpose.**

## Product specifications and availability

Specifications, designs, prices, and availability are subject to change. Medtronic ENT reserves the right to discontinue merchandise or to change specifications without notice or obligation.

## Complaint handling

If you encounter a problem or complaint, please: 1) obtain all relevant facts, including name of user, address, phone number, merchandise number, lot/serial number, date of occurrence, and description of circumstances or nature of problem; and 2) report information immediately to Medtronic ENT, Customer Service at the above phone number. Every effort will be made to address and resolve problems in a manner consistent with our goal of providing quality merchandise and services.

## Returned goods

No merchandise may be returned without the prior authorization of Medtronic ENT. To request a returned goods authorization number (RGA) please contact Customer Service with the following information on the merchandise to be returned: 1) the original invoice number or purchase order number; 2) the reason for return; 3) the catalog number; 4) the quantity; and 5) the lot. Please display the RGA number prominently on the box and include it on all paperwork. All authorized returns must be sent freight prepaid.

- A. All returned merchandise must be in its original package. Any merchandise returned that is obsolete, beyond its expiration date, defaced, altered, damaged or has come in contact with body tissue or fluids will not be eligible for return or credit. Sterile merchandise must be returned in full-box quantities with the outer seals intact.
- B. All merchandise (except battery-powered devices) returned unused within 45 days of date of shipment will be credited at customer's purchase price. Battery-powered devices returned within 45 days of date of purchase will be charged a thirty percent (30%) restocking fee.  
All merchandise returned unused after 45 days and up to 90 days from date of shipment will be credited at customer's purchase price less a fifty percent (50%) restocking fee. Merchandise returned more than 90 days from date of shipment is not eligible for return or credit.
- C. Merchandise returned without prior authorization or not in accordance with A and B above will be returned to customer at customer's expense.

†Terms and conditions not applicable to the purchase of SINUVA™ (mometasone furoate) sinus implants. For terms and conditions related to that product, please contact your Medtronic representative.

# Terms and conditions<sup>†</sup>

## Return of product for repair

To arrange for the return of merchandise for repair, call Customer Service at 800-874-5797 to obtain a Return Goods Authorization number (RGA) and to schedule the repair. Please have the original invoice number or purchase order number available to assist in verifying warranty information.

A purchase order number is required to process all repairs. After obtaining an RGA, customer service will inform you of the repair charge that will be billed against your purchase order.

The RGA number should be prominently displayed on the box and included on all paperwork enclosed with the repair. Merchandise returned for repair should be safely packaged in protective wrapping. Customer must include, with the returned merchandise the 1) purchase order number; 2) correct shipping and billing address; and 3) Repair Order Form or statement of problem/reason for repair.

All repairs should be shipped freight prepaid to:

Medtronic Xomed, Inc.  
4102 Southpoint Boulevard  
Jacksonville, FL 32216  
ATTN: RGA # \_\_\_\_\_

## Important

Due to the potential threat to health, all instruments and devices returned for repair must have been decontaminated and/or sterilized. Any product received in an unsanitary condition will be deemed unrepairable and discarded unless Medtronic ENT is able to sterilize or decontaminate the product. An additional service fee of \$250.00 will be applied to the cost of the repair of any product requiring sterilization or decontamination by Medtronic ENT.

## Loaner policy for repairs

Loaners for merchandise that is being repaired under warranty are provided at no charge and must be returned to Medtronic ENT within ten (10) days after customer's receipt of repaired merchandise. If the loaner is not returned within the allotted time frame, Medtronic ENT reserves the right to charge a daily rental fee. Rental fees will continue to be charged until the loaner is returned to Medtronic ENT.

Customer is responsible for any damage to loaner equipment while in customer's possession or during shipping back to Medtronic ENT. Customer shipping charges associated with the loaner's return.

## Repair warranties

All repairs are warranted to be free of defects in material and workmanship for ninety (90) days after the date of shipment.

## Changes to terms

**Medtronic ENT retains the right to change any of these terms or conditions without notice.**

<sup>†</sup>Terms and conditions not applicable to the purchase of SINUVA™ (mometasone furoate) sinus implants. For terms and conditions related to that product, please contact your Medtronic representative.

# Table of contents

|                                                                         |           |
|-------------------------------------------------------------------------|-----------|
| <b>Catalog directory.....</b>                                           | <b>2</b>  |
| Measurements, standards, and conversions.....                           | 7         |
| <b>Office-based solutions.....</b>                                      | <b>8</b>  |
| XPS Nexus™ ENT power console.....                                       | 9         |
| StraightShot™ M4/M5 microdebrider handpieces.....                       | 9         |
| StealthStation FlexENT™ navigation system.....                          | 10        |
| NuVent™ EM sinus dilation system.....                                   | 10        |
| NuVent™ Eustachian tube dilation balloon.....                           | 11        |
| PROPEL™ Mometasone Furoate Sinus Implants.....                          | 11        |
| <b>Powered surgical equipment.....</b>                                  | <b>13</b> |
| Powered ENT system and components, IPC™.....                            | 14        |
| Powered ENT system and components, XPS Nexus™.....                      | 15        |
| Powered ENT system and components, microdebridors.....                  | 16        |
| Powered ENT system and components, Hydrodebrider™.....                  | 17        |
| 30K burs for the StraightShot™ M4 and M5.....                           | 18        |
| Automated EM tracking blades for the StraightShot™ M4/M5 rotatable..... | 19        |
| Quadcut™ blades for the StraightShot™ M4 and M5.....                    | 20        |
| Straight sinus blades, m4 rotatable.....                                | 21        |
| Curved sinus blades, m4 rotatable.....                                  | 23        |
| Curved sinus blades, non-rotatable.....                                 | 24        |
| Straight sinus burs.....                                                | 25        |
| Straight sinus bur accessories.....                                     | 26        |
| Curved sinus burs, anterior skull base.....                             | 27        |
| Curved sinus burs.....                                                  | 28        |
| Tonsillectomy and adenoidectomy blades.....                             | 29        |
| Airway blades, m4 rotatable.....                                        | 30        |
| Airway blades, non-rotatable.....                                       | 32        |
| <b>Endoscopy products.....</b>                                          | <b>33</b> |
| ENT Slide-On™ EndoSheath™ system.....                                   | 34        |
| Endoscopes.....                                                         | 35        |
| Endo-Scrub™ 2 lens cleaning system.....                                 | 36        |
| <b>ENT navigation solutions.....</b>                                    | <b>37</b> |
| StealthStation™ ENT navigation system.....                              | 38        |
| StealthStation FlexENT™ navigation system.....                          | 38        |
| <b>Nerve monitoring equipment.....</b>                                  | <b>42</b> |
| NIM Vital™ nerve monitoring system.....                                 | 43        |
| EMG endotracheal tubes.....                                             | 44        |
| Probes: monopolar stimulating.....                                      | 45        |
| Probes: bipolar stimulating.....                                        | 46        |
| Other disposables.....                                                  | 46        |
| Electrodes: recording.....                                              | 47        |
| <b>Parathyroid detection equipment.....</b>                             | <b>49</b> |
| PTeye™ parathyroid detection system.....                                | 50        |

# Table of contents

|                                                                        |           |
|------------------------------------------------------------------------|-----------|
| <b>Ventilation tubes and myringotomy accessories .....</b>             | <b>51</b> |
| Ventilation tubes and myringotomy accessories .....                    | 52        |
| Ventilation tubes: temporary.....                                      | 53        |
| Ventilation tubes: permanent and specialty .....                       | 57        |
| <b>Ossicular reconstruction implants and accessories .....</b>         | <b>58</b> |
| Partial ossicular reconstruction implants: titanium.....               | 59        |
| Partial ossicular reconstruction implants: dense hydroxylapatite ..... | 60        |
| Partial ossicular reconstruction implants: Flex H/A™ .....             | 61        |
| Partial ossicular reconstruction implants: Flex H/A™ Necrosis .....    | 62        |
| Partial ossicular reconstruction implants: Polycel™ .....              | 63        |
| Total ossicular reconstruction implants: dense hydroxylapatite .....   | 64        |
| Total ossicular reconstruction implants: Flex H/A.....                 | 64        |
| Total ossicular reconstruction implants: Polycel™ .....                | 65        |
| Universal ossicular reconstruction implants: Flex H/A™ .....           | 65        |
| Middle ear specialty instruments .....                                 | 66        |
| Head supports.....                                                     | 66        |
| <b>Localized Drug Delivery (LDD) .....</b>                             | <b>67</b> |
| PROPEL™ Mometasone Furoate Sinus Implants .....                        | 68        |
| SINUVA™ (mometasone furoate) sinus implant .....                       | 69        |
| <b>Tissue health .....</b>                                             | <b>70</b> |
| Postoperative dressing.....                                            | 71        |
| Nasal packing.....                                                     | 72        |
| Sinus packing .....                                                    | 74        |
| Ear wicks and packs .....                                              | 75        |
| Epistaxis treatment products .....                                     | 76        |
| Splints .....                                                          | 77        |
| Septal buttons .....                                                   | 79        |
| Other products .....                                                   | 80        |

# Measurements, standards, and conversions

## Millimeter equivalents

| MM = INCHES    | MM = INCHES    | MM = INCHES     |
|----------------|----------------|-----------------|
| .1= .00384     | 29.0 = 1.14173 | 65.0 = 2.55905  |
| .2= .00787     | 30.0 = 1.18110 | 66.0 = 2.59842  |
| .3= .01181     | 31.0 = 1.22047 | 67.0 = 2.63779  |
| .4= .01575     | 32.0 = 1.25984 | 68.0 = 2.67716  |
| .5= .01968     | 33.0 = 1.29921 | 69.0 = 2.71653  |
| .6= .02362     | 34.0 = 1.33858 | 70.0 = 2.75590  |
| .7= .02756     | 35.0 = 1.37795 | 71.0 = 2.79527  |
| .8= .03149     | 36.0 = 1.41732 | 72.0 = 2.83464  |
| .9= .03543     | 37.0 = 1.45669 | 73.0 = 2.87401  |
| 1.0 = .03937   | 38.0 = 1.49606 | 74.0 = 2.91338  |
| 2.0 = .07874   | 38.1 = 1.5     | 75.0 = 2.95275  |
| 3.0 = .11811   | 39.0 = 1.53543 | 76.0 = 2.99212  |
| 4.0 = .15748   | 40.0 = 1.57480 | 76.2 = 3.0      |
| 5.0 = .19685   | 41.0 = 1.61417 | 77.0 = 3.03149  |
| 6.0 = .23622   | 42.0 = 1.65354 | 78.0 = 3.07056  |
| 7.0 = .27559   | 43.0 = 1.69291 | 79.0 = 3.11023  |
| 8.0 = .31496   | 44.0 = 1.73228 | 80.0 = 3.14960  |
| 9.0 = .35433   | 45.0 = 1.77165 | 81.0 = 3.18897  |
| 10.0 = .39370  | 46.0 = 1.81102 | 82.0 = 3.22834  |
| 11.0 = .43307  | 47.0 = 1.85039 | 83.0 = 3.26771  |
| 12.0 = .47244  | 48.0 = 1.88976 | 84.0 = 3.30708  |
| 12.7 = .5      | 49.0 = 1.92913 | 85.0 = 3.34645  |
| 13.0 = .51181  | 50.0 = 1.93500 | 86.0 = 3.38582  |
| 14.0 = .55118  | 50.8 = 2.0     | 87.0 = 3.42519  |
| 15.0 = .59055  | 51.0 = 2.00787 | 88.0 = 3.46456  |
| 16.0 = .62992  | 52.0 = 2.04724 | 88.9 = 3.5      |
| 17.0 = .66929  | 53.0 = 2.08661 | 89.0 = 3.50393  |
| 18.0 = .70866  | 54.0 = 2.12598 | 90.0 = 3.54330  |
| 19.0 = .74803  | 55.0 = 2.16535 | 91.0 = 3.58267  |
| 20.0 = .78740  | 56.0 = 2.20472 | 92.0 = 3.62204  |
| 21.0 = .82677  | 57.0 = 2.24409 | 93.0 = 3.66141  |
| 22.0 = .86614  | 58.0 = 2.83460 | 94.0 = 3.70078  |
| 23.0 = .90551  | 59.0 = 2.32283 | 95.0 = 3.74015  |
| 24.0 = .94488  | 67.0 = 2.36220 | 96.0 = 3.77952  |
| 25.0 = .98425  | 61.0 = 2.40157 | 97.0 = 3.81889  |
| 25.4 = 1.0     | 62.0 = 2.44094 | 98.0 = 3.85826  |
| 26.0 = 1.02362 | 63.0 = 2.48031 | 99.0 = 3.89763  |
| 27.0 = 1.06299 | 63.5 = 2.5     | 100.0 = 3.93700 |
| 28.0 = 1.10236 | 64.0 = 2.51968 | 100.6 = 4.0     |

## Needle tubing equivalents

| SIZE<br>GUAGE NO. | APPROX. O.D. |        |
|-------------------|--------------|--------|
|                   | MM           | INCHES |
| 26                | .4           | .018   |
| 24                | .5           | .022   |
| 22                | .7           | .028   |
| 20                | .9           | .035   |
| 18                | 1.3          | .053   |
| 16                | 1.6          | .065   |
| 14                | 2.1          | .083   |
| 12                | 2.8          | .109   |
| 10                | 3.4          | .134   |
| 8                 | 4.2          | .165   |
| 6                 | 5.1          | .203   |

## Suction tubing equivalents

| SIZE<br>FRENCH GUAGE NO. | APPROX. O.D. |        |
|--------------------------|--------------|--------|
|                          | MM           | INCHES |
| 3                        | 1.0          | .039   |
| 4                        | 1.3          | .053   |
| 5                        | 1.7          | .066   |
| 6                        | 2.0          | .079   |
| 7                        | 2.3          | .092   |
| 8                        | 2.7          | .105   |
| 9                        | 3.0          | .118   |
| 10                       | 3.3          | .131   |
| 11                       | 3.7          | .144   |
| 12                       | 4.0          | .158   |
| 13                       | 4.3          | .170   |
| 14                       | 4.7          | .184   |
| 15                       | 5.0          | .197   |
| 16                       | 5.3          | .210   |

# Medtronic

Engineering the extraordinary

ENT product catalog

# Office-based solutions



## Office-based solutions

### XPS Nexus™ ENT power console

- Connects to the StraightShot™ M4 microdebrider
- Handpiece recognition
- Portable
- Built-in irrigation pump
- Foot pedal with non-slip design (sold separately)
- IV pole (sold separately)

| Product | Description       | Qty |
|---------|-------------------|-----|
| 1894510 | XPS Nexus console | 1   |



1894510

### Integrated Power Console (IPC™) system

Our IPC™ system is a powered ENT surgery console with a range of application-specific products.

| Product | Description                                       | Qty |
|---------|---------------------------------------------------|-----|
| 1898001 | Integrated power console                          | 1   |
| 1898430 | Integrated power console multifunction foot pedal | 1   |
| 1897510 | Basket                                            | 1   |
| 1898600 | IPC cart                                          | 1   |



1898001



1898430



1898600

### IPC™ system power cords

| Product | Description                              | Qty |
|---------|------------------------------------------|-----|
| 1895820 | Power cord, standard, North America, 3 M | 1   |
| 1897821 | Power cord, standard, North America, 6 M | 1   |

### StraightShot™ M4 microdebrider handpiece

- Compatible with the IPC™ and the XPS Nexus™
- Rotatable blade options

| Product  | Description                                             | Qty |
|----------|---------------------------------------------------------|-----|
| 1898200T | StraightShot M4 high-speed microdebrider with SureTrack | 1   |



1898200T

### StraightShot™ M5 microdebrider

The StraightShot™ M5 works with:

- Factory-calibrated EM blades for Fusion™ ENT navigation
- Quadcut™ blades
- Cross-compatible with M4 blades and burs

| Product | Description                   | Qty |
|---------|-------------------------------|-----|
| 1899200 | StraightShot M5 microdebrider | 1   |



1899200

# Office-based solutions

## StealthStation FlexENT™ navigation system

### Cart configurations

| Ordering code | FLEXCARTSIDE STAND | FLEXCARTSIDE CLAMP |
|---------------|--------------------|--------------------|
| Side emitter  | ✓                  | ✓                  |
| Emitter stand | ✓                  |                    |
| Bedrail clamp |                    | ✓                  |

### Desktop configurations

| Ordering code | FLEXDESKSIDE STAND | FLEXDESKSIDE CLAMP |
|---------------|--------------------|--------------------|
| Side emitter  | ✓                  | ✓                  |
| Emitter stand | ✓                  |                    |
| Bedrail clamp |                    | ✓                  |

9736226 - StealthStation FlexENT™ software license  
 9736227 - Planning, annotations, and virtual endoscopy  
 9736155 - StealthMerge™ software  
 9735742 - Stealth™ 3D software  
 9733533XOM - Instrument tracker, ENT  
 9734887XOM - AxiEM™ non-invasive patient tracker

## NuVent™ EM sinus dilation system with integrated navigation

- For use with Medtronic ENT electromagnetic (EM) navigation
- Frontal, maxillary, and sphenoid designs
- 5 mm, 6 mm, and 7 mm balloon size options

| Product                   | Description                                           | Qty |
|---------------------------|-------------------------------------------------------|-----|
| <b>70 degree balloons</b> |                                                       |     |
| 1830517FRT70              | Balloon Seeker, EM 70 Degree Balloon 5.0 mm x 17.0 mm | 1   |
| 1830617FRT70              | Balloon Seeker, EM 70 Degree Balloon 6.0 mm x 17.0 mm | 1   |
| 1830717FRT70              | Balloon Seeker, EM 70 Degree Balloon 7.0 mm x 17.0 mm | 1   |
| <b>Frontal balloons</b>   |                                                       |     |
| 1830517FRT                | Balloon seeker, EM frontal balloon 5.0 mm x 17.0 mm   | 1   |
| 1830617FRT                | Balloon seeker, EM frontal balloon 6.0 mm x 17.0 mm   | 1   |
| 1830717FRT                | Balloon seeker, EM frontal balloon 7.0 mm x 17.0 mm   | 1   |
| <b>Maxillary balloons</b> |                                                       |     |
| 1830507MAX                | Balloon seeker, EM maxillary balloon 5.0 mm x 7.0 mm  | 1   |
| 1830607MAX                | Balloon seeker, EM maxillary balloon 6.0 mm x 7.0 mm  | 1   |
| 1830707MAX                | Balloon seeker, EM maxillary balloon 7.0 mm x 7.0 mm  | 1   |
| <b>Sphenoid balloons</b>  |                                                       |     |
| 1830517SPH                | Balloon seeker, EM sphenoid balloon 5.0 mm x 17.0 mm  | 1   |
| 1830617SPH                | Balloon seeker, EM sphenoid balloon 6.0 mm x 17.0 mm  | 1   |
| 1830717SPH                | Balloon seeker, EM sphenoid balloon 7.0 mm x 17.0 mm  | 1   |
| <b>Inflator</b>           |                                                       |     |
| 18INFKIT                  | Kit, inflator                                         | 1   |



NuVent™ balloon package



Frontal balloon (70°)



Frontal balloon (98°)



Maxillary balloon



Sphenoid balloon



Inflator Kit

## Office-based solutions

### NuVent™ Eustachian tube dilation balloon

| Product    | Description                             | Qty |
|------------|-----------------------------------------|-----|
| 1830616EST | Eustachian tube Balloon 6.0mm x 16.0 mm | 1   |

#### Inflator

|          |               |   |
|----------|---------------|---|
| 18INFKIT | Kit, inflator | 1 |
|----------|---------------|---|



NuVent™ Eustachian tube dilation balloon




Included stopcock valve

### Novapak™ nasal sinus packing and stent

- A single-use, dissolvable, nasal packing and stent for use following sinus surgery
- Dimensions 8 cm x 1.8 cm x 1.3 cm



| Product   | Description                      | Qty |
|-----------|----------------------------------|-----|
| CS3600-10 | Novapak™ nasal packing, 10PK STD | 10  |
| CS3900-10 | Novapak™ nasal packing, 10PK FRM | 10  |
| CS3600-2  | Novapak™ nasal packing, 2PK STD  | 2   |
| CS3900-2  | Novapak™ nasal packing, 2PK FRM, | 2   |



### Localized Drug Delivery (LDD)

#### PROPEL™ Mometasone Furoate Sinus Implants


| Product | Description                                                   | Qty |
|---------|---------------------------------------------------------------|-----|
| 70011   | PROPEL sinus implant                                          | 1   |
| 700115  | PROPEL sinus implant                                          | 5   |
| 60011   | PROPEL Mini sinus implant                                     | 1   |
| 600115  | PROPEL Mini sinus implant                                     | 5   |
| 60044-  | PROPEL Mini sinus implant with Straight Delivery System (SDS) | 1   |
| 60044-5 | PROPEL Mini sinus implant with Straight Delivery System (SDS) | 5   |
| 50011-  | PROPEL Contour sinus implant                                  | 1   |
| 50011-5 | PROPEL Contour sinus implant                                  | 5   |



For more information refer to <https://propelopens.com/>

### SINUVA™ (mometasone furoate) sinus implant

| NDC number  | Description                                         | Qty |
|-------------|-----------------------------------------------------|-----|
| 10599000301 | SINUVA™ (mometasone furoate) sinus implant 1350 mcg | 1   |



# Office-based solutions

## SINUVA™ (mometasone furoate) sinus implant

### Ordering

SINUVA sinus implants can be obtained through two accessible avenues:

1. Specialty distributor: providers can purchase SINUVA sinus implants directly from a distributor. Benefits verification and prior authorization from the payer is recommended prior to implanting a patient.
2. Specialty pharmacies: In some cases, SINUVA sinus implants can be filled under the pharmacy benefit by a select specialty pharmacy and delivered directly to the provider. Verifying the patient's benefits will establish the path to fulfillment.

### INDICATION

SINUVA Sinus Implant is a corticosteroid-eluting implant indicated for the treatment of chronic rhinosinusitis with nasal polyps in adult patients  $\geq 18$  years of age who have had ethmoid sinus surgery.

### IMPORTANT SAFETY INFORMATION

#### CONTRAINDICATIONS

Patients with known hypersensitivity to mometasone furoate and any of the ingredients of the SINUVA Sinus Implant.

#### WARNINGS AND PRECAUTIONS

**Local Nasal Adverse Reactions:** Monitor nasal mucosa adjacent to the SINUVA Sinus Implant for any signs of bleeding (epistaxis), irritation, infection, or perforation. Avoid use in patients with nasal ulcers or trauma.

**Glaucoma and Cataracts:** Nasal steroids may result in development of glaucoma and/or cataracts. Glaucoma, cataracts, and clinically significant elevation of intraocular pressure were not observed in patients from the treatment group of one randomized controlled clinical study (N = 53) who underwent bilateral placement of SINUVA Sinus Implants. Close monitoring is warranted in patients with a change in vision or with a history of increased intraocular pressure, glaucoma, and/or cataracts.

**Hypersensitivity Reactions:** Hypersensitivity reactions, including rash, pruritus, angioedema have been reported with the use of corticosteroids.

**Immunosuppression and Risk of Infections:** Persons who are using drugs that suppress the immune system, such as corticosteroids, including SINUVA Sinus Implant are more susceptible to infections than healthy individuals. The safety and effectiveness of SINUVA Sinus Implant have not been established in pediatric patients less than 18 years of age and SINUVA is not indicated for use in this population. Corticosteroids should be used with caution, if at all, in patients with active or quiescent tuberculosis infection of the respiratory tract; untreated systemic fungal, bacterial, viral, or parasitic infections; or ocular herpes simplex.

**Hypercorticism and Adrenal Suppression:** If corticosteroid effects such as hypercorticism and adrenal suppression appear in patients, consider sinus implant removal.

#### ADVERSE REACTIONS

The most common adverse reactions observed ( $> 1\%$  of subjects) in clinical studies were asthma, headache, epistaxis, presyncope, bronchitis, otitis media, and nasopharyngitis.

#### POSTMARKETING EXPERIENCE

The following adverse reactions have been identified during post-approval use of the SINUVA sinus implant. These events include implant migration, lack of efficacy, nasal pain, headache, epistaxis.

**RX only. Please see Full Prescribing Information for SINUVA at [SINUVA.com/hcp](https://www.medtronic.com/hcp)**

# Medtronic

Engineering the extraordinary

ENT product catalog

# Powered surgical equipment

# Powered ENT system and components, IPC™

## Integrated Power Console (IPC™) system

Our IPC™ system is a powered ENT surgery console with a range of application-specific products.

| Product | Description                                       | Qty |
|---------|---------------------------------------------------|-----|
| 1898001 | Integrated power console                          | 1   |
| 1898430 | Integrated power console multifunction foot pedal | 1   |
| 1897510 | Basket                                            | 1   |
| 1898600 | IPC cart                                          | 1   |



1898001

## IPC™ system power cords

| Product | Description                              | Qty |
|---------|------------------------------------------|-----|
| 1895820 | Power cord, standard, North America, 3 M | 1   |
| 1897821 | Power cord, standard, North America, 6 M | 1   |



1898430

## IPC™ system manuals

| Product  | Description                     | Qty |
|----------|---------------------------------|-----|
| 1898851  | IPC manual (English only)       | 1   |
| 1898851A | IPC manual (EL, EN, ES, FR, PT) | 1   |
| 1898851B | IPC manual (DE, EN, FR, IT, NL) | 1   |
| 1898851C | IPC manual (DA, EN, FI, NO, SV) | 1   |
| 1898851D | IPC manual (CS, EN, HU, PL, TR) | 1   |
| 1898851E | IPC manual (EN, RU) – CD only   | 1   |

## IPC™ suction irrigation set with IntelliFlow™ remote control

System includes:

- Coextruded suction-irrigation tubing set
- Multi-size irrigation tip adapters (4) IntelliFlow™ remote control

| Product  | Description                   | Qty |
|----------|-------------------------------|-----|
| 3318604E | IPC suction-irrigation system | 1   |



1898600



3318604E

# Powered ENT system and components, XPS Nexus™

## XPS Nexus™ ENT power console

- Connects to the StraightShot™ M4 microdebrider
- Handpiece recognition
- Portable
- Built-in irrigation pump
- Foot pedal with non-slip design (sold separately)
- IV pole (sold separately)



| Product | Description                        | Qty |
|---------|------------------------------------|-----|
| 1894510 | XPS Nexus™ console                 | 1   |
| 1894520 | XPS Nexus™ foot control            | 1   |
| 1894515 | XPS Nexus™ IV pole                 | 1   |
| 1897821 | Power cord, 6 meter, IEC 320, 115V | 1   |




# Powered ENT system and components, microdebridors

## StraightShot™ M5 Microdebrider

The StraightShot™ M5 works with:

- 30K high-speed burs
- Factory-calibrated EM blades for ENT navigation
- Quadcut™ blades
- Cross-compatible with M4 blades and burs

| Product | Description                      | Qty |
|---------|----------------------------------|-----|
| 1899200 | StraightShot™ M5 microdebrider   | 1   |
| 1899076 | StraightShot™ M5 instrument tray | 1   |



## StraightShot™ M4 microdebrider handpiece

Compatible with the IPC™ and the XPS™ Nexus

- Rotating blade and handpiece

| Product  | Description                                              | Qty |
|----------|----------------------------------------------------------|-----|
| 1898200T | StraightShot™ M4 high-speed microdebrider with SureTrack | 1   |
| 1898075  | StraightShot™ M4 instrument tray                         | 1   |
| 1895522  | StraightShot™ M4 tubing                                  | 5   |





# Powered ENT system and components, microdebridors

## Hydrodebrider™ system

### Power and access

- Powered irrigation system for the paranasal sinuses
- With a rotating spray, the flex-tip incorporates suction and 270° articulation

| Product | Description                       | Qty |
|---------|-----------------------------------|-----|
| 1913001 | Hydrodebrider™ console            | 1   |
| 1914001 | Hydrodebrider™ standard handpiece | 1   |
| 1921002 | Hydrodebrider™ frontal handpiece  | 1   |




1913001

## Hydrodebrider – powered ENT power cords

| Product | Description                     | Qty |
|---------|---------------------------------|-----|
| 1895820 | Power cord, standard, 3 M       | 1   |
| 1897821 | Power cord, standard, 115V, 6 M | 1   |



1914001



Standard handpiece tip close-up



1921002



Frontal handpiece tip close-up

# 30K burs for the StraightShot™ M5

## 30K burs

| Product    | Description                                          | Length | Qty |
|------------|------------------------------------------------------|--------|-----|
| 1884015RTD | 30K tapered diamond choanal atresia bur, 4.0 mm, 15° | 13 cm  | 1   |
| 1883040BLD | 30K bullet diamond bur, 3.0 mm, 40°                  | 13 cm  | 1   |
| 1883040BLC | 30K bullet cutting bur, 3.0 mm, 40°                  | 13 cm  | 1   |
| 1884040RTC | 30K reverse tapered cutting bur, 4.0 mm, 40°         | 13 cm  | 1   |
| 1883040BRC | 30K frontal finesse barrel cutting bur, 3.0 mm, 40°  | 13 cm  | 1   |
| 1883655BRC | 30K barrel cutting bur, 3.6 mm, 55°                  | 13 cm  | 1   |
| 1884070RTD | 30K reverse tapered diamond bur, 4.0 mm, 70°         | 13 cm  | 1   |
| 1884070RTC | 30K reverse tapered cutting bur, 4.0 mm, 70°         | 13 cm  | 1   |
| 1883070BLD | 30K bullet diamond bur, 3.0 mm, 70°                  | 13 cm  | 1   |
| 1883070BLC | 30K bullet cutting bur, 3.0 mm, 70°                  | 13 cm  | 1   |



Maximum speed for 30K burs is 12,000 rpm when used with StraightShot™ M4 microdebrider handpiece

# Automated EM tracking blades for the StraightShot™ M4/M5 rotatable

## QUADCUT™ BLADES



### 4.3 mm Quadcut™ straight rotatable blade with automated EM tracking 1884380EM

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate



### 3.4 mm Quadcut™ straight rotatable blade with automated EM tracking 1883480EM

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate



### 3.0 mm Quadcut™ straight rotatable blade with automated EM tracking 1883080EM

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate

## TRICUT™ BLADE



### 4.0 mm Tricut™ straight rotatable blade with automated EM tracking 1884080EM


- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate

## RAD™ BLADES



### 4.0 mm RAD™ 12 curved rotatable blade with automated EM tracking 1884012EM

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate



### 4.0 mm RAD™ 40 curved rotatable blade with automated EM tracking 1884006EM

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate

## IRRIGATION TUBING

### Irrigation tubing for blades and burs 1895522

- For use with IPC™ blades and burs
- 5 each

Speeds are suggested rpm (revolutions per minute), operated in oscillation mode for blades and (forward) mode for burs.

Measurements are listed in millimeters unless otherwise specified.

# Quadcut™ blades for the StraightShot™ M4 and M5

## QUADCUT™ BLADES

360°



### **4.3 mm Quadcut™ blade** **1884380HR**

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 5 each with irrigation tubing

## QUADCUT™ BLADES

360°




### **3.4 mm Quadcut™ blade** **1883480HRE**

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing

## QUADCUT™ BLADES

360°



### **3.0 mm Quadcut™ blade** **1883080HRE**

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing

*Speeds are suggested rpm (revolutions per minute), operated in oscillation mode for blades and (forward) mode for burs.  
Measurements are listed in millimeters unless otherwise specified.*

# Straight sinus blades, M4 and M5 rotatable

## TRICUT™ BLADES

360°



### **4.0 mm Tricut™ blade** **1884004HR**

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 5 each with irrigation tubing

### **4.0 mm Tricut™ blade** **1884080HR**

- 13.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 5 each with irrigation tubing

## TRICUT™ BLADES

360°



### **3.5 mm Tricut™ blade** **1883504HR**

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 5 each with irrigation tubing



### **2.9 mm Tricut™ blade** **StraightShot™ M5 rotates blade 360°**

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing

## SERRATED BLADES

360°



### **4.0 mm Serrated blade** **1884002HRE**

- 11.0 cm with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing



### **3.5 mm serrated blade** **1883502HRE**

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing



### **2.9 mm serrated blade** **1882902HRE**

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing

# Straight sinus blades, M4 and M5 rotatable

## SILVER BULLET™

360°



### 4.0 mm Silver Bullet™ blade 1884005HRE

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing
- Developed in conjunction with Rodney Lusk, M.D.



### 2.9 mm Silver Bullet™ blade 1882905HRE

- 11.0 cm long with straight shaft
- StraightShot™ M5 rotates blade 360°
- Operating speed: up to 7,500 rpm, oscillate
- 1 each with irrigation tubing
- Developed in conjunction with Rodney Lusk, M.D.

## TURBINATE

360°



### 2.9 mm inferior turbinate blade 1882940HR

- 11.0 cm long
- StraightShot™ M5 rotates blade 360°
- Straight shaft with elevator
- Operating speed: up to 3,000 rpm, oscillate
- 5 each with irrigation tubing
- Developed in conjunction with Laurence O'Halloran, M.D.



### 2.0 mm inferior turbinate blade 1882040HR

- 11.0 cm long
- StraightShot™ M5 rotates blade 360°
- Straight shaft with elevator
- Operating speed: up to 3,000 rpm, oscillate
- 5 each with irrigation tubing
- Developed in conjunction with Laurence O'Halloran, M.D.

# Curved sinus blades, M4 and M5 rotatable

## RAD™ 12 BLADES

360°



### 4 mm RAD™ 12 blade 1884012HR

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 7,500 rpm, oscillate
- 5 each, irrigation tubing separate



### 3.5 mm RAD™ 12 blade 1883512HRE

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate


## IRRIGATION TUBING

### Irrigation tubing for blades and burs 1895522

- For use with XPS™ blades and burs
- 5 each

## RAD™ 40 BLADES

360°



### 4 mm RAD™ 40 blade 1884006HR

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 7,500 rpm, oscillate
- 5 each, irrigation tubing separate




### 3.5 mm RAD™ 40 blade 1883506HRE

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 7,500 rpm, oscillate
- 1 each, irrigation tubing separate

## RAD™ 90 BLADE

360°




### 3.5 mm RAD™ 90 blade 1883519HR

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 3,000 rpm, oscillate
- 3 each, irrigation tubing separate


## RAD™ 60 BLADES

360°



### 4 mm RAD™ 60 blade 1884016HR

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 5,000 rpm, oscillate
- 5 each, irrigation tubing separate




### 3.5 mm RAD™ 60 blade 1883516HRE

- 11.0 cm long with curved shaft
- StraightShot™ M5 rotates blade tip 360° without shaft rotation
- Operating speed: up to 5,000 rpm, oscillate
- 1 each, irrigation tubing separate

## Curved sinus blades, non-rotatable



### RAD™ 60 BLADE



#### **2.9 mm RAD™ 60 blade 1882916**

- 11.0 cm long with curved shaft
- Operating speed: up to 1,500 rpm, oscillate
- 3 each, irrigation tubing separate

### RAD™ 120 BLADE



#### **3.5 mm RAD™ 120 blade 1883517**

- 11.0 cm long with curved shaft
- Operating speed: up to 7,500 rpm, oscillate
- 3 each, irrigation tubing separate



# Straight sinus burs

## OVAL BUR



### **3.2 mm oval bur, high-speed 1883264HS**

- 12.5 cm long with straight shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate

## ROUND BURS



### **4.5 mm round bur, high-speed 1884560HS**

- 12.5 cm long with straight shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate



### **3.2 mm round bur, high-speed 1883262HS**

- 12.5 cm long with straight shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate



### **2.9 mm pediatric round bur 1882960**

- 10.0 cm long with straight shaft
- Operating speed: up to 5,000 rpm (forward)
- 5 each, irrigation tubing separate

## ROUTER BUR




### **4.5 mm aggressive router bur, high-speed 1884562HS**

- 12.5 cm long with straight shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate

# Curved sinus burs, anterior skull base


## ASB CUTTING BUR



**4 mm anterior skull base cutting bur,\* 15°**  
**1884075HSE**


- 15.0 cm long
- Operating speed: up to 12,000 rpm (forward)
- 1 each, irrigation tubing separate
- Developed in conjunction with PJ Wormald, M.D., and Aldo Stamm, M.D.

## ASB DIAMOND BURS



**5 mm anterior skull base diamond bur,\* 15°**  
**1885076HSE**

- 15.0 cm long
- Operating speed: up to 12,000 rpm (forward)
- 1 each, irrigation tubing separate
- Developed in conjunction with PJ Wormald, M.D., and Aldo Stamm, M.D.



**3.2 mm anterior skull base diamond bur,\* 15°**  
**1883274HSE**


- 15.0 cm long
- Operating speed: up to 12,000 rpm (forward)
- 1 each, irrigation tubing separate
- Developed in conjunction with PJ Wormald, M.D., and Aldo Stamm, M.D.

## ASB DIAMOND BURS



**3.2 mm anterior skull base diamond bur,\* 40°**  
**1883277HSE**

- 15.0 cm long
- Operating speed: up to 12,000 rpm (forward)
- 1 each, irrigation tubing separate
- Developed in conjunction with PJ Wormald, M.D., and Aldo Stamm, M.D.




**5 mm anterior skull base diamond bur,\* 70°**  
**1885078HSE**

- 13.0 cm long
- Operating speed: up to 12,000 rpm (forward)
- 1 each, irrigation tubing separate
- Developed in conjunction with PJ Wormald, M.D., and Aldo Stamm, M.D.

# Curved sinus burs

## ROUND DIAMOND BUR



### **5 mm curved round diamond bur, high-speed** **1885061HS**

- 12.5 cm long with curved shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate
- Developed in conjunction with David Kennedy, M.D.

## DCR BURS



### **4 mm curved DCR bur, high-speed** **1884068HS**

- 11.0 cm long with curved shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate
- Developed in conjunction with Michael Mercandetti, M.D.



### **2.5 mm curved diamond DCR bur, high-speed** **1882569HS**

- 11.0 cm long with curved shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate
- Developed in conjunction with PJ Wormald, M.D.

## SEPTOPLASTY BUR




### **3.2 mm septoplasty bur, high-speed** **1883212HS**

- 11.0 cm long with curved shaft
- Operating speed: up to 12,000 rpm (forward)
- 3 each, irrigation tubing separate
- Developed in conjunction with Donald Leopold, M.D., and Eileen Raynor, M.D.

# Tonsillectomy and adenoidectomy blades

## RADENOID™ BLADES



### 4.5 mm RADenoid™ adult blade 1884507

- 13.0 cm long with curved 45° blade
- Operating speed: 1,500 rpm
- 5 each, irrigation tubing separate
- Designed in conjunction with Max April, M.D., and J. Lindhe Guarisco, M.D.



### 4 mm RADenoid™ blade 1884008

- 11.0 cm long with curved 40° blade
- Operating speed: 1,500 rpm
- 5 each, irrigation tubing separate
- Designed in conjunction with Max April, M.D., and J. Lindhe Guarisco, M.D.

## TONSILLECTOMY BLADE



### 4 mm tonsillectomy blade 1884013

- 11.0 cm
- 12° blade
- Operating speed: 1,500 rpm
- 5 each, irrigation tubing separate

## T&A BLADE SET



### Powered T&A blade set 1884008TA

- 13.0 cm
- Removable inner cutting tube
- Operating speed: 1,500 rpm
- 5 each, irrigation tubing separate
- Developed in conjunction with Peter J. Koltai, M.D.

# Airway blades, M4 and M5 rotatable

## SKIMMER™ BLADES



### 2.9 mm Skimmer™ angle-tip blade

**1882979HRE**

- 13.0 cm long double-curved blade
- Operating speed: 60-500 rpm
- 15°
- 1 each with irrigation tubing

### 2.9 mm Skimmer™ angle-tip blade

**1882925HRE**

- 18.0 cm long double-curved blade
- Operating speed: 60-500 rpm
- 15°
- 1 each with irrigation tubing

## SKIMMER™ BLADES



### 2.9 mm Skimmer™ angle-tip blade

**1882923HRE**


- 22.0 cm long double-curved blade
- Operating speed: 60-500 rpm
- 15°
- 1 each with irrigation tubing

### 2.9 mm Skimmer™ angle-tip blade

**1882924HRE**

- 27.0 cm long double-curved blade
- Operating speed: 60-500 rpm
- 15°
- 1 each with irrigation tubing

## TRICUT™ BLADES



### 4 mm Tricut™ angle-tip laryngeal blade

**1884030HRE**

- 22.0 cm long double-curved blade
- Operating speed: 500-1,200 rpm
- 1 each with irrigation tubing
- Developed in conjunction with William Lunn, M.D., and Armin Ernst, M.D.

Speeds are suggested rpm (revolutions per minute), operated in oscillation mode for blades and (forward) mode for burs.

Measurements are listed in millimeters unless otherwise specified.

# Airway blades, rotatable

## TRICUT™ BLADES



**4 mm Tricut™ angle-tip subglottic blade**  
**1884031HRE**

- 27.0 cm long double-curved blade
- Operating speed: 500-1,200 rpm
- 1 each with irrigation tubing
- Developed in conjunction with William Lunn, M.D., and Armin Ernst, M.D.

## TRICUT™ BLADES



**4 mm Tricut™ angle-tip tracheal blade**  
**1884033HRE**

- 37.0 cm long double-curved blade
- Operating speed: 500-1,200 rpm
- 1 each with irrigation tubing
- Developed in conjunction with William Lunn, M.D., and Armin Ernst, M.D.

## TRICUT™ BLADES



**4 mm Tricut™ angle-tip bronchial blade**  
**1884035HRE**

- 45.0 cm long double-curved blade
- Operating speed: 500-1,200 rpm
- 1 each with irrigation tubing
- Developed in conjunction with William Lunn, M.D., and Armin Ernst, M.D.

# Airway blades, non-rotatable

## SKIMMER™ BLADES



### **2.9 mm Skimmer™ angle-tip blade** **1882925**

- 18.0 cm long double-curved blade
- Operating speed: 60-500 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Craig Derkay, M.D., and David Darrow, M.D.

### **2.9 mm Skimmer™ angle-tip blade** **1882923**

- 22.5 cm long double-curved blade
- Operating speed: 60-500 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Craig Derkay, M.D., and David Darrow, M.D.


## SKIMMER™ BLADES



### **3.5 mm Skimmer™ angle-tip subglottic blade** **1883524**

- 27.5 cm long double-curved blade
- Operating speed: 60-500 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Charles Myer, III, M.D.; Paul Wilging, M.D.; Brian Wiatrak, M.D.; Paul Flint, M.D.; David Parsons, M.D.; and John Little, M.D.

## SKIMMER™ BLADES



### **4 mm Skimmer™ angle-tip laryngeal blade** **1884023**

- 22.5 cm long double-curved blade
- Operating speed: 60-500 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Charles Myer, III, M.D.; Paul Wilging, M.D.; Brian Wiatrak, M.D.; Paul Flint, M.D.; David Parsons, M.D.; and John Little, M.D.

### **4 mm Skimmer™ angle-tip subglottic blade** **1884024**

- 27.5 cm long double-curved blade
- Operating speed: 60-500 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Charles Myer, III, M.D.; Paul Wilging, M.D.; Brian Wiatrak, M.D.; Paul Flint, M.D.; David Parsons, M.D.; and John Little, M.D.

# Airway blades, non-rotatable

## TRICUT™ BLADES



**4 mm Tricut™ angle-tip laryngeal blade**  
**1884030**

- 22.5 cm long double-curved blade
- Operating speed: 1,500 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Paul Flint, M.D., and John Little, M.D.

**4 mm Tricut™ angle-tip subglottic blade**  
**1884031**

- 27.5 cm long double-curved blade
- Operating speed: 1,500 rpm
- 3 each with irrigation tubing



**4 mm Tricut™ straight-tip laryngeal blade**  
**1884020**

- 22.5 cm long
- Operating speed: 1,200 rpm
- 3 each with irrigation tubing
- Developed in conjunction with Paul Flint, M.D., and John Little, M.D.

## SERRATED BLADES



**4 mm serrated angle-tip tracheal blade**  
**1884033**

- 37.0 cm long
- Operating speed: 1,200 rpm
- 1 each with irrigation tubing
- Developed in conjunction with Paul Flint, M.D.



**2.9 mm serrated angle-tip blade**  
**1882936E**

- 18.0 cm long double-curved blade
- Operating speed: 500-1,500 rpm
- 1 each with irrigation tubing

**2.9 mm serrated angle-tip blade**  
**1882937E**

- 22.0 cm long double-curved blade
- Operating speed: 500-1,500 rpm
- 1 each with irrigation tubing

## TRACHEAL BLADE



**4 mm straight tracheal blade**  
**1884032**

- 37.0 cm long
- Operating speed: 1,200 rpm
- 1 each with irrigation tubing
- Developed in conjunction with Paul Flint, M.D., and John Little, M.D.



# Medtronic

Engineering the extraordinary

ENT product catalog

# Endoscopy products

# ENT Slide-On™ EndoSheath™ system

## ENT Slide-On™ EndoSheath™ system

The EndoSheath™ system provides a sterile, disposable microbial barrier for use with flexible ENT scopes. It is indicated for endoscopic exams performed with a manufacturer-specific flexible nasopharyngoscope.



## ENT Slide-On™ EndoSheath™ system

### Slide-On™ EndoSheath™ system selection guide

| Make            | Model         | ENT sheath |
|-----------------|---------------|------------|
| Olympus         | ENF-P2        | 224301     |
| Olympus         | ENF-P3        | 224301     |
| Olympus         | ENF-XP        | 224101     |
| Olympus         | ENF-V2        | 226401     |
| Olympus         | ENF-VQ        | 225501     |
| Pentax          | FNL-10S       | 225401     |
| Pentax          | FNL-10BS      | 225401     |
| Pentax          | FNL-7RP3      | 224101     |
| Pentax          | FNL-10RP2     |            |
| Pentax          | FNL-10RP3     | 225401     |
| Pentax          | VNL-1130      | 225501     |
| Pentax          | VNL-1190STK   | 225501     |
| Pentax          | VNL-1170K     | 225501     |
| Pentax          | VNL-1070STK   | 226401     |
| Pentax          | FNL-10RBS     | 225401     |
| Machida         | ENT-3L        | 225301     |
| Machida         | ENT-4L        | 222401     |
| Machida         | ENT-30PIII    | 225301     |
| Medtronic       | 1880500       | 225401     |
| Optim           | ENTity        | 224401     |
| Optim           | ENTity SD     | 225101     |
| Optim           | ENTity XL     | 224401     |
| Optim           | ENTity SDXL   | 225101     |
| Optim           | ENT 3.6F      | 224401     |
| Karl Storz      | 11101RP-1     | 226401     |
| Karl Storz      | 11101RP2      | 226401     |
| Karl Storz      | 11101SK1      | 226101     |
| Welch Allyn     | RL-150        | 225401     |
| Gyrus ACMI      | Duraview OL-1 | 225301     |
| Gyrus ACMI      | 800-600       | 225401     |
| Vision Sciences | ENT-4500      | 225101     |
| Vision Sciences | ENT-4000      | 225601     |
| Vision Sciences | ENT-3000      | 225401     |
| Vision Sciences | ENT-2000      | 225401     |
| Vision Sciences | ENT-1000      | 225101     |
| Vision Sciences | E-F100†       | 225301     |
| Xion            | EV-N          | 225501     |
| Xion            | EF-NS         | 226401     |
| Xion            | EV-NE         | 222401     |
| Richard Wolf    | 7224.001      |            |
| Richard Wolf    | 7224.001 S02  |            |





†Requires use of catalog number 021106, E-F100 Slide-On™ Endosheath™ adapter

# Endoscopes

## Medtronic sinuscope

- Built-in wide-angle field of view of 0°, 30°, 45°, and 70° endoscopes
- Flash autoclaving produces a sterile scope and provides an alternative to soak solutions
- Color-coded for identification
- HD-compatible

| Product  | Diameter | Length   | DOV | FOV  | Color  | Qty |
|----------|----------|----------|-----|------|--------|-----|
| 31.0810a | 2.7 mm   | 187.5 mm | 0°  | 97°  | Green  | 1   |
| 31.0811a | 2.7 mm   | 187.5 mm | 30° | 87°  | Red    | 1   |
| 31.0801a | 4 mm     | 175 mm   | 0°  | 102° | Green  | 1   |
| 31.0802a | 4 mm     | 175 mm   | 30° | 100° | Red    | 1   |
| 31.0803a | 4 mm     | 175 mm   | 70° | 102° | Yellow | 1   |
| 31.0809a | 4 mm     | 175 mm   | 45° | 95°  | Black  | 1   |



# Endo-Scrub™ 2 lens cleaning system

## Endo-Scrub™ 2 lens cleaning system accessories

| Product | Description               | Qty |
|---------|---------------------------|-----|
| 1852000 | Endo-Scrub™ 2 foot switch | 1   |
| 1991015 | Endo-Scrub™ finger switch | 1   |
| 1912030 | Endo-Scrub™ 2 tubing      | 5   |

## Endo-Scrub™ 2 sheaths

| Item no.  | Size   | Angle | Compatible scopes        | Qty |
|-----------|--------|-------|--------------------------|-----|
| 1912023   | 2.7 mm | 0°    | Karl Storz, 27018A       | 5   |
| 1912002   | 4 mm   | 0°    | Karl Storz, 7200A        | 5   |
| 1912000   | 4 mm   | 0°    | Karl Storz, 7230AA       | 5   |
| 197230BVA | 4 mm   | 30°   | Karl Storz, 7230BVA      | 5   |
| 1912010   | 4 mm   | 30°   | Karl Storz, 7230BA       | 5   |
| 1912013   | 4 mm   | 45°   | Karl Storz, 7210FA       | 5   |
| 197230FLA | 4 mm   | 45°   | Karl Storz, 7230FLA      | 5   |
| 1912013   | 4 mm   | 45°   | Karl Storz, 7230FA       | 5   |
| 197230CRA | 4 mm   | 70°   | Karl Storz, 7230CRA      | 5   |
| 197230CVA | 4 mm   | 70°   | Karl Storz, 7230CVA      | 5   |
| 1912020   | 4 mm   | 70°   | Karl Storz, 7230CA       | 5   |
| 1912050   | 4 mm   | 0°    | Lintratec, HD4000R       | 5   |
| 1912051   | 4 mm   | 30°   | Lintratec, HD4300R       | 5   |
| 1912051   | 4 mm   | 30°   | Lintratec, QL4300R       | 5   |
| 1912023   | 2.7 mm | 0°    | Medtronic, 1880160       | 5   |
| 1912000   | 4 mm   | 0°    | Medtronic, 1880152       | 5   |
| 1912010   | 4 mm   | 30°   | Medtronic, 1880153       | 5   |
| 1912015   | 4 mm   | 45°   | Medtronic, 1880155       | 5   |
| 1912020   | 4 mm   | 70°   | Medtronic, 1880154       | 5   |
| 1912032   | 4 mm   | 0°    | Olympus, A7594A          | 5   |
| 1912033   | 4 mm   | 0°    | Olympus, A70940A         | 5   |
| 1912038   | 4 mm   | 0°    | Olympus, WA96200A        | 5   |
| 1912031   | 4 mm   | 30°   | Olympus, A7595A          | 5   |
| 1912034   | 4 mm   | 30°   | Olympus, A70941A         | 5   |
| 1912039   | 4 mm   | 30°   | Olympus, WA96201A        | 5   |
| 1912008   | 4 mm   | 0°    | Richard Wolf, 888.0431   | 5   |
| 1912008   | 4 mm   | 0°    | Richard Wolf, 585412     | 5   |
| 1912018   | 4 mm   | 25°   | Richard Wolf, 8475.42    | 5   |
| 1912019   | 4 mm   | 25°   | Richard Wolf, 8885.433   | 5   |
| 1912019   | 4 mm   | 25°   | Richard Wolf, 8880.433   | 5   |
| 1912028   | 2.7 mm | 30°   | Smith & Nephew, 23-0854  | 5   |
| 1912004   | 4 mm   | 0°    | Smith & Nephew, 23-0850  | 5   |
| 1912004   | 4 mm   | 0°    | Smith & Nephew, 70230910 | 5   |
| 1912014   | 4 mm   | 30°   | Smith & Nephew, 23-0851  | 5   |
| 1912040   | 4 mm   | 0°    | Stryker, 502-477-011     | 5   |
| 1912041   | 4 mm   | 30°   | Stryker, 502-477-031     | 5   |
| 1912043   | 4 mm   | 70°   | Stryker, 502-477-071     | 5   |



Note: All 0° sheaths are packaged with tubing. Tubing may be purchased separately if a 0° sheath is not being used.

# Medtronic

Engineering the extraordinary

ENT product catalog

# Navigation solutions

# ENT navigation solutions

## StealthStation™ ENT navigation system

| Product    | Description                                           | Qty |
|------------|-------------------------------------------------------|-----|
| 9735669KIT | StealthStation™ ENT navigation system with power pack | 1   |

## StealthStation FlexENT™ navigation system

### Cart configurations

| Ordering code | FLEXCARTFLAT | FLEXCARTSIDE STAND | FLEXCARTSIDECLAMP |
|---------------|--------------|--------------------|-------------------|
| Flat emitter  | ✓            |                    |                   |
| Side emitter  |              | ✓                  | ✓                 |
| Emitter stand |              | ✓                  |                   |
| Bedrail clamp |              |                    | ✓                 |

### Desktop configurations

| Ordering code | FLEXDESKFLAT | FLEXDESKSIDE STAND | FLEXDESKSIDECLAMP |
|---------------|--------------|--------------------|-------------------|
| Flat emitter  | ✓            |                    |                   |
| Side emitter  |              | ✓                  | ✓                 |
| Emitter stand |              | ✓                  |                   |
| Bedrail clamp |              |                    | ✓                 |

9736226 - StealthStation FlexENT™ software license

9736227 - Planning, annotations and virtual endoscopy

9736155 - StealthMerge™ software

9735742 - Stealth™ 3D software

NOTE: Product codes may vary if ordered standalone or in a kit. Please confirm the correct product code for your site with your Medtronic representative.

## Reusable standard instrument set


| Product    | Description                      | Qty |
|------------|----------------------------------|-----|
| 9733452XOM | Reusable standard instrument set | 1   |

### Set Includes:

|         |                               |   |
|---------|-------------------------------|---|
| 9733445 | Head frame kit                | 1 |
| 9733446 | Registration probe            | 1 |
| 9733447 | Straight probe                | 1 |
| 9733449 | Straight suction              | 1 |
| 9733450 | Curved suction, 70°           | 1 |
| 9733451 | Curved suction, 90°           | 1 |
| 9733448 | Ostium seeker                 | 1 |
| 9733455 | Instrument tray               | 1 |
| 9732379 | Imaging protocol              | 1 |
| 9730713 | Sterilization reference sheet | 1 |

## Supplemental instrument set

| Product    | Description                 | Qty |
|------------|-----------------------------|-----|
| 9734378XOM | Supplemental instrument set | 1   |



NOTE: The Medtronic family of navigation systems (StealthStation™ ENT, Fusion™, and FUSION™ Compact) will work with the instrument sets, balloons, suctions, and blades listed on these pages.

# ENT navigation solutions

Set includes:

|            |                        |    |
|------------|------------------------|----|
| 9734306    | Frontal suction, 45°   | 1  |
| 9734307    | Frontal suction, 90°   | 1  |
| 9734308    | Small straight suction | 1  |
| 9734311    | Elevator               | 1  |
| 9733446    | Registration probe     | 1  |
| 9733445    | Head frame kit         | 1  |
| 9734320    | Instrument tray        | 1  |
| 960363     | Strap, ENT FESS        | 3  |
| 9732500XOM | Adhesive pad           | 25 |
| 9732501XOM | Silicone pad           | 3  |
| 9732379    | Imaging protocol       | 1  |

## Optional accessories

| Product    | Description                                 | Qty |
|------------|---------------------------------------------|-----|
| 9735960    | Side emitter                                | 1   |
| 9734691XOM | Surgeon mouse kit                           | 1   |
| 9733533XOM | Instrument tracker, ENT                     | 1   |
| 9733534XOM | Patient tracker, ENT                        | 1   |
| 9734887XOM | Axiem™ non-invasive patient tracker         | 1   |
| 9731831XOM | Emitter cover, sterile                      | 1   |
| 9731132    | Dynamic reference (DRF) (EM, bone fixation) | 5   |
| 9732500XOM | Adhesive pads                               | 25  |
| 9734320XOM | Supplemental instrument tray                | 1   |



9734887XOM

## IGS software

| Product | Description                                   | Qty |
|---------|-----------------------------------------------|-----|
| 9735736 | StealthStation™ ENT software                  | 1   |
| 9735920 | StealthStation™ planning station ENT software | 1   |
| 9736155 | StealthMerge™ software                        | 1   |
| 9735742 | Stealth™ 3D software                          | 1   |

## IGS software support program

| Product           | Description                                                                      | Qty |
|-------------------|----------------------------------------------------------------------------------|-----|
| SRVC2-9732560XOM2 | StealthStation™ ENT service and support agreement                                | 1   |
| SERVICEH5-ENTMCT  | StealthStation™ ENT service agreement                                            | 1   |
| SSH5-SRVCFLEXENT  | StealthStation FlexENT™ navigation system service and support agreement          | 1   |
| SSH5-SRVCFLEXENT  | StealthStation FlexENT™ navigation system service and support advanced agreement | 1   |

## Planning station

| Product    | Description                              | Qty |
|------------|------------------------------------------|-----|
| 9735894KIT | StealthStation™ ENT planning station kit | 1   |

Set includes:

NOTE: The Medtronic family of navigation systems (StealthStation™ ENT, Fusion™, and FUSION™ Compact) will work with the instrument sets, balloons, suctions, and blades listed on these pages.


# ENT navigation solutions

|         |                                       |   |
|---------|---------------------------------------|---|
| 9735894 | Stealth S8 planning station           | 1 |
| 9736009 | Power pack PS Stealth S8 English only | 1 |

## NuVent™ EM sinus dilation system balloon dilation with integrated navigation

- For use with Medtronic ENT electromagnetic (EM) navigation
- Frontal, maxillary, and sphenoid designs
- 5 mm, 6 mm, and 7 mm balloon size options

| Product                   | Description                                           | Qty |
|---------------------------|-------------------------------------------------------|-----|
| <b>70 degree balloons</b> |                                                       |     |
| 1830517FRT70              | Balloon Seeker, EM 70 Degree Balloon 5.0 mm x 17.0 mm | 1   |
| 1830617FRT70              | Balloon Seeker, EM 70 Degree Balloon 6.0 mm x 17.0 mm | 1   |
| 1830717FRT70              | Balloon Seeker, EM 70 Degree Balloon 7.0 mm x 17.0 mm | 1   |
| <b>Frontal balloons</b>   |                                                       |     |
| 1830517FRT                | Balloon seeker, EM frontal balloon 5.0 mm x 17.0 mm   | 1   |
| 1830617FRT                | Balloon seeker, EM frontal balloon 6.0 mm x 17.0 mm   | 1   |
| 1830717FRT                | Balloon seeker, EM frontal balloon 7.0 mm x 17.0 mm   | 1   |
| <b>Maxillary balloons</b> |                                                       |     |
| 1830507MAX                | Balloon seeker, EM maxillary balloon 5.0 mm x 7.0 mm  | 1   |
| 1830607MAX                | Balloon seeker, EM maxillary balloon 6.0 mm x 7.0 mm  | 1   |
| 1830707MAX                | Balloon seeker, EM maxillary balloon 7.0 mm x 7.0 mm  | 1   |
| <b>Sphenoid balloons</b>  |                                                       |     |
| 1830517SPH                | Balloon seeker, EM sphenoid balloon 5.0 mm x 17.0 mm  | 1   |
| 1830617SPH                | Balloon seeker, EM sphenoid balloon 6.0 mm x 17.0 mm  | 1   |
| 1830717SPH                | Balloon seeker, EM sphenoid balloon 7.0 mm x 17.0 mm  | 1   |
| <b>Inflator</b>           |                                                       |     |
| 18INFKIT                  | Kit, inflator                                         | 1   |



NOTE: The Medtronic family of Navigation systems (StealthStation ENT, Fusion, and Fusion Compact) will work with the instrument sets, balloons, suctions, and blades listed on these pages.



# ENT navigation solutions

## Malleable suction instruments

- Plug-and-play with Medtronic ENT electromagnetic (EM) navigation
- Multiple lengths, diameters, and tips
- Standard suction design with built-in EM tracker

| Product | Description                  | Qty |
|---------|------------------------------|-----|
| 9735015 | Small standard tip, 14.5 cm  | 1   |
| 9735016 | Medium standard tip, 16.5 cm | 1   |
| 9735017 | Medium angle tip, 16.5 cm    | 1   |
| 9735018 | Medium ball tip, 16.5 cm     | 1   |
| 9735019 | Large standard tip 20 cm     | 1   |



## StraightShot™ automated EM tracking blades

| Product   | Description                                                             | Qty |
|-----------|-------------------------------------------------------------------------|-----|
| 1884012EM | RAD™ 12 curved rotatable blade with automated EM tracking, 4 mm         | 1   |
| 1884006EM | RAD 40 curved rotatable blade with automated EM tracking, 4 mm          | 1   |
| 1884080EM | Tricut™ 13 cm straight rotatable blade with automated EM tracking, 4 mm | 1   |
| 1895522   | Irrigation tubing                                                       | 5   |
| 1883080EM | Quadcut™ straight rotatable blade with automated EM tracking, 3.0 mm    | 1   |
| 1883480EM | Quadcut™ straight rotatable blade with automated EM tracking, 3.4 mm    | 1   |
| 1884380EM | Quadcut™ straight rotatable blade with automated EM tracking, 4.3 mm    | 1   |



## IGS instructional literature

| Product | Description            | Qty |
|---------|------------------------|-----|
| 9732379 | Imaging protocol sheet | 1   |

**Medtronic**

Engineering the extraordinary

ENT product catalog

# Nerve monitoring equipment

# NIM Vital™ nerve monitoring system

## NIM Vital™ nerve monitoring system

### System requirements

| Product  | Description                                | Qty |
|----------|--------------------------------------------|-----|
| NIM4CM01 | NIM Vital™ console                         | 1   |
| 1897821  | Console power cord North America, 6 meters | 1   |
| NIM4CPB1 | NIM Vital™ patient interface 4 channel     | 1   |
| NIM4CPB2 | NIM Vital™ patient interface cable         | 1   |

### Accessories

| Product   | Description                                                              | Qty |
|-----------|--------------------------------------------------------------------------|-----|
| NIM4CAM1  | NIM Vital™ mute probe adaptor                                            | 1   |
| NIMCAD400 | NIM incrementing probe adaptor                                           | 1   |
| NIM4CC01  | NIM Vital™ cart                                                          | 1   |
| NIM4CC02  | NIM Vital™ printer power isolator<br>100-120V power isolator for printer | 1   |
| NIM4CB01  | NIM Vital™ console battery                                               | 1   |
| NIM4CFU1  | NIM Vital™ console fuse kit                                              | 1   |
| 8253075   | NIM 3.0 INTFC 32MA fuse kit                                              | 1   |
| 8220325   | NIM muting probe                                                         | 1   |
| 8253600   | Patient simulator, universal                                             | 1   |

### NIM Vital™ nerve monitoring system compatible printers

- Samsung Xpress m2020w wireless black-and-white laser printer
- HP OfficeJet 200 mobile printer wireless color inkjet printer
- Samsung Xpress C430W



# EMG endotracheal tubes

## NIM TriVantage™ EMG endotracheal tube

- Recording electrode leads are twisted pair
- Packaged sterile with one green and one white subdermal needle
- Single use

| Product | ID     | OD      | Qty |
|---------|--------|---------|-----|
| 8229705 | 5.0 mm | 6.5 mm  | 1   |
| 8229706 | 6.0 mm | 8.2 mm  | 1   |
| 8229707 | 7.0 mm | 9.5 mm  | 1   |
| 8229708 | 8.0 mm | 10.7 mm | 1   |
| 8229735 | 5.0 mm | 6.5 mm  | 3   |
| 8229736 | 6.0 mm | 8.2 mm  | 3   |
| 8229737 | 7.0 mm | 9.5 mm  | 3   |
| 8229738 | 8.0 mm | 10.7 mm | 3   |



All measurements in millimeters.  
ID = Inner Diameter, OD = Outer Diameter

# Probes: monopolar stimulating

## Incrementing probe<sup>†</sup>

- Individual sterile packages
- Single use


| Product  | Description                                               | Qty |
|----------|-----------------------------------------------------------|-----|
| 8225825X | NIM™ surgeon control probe prass tip                      | 3   |
| 8225490X | NIM™ surgeon control probe ball tip                       | 3   |
| NIMDTP35 | NIM™ surgeon control long probe 1 mm ball tip, 35 cm long | 1   |



## Standard prass flush-tip monopolar stimulator probe

- Individual sterile packages
- Single use


| Product | Description               | Qty |
|---------|---------------------------|-----|
| 8225101 | Probes (tips and handles) | 5   |
| 8225110 | Tips only                 | 10  |



## Ball-tip monopolar stimulator probe

Sterile, single use

| Product | Description            | Qty |
|---------|------------------------|-----|
| 8225275 | Tip and handle, 1.0 mm | 1   |
| 8225277 | Tips only, 1.0 mm      | 10  |
| 8225276 | Tip and handle, 2.3 mm | 1   |
| 8225278 | Tips only, 2.3 mm      | 10  |



## Monopolar probe handles

- Individual sterile packages
- Single use

| Product | Description | Qty |
|---------|-------------|-----|
| 8225051 | Handles     | 10  |



## Yingling flex tip monopolar stimulator probe

- Packaged sterile with one green subdermal electrode
- Single use

| Product | Description      | Qty |
|---------|------------------|-----|
| 8225251 | Tips and handles | 5   |



<sup>†</sup>For use with NIM-Response™ and NIM-Neuro™ models 2.0 or later.

## Probes: bipolar stimulating

### Prass bipolar stimulator probe

Sterile, single use

| Product | Description | Qty |
|---------|-------------|-----|
| 8225451 | Probes      | 5   |



### Concentric bipolar stimulating probe

Sterile, single use

| Product | Description | Qty |
|---------|-------------|-----|
| 8225351 | Probes      | 5   |



### Side-by-side bipolar stimulating probe

Sterile, single use

| Product | Description | Qty |
|---------|-------------|-----|
| 8225401 | Probes      | 5   |




## Other disposables

### APS™ (Automatic Periodic Stimulation) electrode†

Sterile, single use

| Product | Description                   | Qty |
|---------|-------------------------------|-----|
| 8228052 | 2 mm, size range 2 mm to 3 mm | 1   |
| 8228053 | 3 mm, Size range 3 mm to 4 mm | 1   |



## Stimulus Instruments

- Individually packed. Single use.
- Connection cable included.

| Product | Description           | Qty |
|---------|-----------------------|-----|
| RIN201  | 3mm Ring Dissector    | 1   |
| CUR201  | Curved Needle         | 1   |
| ELE201  | Elevator              | 1   |
| DIS201  | Stimulating Dissector | 1   |



### Vari-Stim™ III Hand-Held Nerve Locator/Stimulator

- Battery-powered, single use
- Individual sterile packages

| Product | Description            | Qty |
|---------|------------------------|-----|
| 8562010 | Vari-Stim™ Stimulators | 10  |



†For use with NIM Vital, NIM-Response™ 3.0 and NIM-Neuro™ 3.0 models only.


‡For use with NIM™ models 2.0 and later, the Visao™ high-speed otologic drill and Integrated Power Console (IPC™) system.

# Electrodes: recording

## EMG hookwire electrodes

Sterile, single use

| Product | Description              | Qty   |
|---------|--------------------------|-------|
| 8226626 | Single electrodes, 25 mm | 6     |
| 8226326 | Paired electrodes, 25 mm | 1 set |



## Subdermal needle electrodes

- Packaged with one green and one white subdermal electrode
- Sterile, single use

| Product | Description      | Qty |
|---------|------------------|-----|
| 8227103 | 2-channel, 12 mm | 6   |



## Bulk-packed subdermal electrodes

- Packaged in 25 of one color or a mix of colors
- Sterile, lancet-tip needle electrodes
- Available in 10 different colors
- 12.0 mm needle electrodes with 5 ft. wire

| Product | Description                                                                        | Qty |
|---------|------------------------------------------------------------------------------------|-----|
| 8227131 | Needle electrodes, yellow                                                          | 25  |
| 8227133 | Needle electrodes, violet                                                          | 25  |
| 8227134 | Needle electrodes, green                                                           | 25  |
| 8227135 | Needle electrodes, black                                                           | 25  |
| 8227137 | Needle electrodes, white                                                           | 25  |
| 8227139 | Needle electrodes, gray                                                            | 25  |
| 8227140 | Needle electrodes, mixture 1<br>includes 5 red, 5 blue, 5 green, 5 white, 5 yellow | 25  |

# Electrodes: recording

## Subdermal needle twisted pair electrode set

- Packaged with one green and one white subdermal electrode
- Sterile, single use


| Product | Description      | Qty   |
|---------|------------------|-------|
| 8227421 | 4-channel, 22 mm | 1 set |
| 8227422 | 8-channel, 22 mm | 1 set |



## Paired subdermal electrodes

- Packaged with one green and one white subdermal electrode
- Sterile, single use


| Product | Description                             | Qty   |
|---------|-----------------------------------------|-------|
| 8227410 | 2-channel, 12 mm                        | 1 set |
| 8227411 | 4-channel, 12 mm                        | 1 set |
| 8227412 | 8-channel, 12 mm                        | 1 set |
| 8227465 | Mix 1: red, blue, orange, violet, gray  | 25    |
| 8227466 | Mix 2: yellow, green, tan, black, white | 25    |



## Prass paired electrodes

- Packaged with one green and one white subdermal electrode
- Sterile, single use


| Product | Description                 | Qty    |
|---------|-----------------------------|--------|
| 8227304 | 2-channel, 18 mm, 1 M cable | 5 sets |
| 8227301 | 2-channel, 25 mm, 1 M cable | 5 sets |
| 8227307 | 2-channel, 38 mm, 1 M cable | 5 sets |
| 8227500 | 4-channel, 18 mm, 2 M cable | 1 set  |
| 8227510 | 2-channel, 18 mm, 2 M cable | 5 sets |



## Prass paired electrodes with small hub

Sterile, single use

| Product | Description        | Qty    |
|---------|--------------------|--------|
| 8227414 | 2-channel, 12.5 mm | 5 sets |
| 8227415 | 4-channel, 12.5 mm | 5 sets |
| 8227418 | 2-channel, 18.0 mm | 5 sets |
| 8227419 | 4-channel, 18.0 mm | 5 sets |





**Medtronic**

Engineering the extraordinary

ENT product catalog

# PTeye<sup>TM</sup> parathyroid detection system

# PTeye™ parathyroid detection system

## PTeye™ parathyroid detection system

### System requirements

| Item number | Description                  | Qty  |
|-------------|------------------------------|------|
| PT102       | PTeye™ basic bundle          |      |
| PTeye       | PTeye™ console               | 1 ea |
| PTeye - 2   | PTeye™ foot pedal            | 1 ea |
| PTeye - 3   | PTeye™ power supply          | 1 ea |
| PTeye - 4   | PTeye™ power cord 125 V (US) | 1 ea |
| PTeye - 1   | PTeye™ single-use probe      | 1 ea |



### AVERTISSEMENT



INVISIBLE LASER RADIATION  
20mW Max., 0.125ms, 785nm.  
IEC 60825-1:2014, EN 60825-1:2014 + A11:2021.

ANSI/AAMI ES 60601-1:2005/AMD1:2012/AMD2:2021  
CAN/CSA-C22.2 No. 60601-1:14 (R2022)

# Ventilation tubes and myringotomy accessories

# Ventilation tubes and myringotomy accessories

## Dimensions

All dimensions are approximate and subject to normal manufacturing variance.

All measurements are listed as millimeters.

ID = inner diameter


IFD = inner flange diameter

OFD = outer flange diameter

IFDis = interflange distance

L = length

FL = flange length




# Ventilation tubes: temporary

## Armstrong beveled grommets

Designed by Beverly W. Armstrong, M.D., Charlotte, NC


| Material  | Feature | Color | RT  | Qty | ID   | IFD | Product |
|-----------|---------|-------|-----|-----|------|-----|---------|
| Microgel™ | Tab     | Green | INT | 5   | 1.14 | 3.5 | 25501   |
| Silicone  |         | Green | INT | 5   | 1.14 | 3.5 | 1010055 |



## Armstrong beveled grommets, modified

Designed by Beverly W. Armstrong, M.D., Charlotte, NC

| Material | Color | RT  | Qty | ID   | IFD | Product |
|----------|-------|-----|-----|------|-----|---------|
| FLPL     | White | INT | 5   | 1.14 | 3.5 | 1014242 |



## Armstrong beveled straight shanks


Designed by Beverly W. Armstrong, M.D., Charlotte, NC

| Material | Color | RT | Qty | ID   | IFD  | L    | Product |
|----------|-------|----|-----|------|------|------|---------|
| Silicone | Green | ST | 5   | 1.14 | 2.50 | 7.00 | 1010040 |



## Donaldson-type tubes

| Material | Feature | Color | RT | Qty | ID   | IFD | IFDis | Product |
|----------|---------|-------|----|-----|------|-----|-------|---------|
| FLPL     |         | Blue  | ST | 5   | 1.14 | 2.3 | 0.90  | 1015101 |
| Silicone |         | Green | ST | 5   | 1.14 | 2.3 | 0.90  | 1013020 |
| Silicone | Tab     | Green | ST | 5   | 1.14 | 2.3 | 0.90  | 1013015 |




RT = residence time, ST = short term (IFD of 0.0-3.0), INT = intermediate (IFD of 3.1-4.8), LT = long term (IFD of > 4.8)

# Ventilation tubes: temporary

## Goode T-grommets

Endorsed by Richard L. Goode, M.D., Palo Alto, CA


| Material | Color | RT | Qty | ID   | FL  | IFDis | Product |
|----------|-------|----|-----|------|-----|-------|---------|
| Silicone | Green | LT | 5   | 1.27 | 8.5 | 2.50  | 1016020 |



## Goode T-Tubes™


Designed by Richard L. Goode, M.D., Palo Alto, CA

| Material      | Color | RT | Qty     | ID   | FL  | L    | Product |
|---------------|-------|----|---------|------|-----|------|---------|
| Silicone      | Green | LT | 5       | 1.14 | 9.5 | 12.0 | 1016010 |
| Firm silicone | Green | LT | 5       | 1.14 | 9.5 | 12.0 | 1016035 |
| Firm silicone | Green | LT | 5 pairs | 1.14 | 6.0 | 6.0  | 24752   |



## Goode T-Tubes, modified

| Material | Color | RT | Qty | ID   | FL  | L   | Product |
|----------|-------|----|-----|------|-----|-----|---------|
| Silicone | Green | LT | 5   | 1.27 | 5.8 | 6.7 | 24771   |




RT = residence time, ST = short term (IFD of 0.0-3.0), INT = intermediate (IFD of 3.1-4.8), LT = long term (IFD of > 4.8)

# Ventilation tubes: temporary

## Paparella-type vent tubes

Designed by Michael M. Paparella, M.D., Minneapolis, MN


| Material  | Feature     | Color | RT  | Qty | ID   | IFD  | IFDis | Product |
|-----------|-------------|-------|-----|-----|------|------|-------|---------|
| Microgel™ | Notch & tab | Green | ST  | 5   | 1.02 | 2.50 | 1.00  | 25441   |
| Silicone  | Notch & tab | Green | ST  | 5   | 1.14 | 2.40 | 1.10  | 1025001 |
| Silicone  | Notch       | Green | ST  | 5   | 1.14 | 2.40 | 1.10  | 1025005 |
| Microgel™ | No feature  | Blue  | INT | 5   | 1.14 | —    | —     | 1082201 |



## Reuter bobbins

Designed by S. Harold Reuter, M.D., Houston, TX


| Material        | Feature  | Color    | RT | Qty | ID   | IFD | IFDis | Product |
|-----------------|----------|----------|----|-----|------|-----|-------|---------|
| FLPL            | No holes | Blue     | ST | 5   | 1.14 | 2.6 | 1.0   | 1010201 |
| Titanium        | No holes | Titanium | ST | 5   | 1.27 | 2.6 | 1.0   | 1027070 |
| Stainless steel | Holes    | Steel    | ST | 5   | 1.14 | 2.6 | 1.0   | 1027015 |



## Shah-type ventilation tubes

Designed by N. Shah, FRCS, DLO, London, England


| Material | Feature | Color | RT | Qty | ID   | IFD | IFDis | Product |
|----------|---------|-------|----|-----|------|-----|-------|---------|
| FLPL     | Mini    | White | ST | 5   | 0.76 |     | 1.00  | 14311   |
| FLPL     |         | Blue  | ST | 5   | 1.14 |     | 1.58  | 1018100 |
| FLPL     | Wire    | Blue  | ST | 5   | 1.14 |     | 1.58  | 1018101 |



## Sheehy-type collar buttons

Designed by James L. Sheehy, M.D., Los Angeles, CA

| Material  | Color | RT  | Qty | ID   | IFD | IFDis | Product |
|-----------|-------|-----|-----|------|-----|-------|---------|
| Microgel™ | Blue  | ST  | 5   | 1.27 | 2.8 | 1.32  | 1083301 |
| FLPL      | Blue  | INT | 5   | 1.27 | 3.0 | 1.55  | 1013301 |
| Silicone  | Green | ST  | 5   | 1.27 | 2.8 | 1.32  | 1028145 |




RT = residence time, ST = short term (IFD of 0.0-3.0), INT = intermediate (IFD of 3.1-4.8), LT = long term (IFD of > 4.8)

# Ventilation tubes: temporary

## Shepard grommets

Designed by Marvin G. Shepard, M.D., Dallas, TX


| Material | Feature | Color | RT | Qty | ID   | IFD  | Product |
|----------|---------|-------|----|-----|------|------|---------|
| FLPL     |         | Blue  | ST | 5   | 1.02 | 2.40 | 1016301 |
| FLPL     | Wire    | Blue  | ST | 5   | 1.02 | 2.40 | 1016501 |
| FLPL     | Wire    | Blue  | ST | 5   | 1.14 | 2.40 | 1016201 |
| FLPL     | Tab     | Blue  | ST | 5   | 1.14 | 2.40 | 1016101 |
| Silicone | Tail    | Green | ST | 5   | 1.14 | 2.30 | 1028025 |



## Soileau Tytan™ titanium ventilation tubes


Endorsed by James S. Soileau, M.D., Gonzales, LA

| Material       | Color    | RT | Qty     | ID   | IFD  | Product |
|----------------|----------|----|---------|------|------|---------|
| Tytan titanium | Titanium | ST | 5 pairs | 1.27 | 2.75 | 1066111 |




## Umbrella ventilation tubes

| Material | Feature | Color | RT  | Qty | ID   | IFD  | Product |
|----------|---------|-------|-----|-----|------|------|---------|
| Silicone | Tab     | Green | INT | 5   | 1.02 | 2.26 | 1016030 |



## V-T grommet ventilation tubes

| Material  | Color | RT  | Qty | ID   | IFD | IFDis | Product |
|-----------|-------|-----|-----|------|-----|-------|---------|
| Microgel™ | Blue  | INT | 5   | 1.27 | 4.5 | 1.3   | 1084401 |



RT = residence time, ST = short term (IFD of 0.0-3.0), INT = intermediate (IFD of 3.1-4.8), LT = long term (IFD of > 4.8)




# Ventilation tubes: permanent and specialty

## Pediatric ventilation tubes

1056101, 1056102, and 1056159 Endorsed by M. Coyle Shea, M.D., Memphis, TN  
1056108 Endorsed by William H. Moretz, Jr., M.D., Augusta, GA

| Feature   | Material             | Color    | RT | Qty | ID   | IFD | Product |
|-----------|----------------------|----------|----|-----|------|-----|---------|
|           | Tiny Tytan™ Titanium | Titanium | ST | 5   | 0.76 | 1.5 | 1056101 |
| Mini shah | FLPL                 | White    | ST | 5   | 0.76 | 2.0 | 14311   |



RT = residence time, ST = short term (IFD of 0.0-3.0), INT = intermediate (IFD of 3.1-4.8), LT = long term (IFD of > 4.8), P = permanent

ENT product catalog

# Ossicular reconstruction implants and accessories

# Partial ossicular reconstruction implants: titanium

## Scheer piston and wire

Fluoroplastic piston, stainless steel wire  
Endorsed by Harold F. Schuknecht, M.D.



| Product | Description                       | SD     | SL     | L      | Qty |
|---------|-----------------------------------|--------|--------|--------|-----|
| 1128190 | Stainless steel and fluoroplastic | 0.6 mm | 1.8 mm | 5.0 mm | 1   |



## Partial titanium prosthesis

Fully cannulated design

| Product | Description | HD     | L    | ID     | Qty |
|---------|-------------|--------|------|--------|-----|
| 1160000 | Titanium    | 2.5 mm | 2 mm | 1.1 mm | 1   |



All dimensions are approximate.  
LID = loop inner diameter    SD = shaft diameter    SL = shaft length    L = length

# Partial ossicular reconstruction implants: dense hydroxylapatite

## Incus partial prosthesis

- Grooved top
- Partially cannulated shaft


| Product | Description     | SID    | L       | Qty |
|---------|-----------------|--------|---------|-----|
| 1112075 | Hydroxylapatite | 1.1 mm | 3.75 mm | 1   |



## Notched offset partial

- Dense H/A cap
- Offset cap


| Product | Description     | SID     | HD   | L    | Qty |
|---------|-----------------|---------|------|------|-----|
| 0362    | Hydroxylapatite | 1.14 mm | 4 mm | 5 mm | 1   |



## Partial 90° off-center round

- Off-centered round head
- Fully cannulated shaft


| Product | Description     | SID    | HD   | L    | Qty |
|---------|-----------------|--------|------|------|-----|
| 1112095 | Hydroxylapatite | 1.1 mm | 4 mm | 5 mm | 1   |



## Partial 90° round

- Round head
- Fully cannulated shaft

| Product | Description     | SID    | HD   | L    | Qty |
|---------|-----------------|--------|------|------|-----|
| 1112090 | Hydroxylapatite | 1.1 mm | 4 mm | 5 mm | 1   |



All dimensions are approximate.


SID = shaft inner diameter    H = head    HD = head diameter    L = length

# Partial ossicular reconstruction implants: Flex H/A™

## Silverstein incus-stapes connector

Endorsed by Herbert Silverstein, M.D., Sarasota, FL

| Product | Description     | ID     | OD     | H    | Qty |
|---------|-----------------|--------|--------|------|-----|
| 1112120 | Hydroxylapatite | 1.1 mm | 2.1 mm | 1 mm | 1   |




## Causse notched offset partial with malleable link

- Offset, dense H/A cap with notch
- Flex H/A shaft may be trimmed with a knife
- Malleable titanium link

Endorsed by Jean-Bernard Causse, M.D., Béziers, France

| Product | Description                                | SID     | HD      | L    | Qty |
|---------|--------------------------------------------|---------|---------|------|-----|
| 0526    | Hydroxylapatite with Flex H/A and titanium | 1.14 mm | 3.25 mm | 5 mm | 1   |




## Causse-vincent notched narrow partial with short malleable link

- Offset, dense H/A cap with notch
- Flex H/A shaft may be trimmed with a knife
- Malleable titanium link

Endorsed by Jean-Bernard Causse, M.D., and Robert Vincent, M.D., Béziers, France

| Product | Description                                | SID     | HD      | L     | Qty |
|---------|--------------------------------------------|---------|---------|-------|-----|
| 0543    | Hydroxylapatite with Flex H/A and titanium | 1.14 mm | 3.25 mm | 14 mm | 1   |



## Centered partial, non-malleable

- Prosthesis is fully cannulated
- Flex H/A shaft may be trimmed with a knife

| Product | Description                  | SID     | HD   | L    | Qty |
|---------|------------------------------|---------|------|------|-----|
| 0524    | Hydroxylapatite and Flex H/A | 1.14 mm | 3 mm | 5 mm | 1   |



All dimensions are approximate.


ID = inner diameter    OD = outer diameter    SID = shaft inner diameter    H = head    HD = head diameter    L = length

## Partial ossicular reconstruction implants: Flex H/A™

### Malleus cradle partial with malleable link

- Dense H/A cap
- Titanium link
- Trimmable Flex H/A shaft


| Product | Description                                | SID     | H          | L    | Qty |
|---------|--------------------------------------------|---------|------------|------|-----|
| 0541    | Hydroxylapatite with Flex H/A and titanium | 1.14 mm | 2 x 4.4 mm | 5 mm | 1   |



### Notched offset partial with malleable link

- Dense H/A cap
- Cap is offset
- Trimmable Flex H/A shaft

| Product | Description                                | SID     | HD   | L    | Qty |
|---------|--------------------------------------------|---------|------|------|-----|
| 0570    | Hydroxylapatite with Flex H/A and titanium | 1.14 mm | 4 mm | 5 mm | 1   |




## Partial ossicular reconstruction implants: Flex H/A™ necrosis

### Flex H/A incus necrosis prostheses

Flex H/A material can be modified with knife to accommodate anatomical requirements

| Product | Description      | IID    | STID    | Qty |
|---------|------------------|--------|---------|-----|
| 0580    | Flex H/A, medium | 0.8 mm | 1.17 mm | 1   |



All dimensions are approximate.

SID = shaft inner diameter    H = head    HD = head diameter    L = length

# Partial ossicular reconstruction implants: Polycel™

## Austin modified partial

- Head is off-center
- Fully cannulated

Endorsed by David F. Austin, M.D., Idaho Falls, ID

| Product | Description | SID     | HD   | L       | Qty |
|---------|-------------|---------|------|---------|-----|
| 1156361 | Polycel™    | 1.17 mm | 3 mm | 4.75 mm | 1   |




## Causse fluoroplastic partial

- Polycel head
- Fluoroplastic shaft

Endorsed by Jean-Bernard Causse, M.D., Béziers, France

| Product | Description                | SID     | HD   | L       | Qty |
|---------|----------------------------|---------|------|---------|-----|
| 1156376 | Polycel™ and Fluoroplastic | 1.17 mm | 4 mm | 4.75 mm | 1   |



## Sheehy partial ossicular prosthesis (POP)

Endorsed by James L. Sheehy, M.D., Los Angeles, CA

| Product | Description | SID     | HD   | L       | Qty |
|---------|-------------|---------|------|---------|-----|
| 1156362 | Polycel™    | 1.17 mm | 3 mm | 4.75 mm | 1   |



All dimensions are approximate.

SID = shaft inner diameter    H = head    HD = head diameter    L = length

## Total ossicular reconstruction implants: dense hydroxylapatite

### Notched offset total with malleable link

- Flexible titanium link
- Notched cap
- Offset design

| Product | Description                  | SID    | HD   | L    | Qty |
|---------|------------------------------|--------|------|------|-----|
| 0363    | Hydroxylapatite and titanium | 0.8 mm | 4 mm | 7 mm | 1   |



## Total ossicular reconstruction implants: Flex H/A™

### Cause notched offset total with short malleable link

- Shortened malleable link allows trimming to short lengths
- Notched, dense H/A cap
- Trimmable Flex H/A shaft

Endorsed by Jean-Bernard Causse, M.D., Béziers, France

| Product | Description                                | SD   | HD      | L    | Qty |
|---------|--------------------------------------------|------|---------|------|-----|
| 0518    | Hydroxylapatite with Flex H/A and titanium | 1 mm | 3.25 mm | 9 mm | 1   |



All dimensions are approximate.

SID = shaft inner diameter    H = head    HD = head diameter    L = length



# Total ossicular reconstruction implants: Flex H/A™

## Centered total, non-malleable

- Dense H/A cap
- Non-malleable shaft, made totally of trimmable Flex H/A

| Product | Description                   | SD   | HD   | L    | Qty |
|---------|-------------------------------|------|------|------|-----|
| 0523    | Hydroxylapatite with Flex H/A | 1 mm | 3 mm | 9 mm | 1   |



## Lesinski offset total with malleable link

- Trimmable Flex H/A shaft

Designed by S. George Lesinski, M.D., Cincinnati, OH

| Product | Description                                | SD   | HD      | L    | Qty |
|---------|--------------------------------------------|------|---------|------|-----|
| 0530    | Hydroxylapatite with Flex H/A and titanium | 1 mm | 3.25 mm | 8 mm | 1   |



All dimensions are approximate.


SID = shaft inner diameter    H = head    HD = head diameter    L = length

# Universal ossicular reconstruction implants: Flex H/A™

## Bojrab universal prosthesis

Designed by Dennis Bojrab, M.D., Farmington Hills, MI

| Product | Description                         | PID     | HD   | PL      | TL   | Qty |
|---------|-------------------------------------|---------|------|---------|------|-----|
| 0506    | Hydroxylapatite with Flex H/A shaft | 1.14 mm | 3 mm | 3.65 mm | 7 mm | 1   |




## Head supports

### Bio-flek™ surgical head supports

Offset design endorsed by William H. Moretz, Jr., M.D., Augusta, GA

| Product | Size                                  | Qty |
|---------|---------------------------------------|-----|
| 1358225 | Standard, 50.8 mm x 22.9 cm (2" x 9") | 1   |



All dimensions are approximate.  
SID = shaft inner diameter    H = head    HD = head diameter    L = length

# Medtronic

Engineering the extraordinary



ENT product catalog

# Localized Drug Delivery (LDD)

# Localized Drug Delivery (LDD)


## PROPEL™ Mometasone Furoate Sinus Implants

| Product | Description                                                   | Qty |
|---------|---------------------------------------------------------------|-----|
| 70011   | PROPEL sinus implant                                          | 1   |
| 700115  | PROPEL sinus implant                                          | 5   |
| 60011   | PROPEL Mini sinus implant                                     | 1   |
| 600115  | PROPEL Mini sinus implant                                     | 5   |
| 60044-  | PROPEL Mini sinus implant with Straight Delivery System (SDS) | 1   |
| 60044-5 | PROPEL Mini sinus implant with Straight Delivery System (SDS) | 5   |
| 50011-  | PROPEL Contour sinus implant                                  | 1   |
| 50011-5 | PROPEL Contour sinus implant                                  | 5   |



## SINUVA™ (mometasone furoate) sinus implant

| NDC number  | Descriptor                           | Qty |
|-------------|--------------------------------------|-----|
| 10599000301 | SINUVA™ (mometasone furoate) 1350mcg | 1   |



### Ordering

SINUVA sinus implants can be obtained through two accessible avenues:

1. Specialty distributor: Providers can purchase SINUVA sinus implants directly from a distributor. Benefits verification and prior authorization from the payer is recommended prior to implanting a patient.
2. Specialty pharmacies: In some cases, SINUVA sinus implants can be filled under the pharmacy benefit by a select specialty pharmacy and delivered directly to the provider. Verifying the patient's benefits will establish the path to fulfillment.
3. For more information about procurement pathway, please reach out to: [ent.us.reimbursement@medtronic.com](mailto:ent.us.reimbursement@medtronic.com)

# Localized Drug Delivery (LDD)

## SINUVA™ (mometasone furoate) sinus implant

### INDICATION

SINUVA Sinus Implant is a corticosteroid-eluting implant indicated for the treatment of chronic rhinosinusitis with nasal polyps in adult patients  $\geq 18$  years of age who have had ethmoid sinus surgery.

### IMPORTANT SAFETY INFORMATION

#### CONTRAINDICATIONS

Patients with known hypersensitivity to mometasone furoate and any of the ingredients of the SINUVA Sinus Implant.

#### WARNINGS AND PRECAUTIONS

**Local Nasal Adverse Reactions:** Monitor nasal mucosa adjacent to the SINUVA Sinus Implant for any signs of bleeding (epistaxis), irritation, infection, or perforation. Avoid use in patients with nasal ulcers or trauma.

**Glaucoma and Cataracts:** Nasal steroids may result in development of glaucoma and/or cataracts. Glaucoma, cataracts, and clinically significant elevation of intraocular pressure were not observed in patients from the treatment group of one randomized controlled clinical study (N = 53) who underwent bilateral placement of SINUVA Sinus Implants. Close monitoring is warranted in patients with a change in vision or with a history of increased intraocular pressure, glaucoma, and/or cataracts.

**Hypersensitivity Reactions:** Hypersensitivity reactions, including rash, pruritus, angioedema have been reported with the use of corticosteroids.

**Immunosuppression and Risk of Infections:** Persons who are using drugs that suppress the immune system, such as corticosteroids, including SINUVA Sinus Implant are more susceptible to infections than healthy individuals. The safety and effectiveness of SINUVA Sinus Implant have not been established in pediatric patients less than 18 years of age and SINUVA is not indicated for use in this population. Corticosteroids should be used with caution, if at all, in patients with active or quiescent tuberculosis infection of the respiratory tract; untreated systemic fungal, bacterial, viral, or parasitic infections; or ocular herpes simplex.

**Hypercorticism and Adrenal Suppression:** If corticosteroid effects such as hypercorticism and adrenal suppression appear in patients, consider sinus implant removal.

#### ADVERSE REACTIONS

The most common adverse reactions observed ( $> 1\%$  of subjects) in clinical studies were asthma, headache, epistaxis, presyncope, bronchitis, otitis media, and nasopharyngitis.

#### POSTMARKETING EXPERIENCE

The following adverse reactions have been identified during post-approval use of the SINUVA sinus implant. These events include implant migration, lack of efficacy, nasal pain, headache, epistaxis.

**RX only.** Please see Full Prescribing Information for SINUVA at [SINUVA.com/hcp](https://www.medtronic.com/hcp)

# Medtronic

Engineering the extraordinary

ENT product catalog


# Tissue health

# Postoperative dressing

## Novapak™ nasal sinus packing and stent

- A single-use, dissolvable nasal packing and stent for use following sinus surgery
- Dimensions 8 cm x 1.8 cm x 1.3 cm


| Product   | Description                      | Qty |
|-----------|----------------------------------|-----|
| CS3600-10 | Novapak™ nasal packing, 10PK STD | 10  |
| CS3900-10 | Novapak™ nasal packing, 10PK FRM | 10  |
| CS3600-2  | Novapak™ nasal packing, 2PK STD  | 2   |
| CS3900-2  | Novapak™ nasal packing, 2PK FRM  | 2   |



# Nasal packing

## Merocel™ Doyle nasal dressing


| Product | Description                 | Length | Thickness | Height | Qty |
|---------|-----------------------------|--------|-----------|--------|-----|
| 400411  | Without string, with airway | 8 cm   | 1.5 cm    | 3 cm   | 10  |
| 440411  | With string, with airway    | 8 cm   | 1.5 cm    | 3 cm   | 10  |
| 400424  | Without string              | 8 cm   | 1.5 cm    | 3 cm   | 10  |
| 440424  | With string                 | 8 cm   | 1.5 cm    | 3 cm   | 10  |
| 400413  | Without string              | 8 cm   | 1 cm      | 3 cm   | 10  |
| 440413  | With string                 | 8 cm   | 1 cm      | 3 cm   | 10  |
| 440412  | With string, Pediatric      | 5.5 cm | 1 cm      | 2.5cm  | 10  |



## Merocel™ Kennedy intranasal surgical sponge (KISS) intra-op throat pack

Designed by David Kennedy, M.D., Philadelphia, PA



| Product | Description | Qty |
|---------|-------------|-----|
| 400515  | With string | 10  |



# Nasal packing


## Merocel™ 2000 nasal packing

| Product | Description              | Length | Thickness | Height | Qty |
|---------|--------------------------|--------|-----------|--------|-----|
| 470401  | With string              | 4.5 cm | 1.5 cm    | 2 cm   | 10  |
| 470404  | With string              | 8 cm   | 1.5 cm    | 2 cm   | 10  |
| 470408  | With string, with airway | 4.5 cm | 1.5 cm    | 2 cm   | 10  |
| 470405  | With String, with airway | 8 cm   | 1.5 cm    | 2 cm   | 10  |



## Merocel™ standard dressing

| Product | Description                                     | Length | Thickness | Height | Qty |
|---------|-------------------------------------------------|--------|-----------|--------|-----|
| 400400  | Without string                                  | 4.5 cm | 1.5 cm    | 2 cm   | 20  |
| 440400  | With string                                     | 4.5 cm | 1.5 cm    | 2 cm   | 10  |
| 400402  | Without string                                  | 8 cm   | 1.5 cm    | 2 cm   | 20  |
| 440402  | With string                                     | 8 cm   | 1.5 cm    | 2 cm   | 10  |
| 400403  | Without string, Single Compression              | 8 cm   | 1.5 cm    | 2 cm   | 20  |
| 400409  | Without string, with airway                     | 4.5 cm | 1.5 cm    | 2 cm   | 20  |
| 440409  | With string, with airway                        | 4.5 cm | 1.5 cm    | 2 cm   | 10  |
| 400410  | Without string, with airway                     | 8 cm   | 1.5 cm    | 2 cm   | 20  |
| 400405  | Without string, single compression, with airway | 8 cm   | 1.5 cm    | 2 cm   | 20  |
| 440410  | With string, with airway                        | 8 cm   | 1.5 cm    | 2 cm   | 10  |



## Merocel™ Pope epistaxis packing

Available with or without drawstring


| Product | Description    | Length | Thickness | Height | Qty |
|---------|----------------|--------|-----------|--------|-----|
| 400406  | Without string | 10 cm  | 1.5 cm    | 2.5 cm | 10  |
| 440406  | With string    | 10 cm  | 1.5 cm    | 2.5 cm | 10  |



## Merocel™ Pope Flex-Pak™

Flexible epistaxis packing with pinholes and drawstring

| Product | Description | Length | Thickness | Height | Qty |
|---------|-------------|--------|-----------|--------|-----|
| 460406  | With string | 10 cm  | 1.5 cm    | 2.5 cm | 10  |



All dimensions are approximate.



## Nasal packing

### Merocel™ Staxi-Stat™ pack

| Product | Description    | Qty |
|---------|----------------|-----|
| 400416  | Without string | 10  |

### Merocel™ Weimert epistaxis packing

| Product | Description | Length | Thickness | Height | Qty |
|---------|-------------|--------|-----------|--------|-----|
| 400421  | With string | 5.5 cm | 1.5 cm    | 2.5 cm | 20  |



400421


## Sinus packing

### Merocel™ 2000 Laminated Kennedy SinuStent™ dressing

| Product | Description                     | Length | Thickness | Height | Qty      |
|---------|---------------------------------|--------|-----------|--------|----------|
| 470530  | With string, ultra-slim profile | 3.5 cm | 0.6 cm    | 1.2 cm | 10 pairs |

### Merocel™ Baron sinus-pak

| Product | Description          | Length | Thickness | Height | Qty      |
|---------|----------------------|--------|-----------|--------|----------|
| 400407  | With string, regular | 2.5 cm | 1.2 cm    | 2 cm   | 10 pairs |




400407

### Merocel™ Kennedy sinus-pak

| Product | Description                     | Length | Thickness | Height | Qty      |
|---------|---------------------------------|--------|-----------|--------|----------|
| 400422  | With string, regular            | 3.5 cm | 1.2 cm    | 1.2 cm | 10 pairs |
| 400426  | With string, slim-profile       | 3.5 cm | 0.9 cm    | 1.2 cm | 10 pairs |
| 400530  | With string, ultra-slim profile | 3.5 cm | 0.6 cm    | 1.2 cm | 10 pairs |

### Merocel™ Oro-Pak intra-op throat pack

| Product | Description | Qty |
|---------|-------------|-----|
| 400520  | With string | 10  |



400520


All dimensions are approximate.

## Sinus packing

### Merocel™ Toffel sinus-pak kit with airway tube

- Sinus-pak (4.5 cm L x 1.3 cm T x 2.0 cm H)
- Nasal dressing (7.5 cm L x 1.5 cm T x 3.5 cm H)
- Airway tube (8 cm L)


| Product | Description | Qty     |
|---------|-------------|---------|
| 400423  | With string | 10 kits |



400423

### Merocel™ Toffel slim-profile sinus-pak

| Product | Description  | Length | Thickness | Height | Qty      |
|---------|--------------|--------|-----------|--------|----------|
| 400428  | Slim-profile | 4.5 cm | 0.9 cm    | 2 cm   | 10 pairs |




400428

## Ear wicks and packs

### Merocel™ ear packs, ambrus


| Product | Diameter | Length | Qty |
|---------|----------|--------|-----|
| 400354  | 15 mm    | 20 mm  | 20  |



400354

### Merocel™ micro spheres

| Product | Diameter | Thickness | Qty |
|---------|----------|-----------|-----|
| 400151  | 4 cm     | 0.7 cm    | 100 |
| 400152  | 5 cm     | 0.7 cm    | 100 |



### Merocel™ ear wicks, pope

| Product | Diameter | Length | Qty |
|---------|----------|--------|-----|
| 400141  | 9 mm     | 15 mm  | 50  |
| 425141  | 9 mm     | 15 mm  | 25  |
| 400146  | 9 mm     | 24 mm  | 10  |
| 400144  | 12 mm    | 24 mm  | 10  |




400141

All dimensions are approximate.

# Ear wicks and packs

## Merocel™ ear packs, Schindler

| Product | Description     | Diameter | Length | Qty |
|---------|-----------------|----------|--------|-----|
| 400355  | With drawstring | 12 mm    | 15 mm  | 10  |
| 400356  | With drawstring | 9 mm     | 13 mm  | 10  |



## Merocel™ Pope Oto-Wick™

| Product | Description | Size  | Qty |
|---------|-------------|-------|-----|
| 1425008 | Fenestrated | 15 mm | 50  |



# Epistaxis treatment products

## Epistat™ nasal catheter

- Two independently inflatable cuffs
- Hollow catheter
- Silicone construction
- Beveled posterior tip

| Product | Description             | Qty |
|---------|-------------------------|-----|
| 1527031 | Epistat™ nasal catheter | 1   |



All dimensions are approximate.

# Splints

## Doyle bi-valve septal splint

| Product | Qty    |
|---------|--------|
| 1524075 | 1 pair |



1524075

## Doyle II intranasal airway splints

- Integrated airway
- Medical grade silicone
- Hole in anterior tip

| Product | Qty     |
|---------|---------|
| 1524050 | 5 pairs |
| 1524055 | 1 pair  |



1524050

## External nasal splints

Lightweight aluminum lined with a thin layer of foam rubber

| Product | Description    | Qty |
|---------|----------------|-----|
| 1528110 | Single, small  | 25  |
| 1528111 | Single, small  | 5   |
| 1528120 | Single, medium | 25  |
| 1528121 | Single, medium | 5   |
| 1528137 | Single, large  | 25  |
| 1528138 | Single, large  | 5   |



All dimensions are approximate.


# Splints

## Reuter splint

- Constructed of fluoroplastic
- Pre-punched holes for easier suturing

Designed by S. Harold Reuter, M.D., Houston, TX

| Product | Description           | Qty     |
|---------|-----------------------|---------|
| 1527010 | Thin 0.25 mm          | 5 pairs |
| 1527015 | Thick 0.5 mm          | 5 pairs |
| 1527020 | Thick 0.5 mm, no slit | 5 pairs |
| 1527025 | Thin 0.25 mm, no slit | 5 pairs |
| 1527000 | Thin 0.25 mm          | 5 pairs |
| 1527005 | Thick 0.5 mm          | 5 pairs |




## Merocel™ silicone splint

- Septal splint folds
- Firm, yet flexible silicone
- Available to order individually or as kit

Kit includes:

- 2 Silicone septal splints (400561)
- 2 Merocel™ nasal packs

| Product | Description | Length | Thickness | Height | Qty     |
|---------|-------------|--------|-----------|--------|---------|
| 400561  | Splint      | 75 mm  | 25 mm     | 1 mm   | 5 pairs |
| 400562  | Kit         |        |           |        | 5 kits  |



All dimensions are approximate.

# Splints


## Simple splint

- Simple, inexpensive fluoroplastic splint
- Easily trimmed to fit

| Product | Length | Thickness | Height | Qty      |
|---------|--------|-----------|--------|----------|
| 1522000 | 70 mm  | 0.5 mm    | 44 mm  | 10 pairs |

## Thermaspilnts

| Product | Description | Length  | Height | Qty |
|---------|-------------|---------|--------|-----|
| 1529000 | Small       | 3.18 cm | 5.7 cm | 10  |
| 1529010 | Medium      | 3.9 cm  | 6.0 cm | 10  |
| 1529020 | Large       | 4.4 cm  | 6.7 cm | 10  |



1529010

# Septal buttons

## Nasal septal button

Medical silicone construction

| Product | Diameter | Qty |
|---------|----------|-----|
| 1524105 | 3 cm     | 1   |




1524105

## Two-Part adjustable septal button

- For non-surgical closure of nasal septal perforations
- Two-part button
- Medical silicone construction

Designed by Manning Goldsmith, III, M.D., Savannah, GA

| Product | Diameter | Qty |
|---------|----------|-----|
| 1524110 | 3 cm     | 1   |



1524110


All dimensions are approximate.

# Other products

## Instrument wipes

Removes fluid and debris from delicate instruments


| Product | Description                   | Qty |
|---------|-------------------------------|-----|
| 2218525 | 3.5" x 3.5" (8.9 cm x 8.9 cm) | 25  |



## Silicone sheeting

- Sheets are made of biocompatible silicone elastomer
- Sheets can be cut and trimmed as needed during procedures

| Product | Thickness                              | Qty |
|---------|----------------------------------------|-----|
| 1532500 | 5 cm x 5 cm x .25 mm (2" x 2" x .010") | 10  |
| 1532505 | 5 cm x 5 cm x .13 mm (2" x 2" x .005") | 10  |
| 1532510 | 5 cm x 5 cm x .5 mm (2" x 2" x .020")  | 10  |
| 1532520 | 5 cm x 5 cm x 1 mm (2" x 2" x .040")   | 10  |
| 1532530 | 5 cm x 5 cm x 1.6 mm (2" x 2" x .062") | 10  |



All dimensions are approximate.

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800-874-5797 or consult Medtronic's website at **medtronicent.com**.

# Medtronic

## ENT

6743 Southpoint Drive N  
Jacksonville, FL 32216  
USA  
Toll free: 800-874-5797  
Telephone: 904-296-9600  
Fax: 800-678-3995

©2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.™ Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.  
03/2025 - US-ENT-XXXX - [WF#13818251]