

Medtronic navigation comparisons



Hardware	StealthStation FlexENT™ system	Fusion Compact™ system	Fusion™ system (cart)
Cart	Optional Equipment	No	Integrated
Screen	24", multi touch	19", single touch	19", non-touch
Emitter	Under the head, or side emitter	Side emitter	Side emitter
Instrument ports	6	8	4
Instruments reusable	Yes	Yes	Yes
Footswitch	Optional Equipment	No	No
Wireless connectivity	No	No	No
Scan input	CD, USB, Wired Ethernet	CD, USB, Wired Ethernet	CD, USB, Wired Ethernet
Support balloon	Yes	Yes	Yes

Software	StealthStation FlexENT™ system	Fusion Compact™ system	Fusion™ system (cart)
Registration type	Touch, Trace, 3-Point Touch	Tracer, point merge	Tracer, point merge
Video recording	Yes	No	No
Surgical planning	Optional license required	No	No
Merge exams	Optional license required	Yes	Optional License - 4 exams
3D modelling	Optional license required	No	No
Virtual endoscopic views	Optional license required	No	No
DICOM Q/R	Included	Optional	Optional

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at **medtronicent.com**.

Medtronic

ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
Toll free: (800) 874-5797
Telephone: (904) 296-9600
Fax: (800) 678-3995

International telephone numbers

Europe HQ Switzerland 41-21-802-7000
Latin America HQ 305-500-9328
Adriatic Region 385-1-488-1120
Argentina 54-11-4898-5700
Australia 1-800-668-670
Baltic Region 37-1-67560226
Belgium 32-2456-09-09
Brazil 55-11-2182-9200
Canada 1-800-268-5346
Chile 56-2-2655-5110
China 86-10-5869-8989
Colombia 57-1-742-7300
Czech Republic 420-233-059-111
France 33-470-679-800
Germany 49-2159-8149-353
Greece 30-210-67-79-099
Hong Kong 852-2919-1300
Hungary 36-188-90600

India 91-22-33074700
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-3-6774-4611
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09
Malaysia 60-37-953-4800
Mexico 52-55-11-02-90-30
Netherlands 31-45-566-8800
Poland 48-22-4656900
Russian Fed. 7-495-580-73-77
Singapore 65-6436-5000
South Africa 27-11-260-9300
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
Thailand 662-232-7400
UK 44-1923-205-166

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC202103344a EN