


Medtronic

Malleable Suction Instruments

Bendable suction with built-in navigation



Re-shape a suction instrument
during image-guided surgery
without having to re-verify it.



Malleable Suction Instruments

Malleable Suction™ Instruments feature built-in AxiEM™ navigation technology within the instrument tip. Simply bend the instrument as necessary, and continue operating. There's no need to re-verify.

- Plug-and-play with Medtronic ENT electro-magnetic (EM) navigation
- Multiple lengths, diameters, and tips for precise localization and access
- Standard suction design with built-in EM tracker

Part Number	Description	Diameter	Length
9735015	Small standard tip	Approx. 7fr	14.5 cm
9735016	Medium standard tip	Approx. 9fr	16.5 cm
9735017	Medium angle tip	Approx. 9fr	16.5 cm
9735018	Medium ball tip	Approx. 9fr	16.5 cm
9735019	Large standard tip	Approx. 12fr	20 cm



Standard tip



Angle tip



Ball tip

Medtronic

ENT

6743 Southpoint Drive N
Jacksonville, FL 32216, USA
Toll free: (800) 874-5797
Telephone: (904) 296-9600
Fax: (800) 678-3995

International telephone numbers

Europe HQ Switzerland 41-21-802-7000
Latin America HQ 305-500-9328
Adriatic Region 385-1-488-1120
Argentina 54-11-4898-5700
Australia 1-800-668-670
Baltic Region 37-1-67560226
Belgium 32-2456-09-09
Brazil 55-11-2182-9200
Canada 1-800-268-5346
Chile 56-2-2655-5110
China 86-10-5869-8989
Colombia 57-1-742-7300

Czech Republic 420-233-059-111
France 33-470-679-800
Germany 49-2159-8149-209
Greece 30-210-67-79-099
Hong Kong 852-2919-1300
Hungary 36-188-90600
India 91-22-33074700
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-3-6774-4611
Korea 82-2-3404-3600
Lebanon 961-1-370-670

Luxembourg 32-2456-09-09
Malaysia 60-37-953-4800
Mexico 52-55-11-02-90-30
Netherlands 31-45-566-8800
Poland 48-22-4656900
Russian Fed. 7-495-580-73-77
Singapore 65-6436-5000
South Africa 27-11-260-9300
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
Thailand 662-232-7400
UK 44-1923-205-166



Plug-and-play with the StealthStation™ ENT Navigation System

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.



For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at www.medtronicent.com.