

# Medtronic navigation comparisons



| Hardware              | StealthStation FlexENT™ system  | Fusion Compact™ system  | Fusion™ system (cart)   |
|-----------------------|---------------------------------|-------------------------|-------------------------|
| Cart                  | Optional Equipment              | No                      | Integrated              |
| Screen                | 24", multi touch                | 19", single touch       | 19", non-touch          |
| Emitter               | Under the head, or side emitter | Side emitter            | Side emitter            |
| Instrument ports      | 6                               | 8                       | 4                       |
| Instruments reusable  | Yes                             | Yes                     | Yes                     |
| Footswitch            | Optional Equipment              | No                      | No                      |
| Wireless connectivity | No                              | No                      | No                      |
| Scan input            | CD, USB, Wired Ethernet         | CD, USB, Wired Ethernet | CD, USB, Wired Ethernet |
| Support balloon       | Yes                             | Yes                     | Yes                     |

| Software                    | StealthStation<br>FlexENT™ system | Fusion Compact™<br>system | Fusion™ system (cart)      |
|-----------------------------|-----------------------------------|---------------------------|----------------------------|
| Registration type           | Touch, Trace, 3-Point Touch       | Tracer, point merge       | Tracer, point merge        |
| Video recording             | Yes                               | No                        | No                         |
| Surgical planning           | Optional license required         | No                        | No                         |
| Merge exams                 | Optional license required         | Yes                       | Optional License - 4 exams |
| 3D modelling                | Optional license required         | No                        | No                         |
| Virtual endoscopic<br>views | Optional license required         | No                        | No                         |
| DICOM Q/R                   | Included                          | Optional                  | Optional                   |

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at **medtronicent.com**.

## Medtronic

### ENT

6743 Southpoint Drive N  
Jacksonville, FL 32216  
USA  
Toll free: (800) 874-5797  
Telephone: (904) 296-9600  
Fax: (800) 678-3995

### International telephone numbers

Europe HQ Switzerland 41-21-802-7000  
Latin America HQ 305-500-9328  
Adriatic Region 385-1-488-1120  
Argentina 54-11-4898-5700  
Australia 1-800-668-670  
Baltic Region 37-1-67560226  
Belgium 32-2456-09-09  
Brazil 55-11-2182-9200  
Canada 1-800-268-5346  
Chile 56-2-2655-5110  
China 86-10-5869-8989  
Colombia 57-1-742-7300  
Czech Republic 420-233-059-111  
France 33-470-679-800  
Germany 49-2159-8149-353  
Greece 30-210-67-79-099  
Hong Kong 852-2919-1300  
Hungary 36-188-90600

India 91-22-33074700  
Israel 972-9-972-4400  
Italy 39-02-24137-324  
Japan 81-3-6774-4611  
Korea 82-2-3404-3600  
Lebanon 961-1-370-670  
Luxembourg 32-2456-09-09  
Malaysia 60-37-953-4800  
Mexico 52-55-11-02-90-30  
Netherlands 31-45-566-8800  
Poland 48-22-4656900  
Russian Fed. 7-495-580-73-77  
Singapore 65-6436-5000  
South Africa 27-11-260-9300  
Spain 34-91-625-05-40  
Taiwan 886-2-2183-6000  
Thailand 662-232-7400  
UK 44-1923-205-166

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC202103344a EN