




TruClear™ hysteroscope set duo

# Integrate hysteroscopy into your office.



Perform see-and-treat hysteroscopic procedures with our expanded portfolio of TruClear™ hysteroscopes.

The latest to join the TruClear™ system family of products, the TruClear™ hysteroscope set duo provides the option of a diagnostic or operative sheath for the office setting.

It's the complete hysteroscopic offering you need to expand diagnostic and therapeutic procedures for your patients in the comfort and convenience of your office.

## Diagnostic benefits

- Easy insertion into the cervical canal and uterine cavity with slanted tip design<sup>1</sup>
- Angled view you need to see the uterine cavity with a 30° scope
- Evaluate the uterine cavity with a small 4-mm diagnostic sheath
- Little-to-no dilation with diagnostic sheath<sup>1,2</sup>

## Operative benefits


- Maintains clear operative field with built-in continuous flow<sup>2,3</sup>
- Access to the uterine cavity with long scope length
- Can be used with your preferred hand instrument†

## Practice benefits


- Perform see-and-treat hysteroscopic procedures with the flexibility of using either a diagnostic and operative sheath

†Hand instruments must be 5Fr diameter and 18 inches in length

## TruClear™ operative sheath duo




## TruClear™ diagnostic sheath duo



# Ordering information

| Product number | Description  | Unit of measure |
|----------------|--|-----------------|
| 8300057945     | TruClear™ Hysteroscope Duo   | Each            |
| 8300057946     | TruClear™ Diagnostic Sheath Duo  | Each            |
| 8300057947     | TruClear™ Operative Sheath Duo   | Each            |
| 8300055246     | Rubber Replacement Seal  | Each            |
| 8300059605     | TruClear™ Stopcock Replacement Set, Diagnostic (for use with 8300057946) | Each            |
| 8300059594     | TruClear™ Threaded Ring, Duo (for use with 8300057946)                   | Each            |
| 8300059604     | TruClear™ O-Ring, Duo (for use with 8300057946)                          | Pack of 5       |
| 8300054627     | Polishing Paste (for use with 8300057945)                                | Each            |
| 8300055288     | Small Stopcock Replacement Set (for use with 8300057947)                 | Each            |



Discover more ways to streamline your in-office procedures with the TruClear™ hysteroscopic tissue removal system.

For more information, call your Medtronic sales representative or visit [medtronic.com/truclear](https://medtronic.com/truclear).

1. Internal test report RE00440188 Rev Duo Insertion Claim Test Report, 26 January 2023  
2. Internal test report RE00192361 Dynamite Hysteroscope Evaluation Testing Report, 29 March 2019  
3. Internal test report RE00192361-1 Appendix A - OEM-Private Label Screening Test Results, 29 March 2019

©2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 02/2023 - US-UW-2300007- [WF#7114153]

555 Long Wharf Drive  
New Haven, CT 06511  
800.722.8772  
[medtronic.com/covidien](https://medtronic.com/covidien)

