

MR Statement for Axiom™ Detachable Coils

Axiom Detachable Coils can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss or less
- Maximum average specific absorption rate (SAR) of 3W/kg for 15 minutes of scanning

Test Report #TR07-198ADD1

CAUTION: Federal (USA) law restricts this device to sale distribution and use by or on order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

www.medtronic.com

ev3 Neurovascular

9775 Toledo Way
Irvine, CA 92618

Customer Service:
+ 1.800.684.6733



Patient Name: _____

Date of Procedure: _____

Treating Physician: _____

Physician Phone #: _____

Hospital: _____

Hospital Phone #: _____