





Specimen retrieval bags

Bag dimensions and volume	Product name	Order code	Bag material	Trocar size
 <p>Volume 220 ml</p> <p>2.28" 5.9"</p>	Endo Catch™ Gold device 10 mm	173050G	Polyurethane	11 mm or 12 mm
 <p>Volume 1500 ml</p> <p>3.83" 9.3"</p>	Endo Catch™ II device 15 mm	173049	Polyurethane	15 mm
 <p>Volume 150 ml</p> <p>3" 6"</p>	Endo Bag™ device	250-30	Polyurethane	10 mm, 11 mm, or 12 mm
 <p>Volume 750 ml</p> <p>5" 8"</p>	Endo Bag™ device	250-40	Polyurethane	10 mm, 11 mm, or 12 mm

 For more information, call 800-772-8772 or visit us at [medtronic.com/covidien](https://www.medtronic.com/covidien)