

**Medtronic**

McGRATH™ MAC video laryngoscope

# The right fit for your routine



# Inspired by insight – invested in safety

Better airway visibility drives better first-pass intubation success<sup>1</sup>. The McGRATH™ MAC video laryngoscope combines your core laryngoscopy skills with our evolved technology – so your first intubation attempt is your best.

The McGRATH™ MAC video laryngoscope provides<sup>1,2</sup>:

- Better subglottic visualization than traditional direct laryngoscopy
- Reduced airway complications
- Less hemodynamic stimulation than DL and certain VL

**The McGRATH™ MAC video laryngoscope delivers confidence, simplicity, and convenience.**

## Confidence in a familiar technique


- Combines line of sight video from its portrait display to retain your traditional laryngoscopy skills
- Direct line of sight for the healthcare provider while remaining upright increasing the distance from the patient's mouth<sup>3</sup>

## Simplicity in a versatile solution

- Single-button functionality with minute-by-minute battery indication and auto shut off feature to optimize battery life
- Multiple blade size options for pediatric to adult patients and routine extreme airways all with a single handle

## Convenience in one system

- Video laryngoscopy is cost-effective and recommended by The Society of Airway Management and Difficult Airway Society<sup>4</sup>
- Sterile, disposable blades and flexible handle cleaning options





## Ordering information

SKU No.	Description	Quantity
301-000-000	McGRATH™ MAC handle	1
350-082-000	McGRATH™ MAC blade size 1	10/box
350-072-000	McGRATH™ MAC blade size 1	50/box
350-084-000	McGRATH™ MAC blade size 2	10/box
350-017-000	McGRATH™ MAC blade size 2	50/box
350-086-000	McGRATH™ MAC blade size 3	10/box
350-005-000	McGRATH™ MAC blade size 3	50/box
350-088-000	McGRATH™ MAC blade size 4	10/box
350-013-000	McGRATH™ MAC blade size 4	50/box
X3-003-000	McGRATH™ MAC X blade size 3	10/box
340-000-000	Battery pack, 250-minute	1
300-000-100	Carry case	1

## Technical specifications

### McGRATH™ MAC handle

Size	180 mm x 68 mm x 110 mm
Weight	180 g
Power	3.6V lithium battery pack (c. 250 minutes)
Light source	High-intensity LED
Display	2.5" LCD color display
Camera	CMOS
Materials	Durable medical-grade thermoplastics composite with structural alloy core. The device and packaging are latex free.
Warranty	36 months

### McGRATH™ MAC disposable blades

Material	Fog-free, medical-grade, optical polymer
Packaging	Packaged sterile for single use

For demonstrations and quotes,  
please contact your local representative.  
Find out more at [medtronic.com/mac-first](http://medtronic.com/mac-first)



McGRATH™  
MAC X  
blade size 3



McGRATH™  
MAC  
blade size 4



McGRATH™  
MAC  
blade size 3



McGRATH™  
MAC  
blade size 2



McGRATH™  
MAC  
blade size 1



Battery pack



Carry case



For trained personnel only. For specific indications and instructions for use, please refer to the product manual.

#### References

1. Cook TM. Airway management guidance for the endemic phase of COVID-19. *Anaesthesia*. 2021 Feb;76(2):251-260. doi:10.1111/anese.15253.
2. El-Boghdady K, Wong DJN, Owen R, et al. Risks to healthcare workers following tracheal intubation of patients with COVID-19: a prospective international multicentre cohort study. *Anaesthesia*. 2020 Nov;75(11):1437-1447. doi: 10.1111/anae.15170. Epub 2020 Jul 9. PMID: 32516833; PMCID: PMC7300828.
3. Foley LJ, Urdaneta F, Berkow L, et al. Difficult Airway Management in Adult COVID-19 Patients: Statement by the Society of Airway Management. *Anesth Analg*. 2021 Mar 12. doi: 10.1213/ANE.0000000000005554. Epub ahead of print. PMID: 33711004.
4. Frerk C, Mitchell VS, McNarry AF, et al. Difficult Airway Society intubation guidelines working group. Difficult Airway Society 2015 guidelines for management of unanticipated difficult intubation in adults. *Br J Anaesth*. 2015 Dec;115(6):827-48. doi: 10.1093/bja/aev371. Epub 2015 Nov 10. PMID: 26556848; PMCID: PMC4650961.

[medtronic.com/mac-first](https://www.medtronic.com/mac-first)

©2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 03/2023 - US-RE-2200002 - [WF#2927481]

6135 Gunbarrel Avenue  
Boulder, CO 80301

800.635.5267

**Medtronic**