

Medtronic

Enhanced visualization. Every intubation.

The next generation McGRATH™ MAC video laryngoscope† features enhanced optics for clearer visibility and increased durability.

Inspired by you

Based on your feedback, we designed the McGRATH™ MAC video laryngoscope so your first attempt will be your best – to help avoid airway complications.

Designed for routine use

The McGRATH™ MAC video laryngoscope provides a simple and convenient solution for confident intubation.

New auto-off feature, decreases unnecessary battery drainage

Enhanced visualization with brighter, warmer LED

New camera to see more with the same size screen

Intelligent battery management – minute-by-minute power indicator for the exact remainder of battery capacity



Learn more about the next generation McGRATH™ MAC video laryngoscope.
[medtronic.com/mac-first](https://www.medtronic.com/mac-first)

Brighter and smarter

2x
light spread[†]

3x
brighter[†]

4x
resolution[†]



Enhanced design

Legacy design

Visualization		
Light color	Natural, warmer light	Bright white light
Light brightness	3X brighter light source ¹	Standard light source
White balance	Automatically adjusts	Fixed
Camera detail	Enhanced resolution	Standard resolution
Durability		
McGRATH™ MAC Handle	Reinforced screen hinge	Standard screen hinge
IP rating	IPX8	IPX7
Battery management		
	Auto-off feature	Manual power on/off feature
Ordering information		
SKU number	301-000-000	300-000-000

For trained personnel only. For specific indications and instructions for use, please refer to the product manual.

[†] May not be available for sale in all markets

[‡] As compared to the previous version of the McGrath™ MAC video laryngoscope

1. Based on internal test report #RE001500188, V&V plan.

©2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 03/2023 -US-RE-2200005-[WF#3429769]

[medtronic.com/mac-first](https://www.medtronic.com/mac-first)

Medtronic