



Enhance your ventilator experience with our service agreements software options add-on

Upgrade your Puritan Bennett™ 980 ventilators with enhanced capabilities tailored to your specific areas of care.

Our software option upgrades offer powerful tools that may help to improve ventilation and patient care. This software option add-on includes our four software options that do not require you to upgrade your hardware:

- **High-Flow Oxygen Therapy software:** enables effective high-flow oxygen delivery. Software provides flows of up to 50 L/min in neonatal and pediatric patients and up to 80 L/min in adult patients.
- **NeoMode 2.0 software[†]:** delivers specialized neonatal ventilation to patients weighing 300 grams to 7.5 kg.
- **NIV+ software:** monitors pressure at the patient interface to help reduce the uncertainty around pressure delivery to neonatal patients.
- **IE Sync™ software:** enhances patient-ventilator synchronization, that may help in improving treatment outcomes.

Now, access all four of these options for a special price of \$600, delivering a value of \$5,350.

This offer is available when you choose the software options add-on as part of the service agreement of your choice for the Puritan Bennett™ 980 ventilator.

To further enhance your PB980 ventilator systems, add the Integrated Nebulizer software option[‡] for an additional \$600, normally priced at \$1,800.

Integrated Nebulizer is designed to power and control the Aerogen™ nebulizer to deliver physician-prescribed medication for the patient.

[†] NeoMode 2.0 software does not include Proximal flow sensor. The proximal flow upgrade includes software and hardware, so it is excluded from this offer.

[‡] If you opt for the Integrated Nebulizer software option, it includes some hardware that will be installed during the software upgrade. Only available while supplies last.





Learn more

Start a new service plan or enhance your current one – contact us to find the right fit and unlock software upgrades for your PB980.

Contact our Service Agreements team at **800-255-6774 option 3** or **RS.ServiceAgreements@Medtronic.com**.

Clinical questions?

Contact our Clinical Support team at **ventilationclinicalsupport@medtronic.com**

Additional terms and conditions

- Choose from quarterly or monthly payment options to pay over time, and we'll upgrade your fleet with the selected software suite during your coverage term.†
- This add-on option can be added to service agreements of any duration, but payment for this add-on must be completed in full within the first year of coverage.
- The pricing for this offering is fixed and cannot be adjusted, even if any part of your fleet already includes some of the options listed above.
- Should your agreement be terminated early, the remaining balance of this add-on will become payable in full.

† Additional travel fees apply when not partnered with an existing contract-included preventive maintenance service event.

©2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.™™ Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 05/2025 - US-RE-2500046 - [WF#17306475]

[medtronic.com](https://www.medtronic.com)

Medtronic