

# Nebulization simplified

## Integrated nebulizer option for the Puritan Bennett™ 980 ventilator

The Puritan Bennett™ 980 ventilator includes an integrated nebulizer option – designed to power and control the Aerogen™ nebulizer, which delivers prescribed medication to patients.

It allows you to set, start, extend, and stop timed nebulization cycles directly on the ventilator's user interface – eliminating the need for the external controller and power supply.

On-screen options:

- Nebulization cycle time
  - Solo: 5-60 minutes – five minute increments, continuous
- Start, extend, or stop cycle



## Ordering information

Description	Part number
Integrated nebulizer option for Puritan Bennett™ 980 ventilator	980NEB
Aerogen™ controller cable	AG-AP1085
Aerogen™ Solo nebulizer, disposable (5/pack)	AG-AS3100
Aerogen™ Solo nebulizer, disposable (10/pack)	AG-AS3200

Please refer to instructions for use for the indications, contraindications, warnings, risks, and precautions associated with the device and the nebulizer option.

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. TM\* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 08/2022 – US-RE-2100159 – [WF#5196388]

6135 Gunbarrel Avenue  
Boulder, CO 80301

800.635.5267

