


Understanding tracheostomy tube sizing for optimal patient fit

Shiley™ flexible adult tracheostomy tubes

Meet the unique needs of your tracheostomy patients. Explore expanded Shiley™ flexible adult tracheostomy tube sizes and the updated visual design of the flange. The visual upgrades to the flange include: product information that’s organized in a more streamlined manner, increased font size that enhances readability, and etched markings to prevent information from rubbing off over time.



4	CN	65	A
Jackson Size	Cuff Type	ISO Size	Inner Cannula Type




Scan the QR code for product support resources.

Product code logic and sizing for Shiley™ flexible tracheostomy tubes

Shiley™ flexible tracheostomy tube with TaperGuard™ cuff

Disposable inner cannula product		Reusable inner cannula product†	Dimensions			
Tracheostomy tube SKU#	Inner cannula SKU#	Tracheostomy tube SKU#	I.D.	O.D.	Length	Inner cannula I.D.
4CN65A	4IC65	4CN65R	6.5 mm	9.4 mm	62 mm	5.5 mm
5CN70A	5IC70	5CN70R	7.0 mm	10.1 mm	68 mm	6.0 mm
6CN75A	6IC75	6CN75R	7.5 mm	10.8 mm	74 mm	6.5 mm
7CN80A	7IC80	7CN80R	8.0 mm	11.4 mm	77 mm	7.0 mm
8CN85A	8IC85	8CN85R	8.5 mm	12.2 mm	79 mm	7.5 mm
9CN90A	9IC90	9CN90R	9.0 mm	12.7 mm	79 mm	8.0 mm
10CN10A	10IC10	10CN10R	10.0 mm	13.8 mm	79 mm	9.0 mm



Shiley™ flexible tracheostomy tube cuffless

Disposable inner cannula product		Reusable inner cannula product†	Dimensions			
Tracheostomy tube SKU#	Inner cannula SKU#	Tracheostomy tube SKU#	I.D.	O.D.	Length	Inner cannula I.D.
4UN65A	4IC65	4UN65R	6.5 mm	9.4 mm	62 mm	5.5 mm
5UN70A	5IC70	5UN70R	7.0 mm	10.1 mm	68 mm	6.0 mm
6UN75A	6IC75	6UN75R	7.5 mm	10.8 mm	74 mm	6.5 mm
7UN80A	7IC80	7UN80R	8.0 mm	11.4 mm	77 mm	7.0 mm
8UN85A	8IC85	8UN85R	8.5 mm	12.2 mm	79 mm	7.5 mm
9UN90A	9IC90	9UN90R	9.0 mm	12.7 mm	79 mm	8.0 mm
10UN10A	10IC10	10UN10R	10.0 mm	13.8 mm	79 mm	9.0 mm



Scan the QR code for product support resources.

If you can't find the right size tube for your patient, please call our Customer Service department at 1-800-635-5267.

For trained personnel only. For specific indications and instructions for use, please refer to the IFU: PT00148006, PT00148007, PT00148008, PT00164499, PT00164528. Duration of use should not exceed twenty-nine (29) days. Please refer to the product manual for detailed usage and troubleshooting instructions.

©2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 07/2025 - US-RE-2500088 - [WF#18288141]

[medtronic.com/covidien](https://www.medtronic.com/covidien)

Medtronic