


A new look for Shiley™ flexible adult tracheostomy tubes

Same trusted product, reimagined for enhanced readability

Based on customer feedback, the visual design of the flange for Shiley™ flexible adult tracheostomy tubes has been updated. The updated flange gives you everything you need to get a better read on critical product information.




Explore the visual upgrades:

- Product information that's organized in a more streamlined manner
- Increased font size that enhances readability and is more discreet for patients
- Etched markings – instead of printed – to prevent information from rubbing off over time



Interpreting the REF

The product reference code, or SKU, featured on the Shiley™ flexible tracheostomy tube flange indicates the product size, whether the tube is cuffed or uncuffed, and if the inner cannula is disposable or reusable.



Product sizing

While the visual design of the flange changed, the product design, clinical product use, and product ordering information remain the same.

Shiley™ Flexible Adult Tracheostomy Tubes with Disposable Inner Cannula


Size		SKU#			Specifications				
Jackson	ISO	Cuffed	Cuffless	Inner Cannula	Outer Cannula Diameter (mm O.D.)	Outer Cannula Inner Diameter (mm I.D.)	Length (mm)	Cuff Resting Diameter (mm)	Inner Cannula Inner Diameter (mm I.D.)
4	6.5	4CN65H	4UN65H	4IC65	9.4	6.5	62	20.6	5.5
5	7.0	5CN70H	5UN70H	5IC70	10.1	7.0	68	23.0	6.0
6	7.5	6CN75H	6UN75H	6IC75	10.8	7.5	74	25.4	6.5
7	8.0	7CN80H	7UN80H	7IC80	11.4	8.0	77	25.4	7.0
8	8.5	8CN85H	8UN85H	8IC85	12.2	8.5	79	26.6	7.5
9	9.0	9CN90H	9UN90H	9IC90	12.7	9.0	79	27.6	8.0
10	10.0	10CN10H	10UN10H	10IC10	13.8	10.0	79	28.6	9.0

Shiley™ Flexible Adult Tracheostomy Tubes with Reusable Inner Cannula

4	6.5	4CN65R	4UN65R	9.4	6.5	62	20.6	5.5
5	7.0	5CN70R	5UN70R	10.1	7.0	68	23.0	6.0
6	7.5	6CN75R	6UN75R	10.8	7.5	74	25.4	6.5
7	8.0	7CN80R	7UN80R	11.4	8.0	77	25.4	7.0
8	8.5	8CN85R	8UN85R	12.2	8.5	79	26.6	7.5
9	9.0	9CN90R	9UN90R	12.7	9.0	79	27.6	8.0
10	10.0	10CN10R	10UN10R	13.8	10.0	79	28.6	9.0

Ready for release in a phased rollout

Over the coming year, we'll begin a phased rollout of Shiley™ flexible adult tracheostomy products featuring the updated flange. Fulfillment from our two manufacturing facilities, existing product in stock with distributors, and each facility's existing inventory may dictate when the new products will reach clinical areas. All facilities should receive products with the updated flange by early 2023.



Visit [medtronic.com/ShileyFlex](https://www.medtronic.com/ShileyFlex) to learn more.
Or scan the QR code.

©2022 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 03/2022 - US-ENT-2200002 - [WF#6383262]

6135 Gunbarrel Avenue
Boulder, CO 80301
800.635.5267
[medtronic.com/covidien](https://www.medtronic.com/covidien)

Medtronic