

Trusted brand. Advanced design.



Shiley™ neonatal and pediatric tracheostomy tubes

Shiley™ neonatal and pediatric tracheostomy solutions incorporate a range of clinician-inspired design advances – tailored to meet the needs of neonatal and pediatric patients.

Flange distinction for neonatal or pediatric patients



Comfort recess on the bottom of the connector to help reduce inferior stoma irritation¹



TaperGuard™ cuff may reduce removal force¹



Transparent soft flange for visualization of skin below



Neonatal tracheostomy tubes

Product Sku	Description	I.D.	O.D.	Length
2.5NEF	Cuffless	2.5 mm	4.2 mm	28 mm
3.0NEF		3.0 mm	4.8 mm	30 mm
3.5NEF		3.5 mm	5.4 mm	32 mm
4.0NEF		4.0 mm	6.0 mm	34 mm
4.5NEF		4.5 mm	6.7 mm	36 mm
2.5NCF	With TaperGuard™ cuff	2.5 mm	4.2 mm	28 mm
3.0NCF		3.0 mm	4.8 mm	30 mm
3.5NCF		3.5 mm	5.4 mm	32 mm
4.0NCF		4.0 mm	6.0 mm	34 mm
4.5NCF		4.5 mm	6.7 mm	36 mm

Pediatric tracheostomy tubes

Product Sku	Description	I.D.	O.D.	Length
2.5PEF	Cuffless	2.5 mm	4.2 mm	38 mm
3.0PEF		3.0 mm	4.8 mm	39 mm
3.5PEF		3.5 mm	5.4 mm	40 mm
4.0PEF		4.0 mm	6.0 mm	41 mm
4.5PEF		4.5 mm	6.7 mm	42 mm
5.0PEF		5.0 mm	7.3 mm	44 mm
5.5PEF	Extra long, cuffless	5.5 mm	7.9 mm	46 mm
5.0PELF		5.0 mm	7.3 mm	50 mm
5.5PELF		5.5 mm	7.9 mm	52 mm
6.0PELF		6.0 mm	8.5 mm	54 mm
6.5PELF	With TaperGuard™ cuff	6.5 mm	9.0 mm	56 mm
2.5PCF		2.5 mm	4.2 mm	38 mm
3.0PCF		3.0 mm	4.8 mm	39 mm
3.5PCF		3.5 mm	5.4 mm	40 mm
4.0PCF		4.0 mm	6.0 mm	41 mm
4.5PCF		4.5 mm	6.7 mm	42 mm
5.0PCF	Extra long with TaperGuard™ cuff	5.0 mm	7.3 mm	44 mm
5.5PCF		5.5 mm	7.9 mm	46 mm
5.0PLCF		5.0 mm	7.3 mm	50 mm
5.5PLCF		5.5 mm	7.9 mm	52 mm
6.0PLCF	Extra long with TaperGuard™ cuff	6.0 mm	8.5 mm	54 mm
6.5PLCF		6.5 mm	9.0 mm	56 mm

For information or to place an order, contact your local Medtronic representative or call 800-722-8772.



[LivingwithaTrach.com](https://www.livingwithatrach.com)

Access educational videos and home care support for tracheostomy patients



[Medtronic.com/Covidien](https://www.Medtronic.com/Covidien)

View product information and clinical education opportunities

For trained personnel only. For specific indications and instructions for use, please refer to the product manual.

1. Based on internal testing. Comparative ventilator air leak test performed using the Shiley™ 4.0PCF and 6.5PLCF pediatric tracheostomy tubes with taper-shaped cuff vs. predicate Shiley™ 4.0PDC and 6.5PLC pediatric tracheostomy tubes with barrel-shaped cuff.

©2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 01/2023-US-RE-2300091-[WF#3901073]

6135 Gunbarrel Avenue Boulder, CO 80301 800.635.5267 [medtronic.com/covidien](https://www.Medtronic.com/covidien)

Medtronic