

Strong sutures. Confident closure.

By combining greater pliability^{†,1} with 3-week-long tensile strength and absorption in 90 to 110 days,^{1,2} Trexon[™] monofilament synthetic absorbable sutures embody the next evolution in Medtronic wound closure offerings.

Strength when your patients need it most.

- Superior strength during the critical early stages of wound healing compared to the competition^{‡,1}
- Equivalent or superior tensile strength out of the package than legacy Biosyn[™] sutures^{§,1}
- Superior tensile strength out of the package compared to Monocryl^{™*} and Monocryl^{™*} Plus sutures^{¶,1}

[†] Compared to legacy Biosyn[™] sutures. [‡] Benchmark testing may not be indicative of clinical performance. As compared to Monocryl^{™*} Plus sutures and Monocryl^{™*} sutures. [§] For sizes 2/0, 4/0, and 6/0. [¶] For sizes 1, 2/0, and 4/0. ^{||} For sizes 1, 4/0, and 6/0.

Feel the difference of superior pliability.^{†,1}

- More pliable than legacy Biosyn[™] sutures¹
- Equivalent or superior pliability compared to Monocryl^{™*} sutures^{||,1} and Monocryl^{™*} Plus sutures¹

Unpack a memorable suture performance.

Trexon[™] sutures offer less memory out of the package than Biosyn[™] sutures,¹ with an updated packaging design aligned to industry standards and surgeons' needs.



Medtronic

Ordering information

Choose a suture made for high standards.

Order code	Product description	Size	Color	Strands	Length (in)	Box quantity	Needle type
CY811P	Trexon™ 2/0 30" VL GS-21	2/0	Violet	1	30	36	GS-21
CY843P	Trexon™ 2/0 30" UD GS-21	2/0	Undyed	1	30	36	GS-21
CY883P	Trexon™ 2/0 30" VL GS-22	2/0	Violet	1	30	36	GS-22
CY904P	Trexon™ 0 36" VL GS-25	0	Violet	1	36	36	GS-25
CY953P	Trexon™ 2/0 36" UD GS-21	2/0	Undyed	1	36	36	GS-21
CY954P	Trexon™ 0 36" UD GS-21	0	Undyed	1	36	36	GS-21
GY121P	Trexon™ 4/0 30" VL V-20	4/0	Violet	1	30	36	V-20
GY122P	Trexon™ 3/0 30" VL V-20	3/0	Violet	1	30	36	V-20
GY123P	Trexon™ 2/0 30" VL V-20	2/0	Violet	1	30	36	V-20
GY321P	Trexon™ 4/0 30" UD V-20	4/0	Undyed	1	30	36	V-20
GY322P	Trexon™ 3/0 30" UD V-20	3/0	Undyed	1	30	36	V-20
GY323P	Trexon™ 2/0 30" UD V-20	2/0	Undyed	1	30	36	V-20
SY5627P	Trexon™ 4/0 18" UD P-12	4/0	Undyed	1	18	12	P-12
SY5628P	Trexon™ 3/0 18" UD P-12	3/0	Undyed	1	18	12	P-12
SY5637P	Trexon™ 4/0 30" UD P-12	4/0	Undyed	1	30	12	P-12
SY5638P	Trexon™ 3/0 30" UD P-12	3/0	Undyed	1	30	12	P-12
SY5678P	Trexon™ 4/0 30" UD P-14	4/0	Undyed	1	30	36	P-14
SY5679P	Trexon™ 4/0 18" UD P-14	4/0	Undyed	1	18	24	P-14
SY5679GP	Trexon™ 3/0 30" UD P-14	3/0	Undyed	1	30	12	P-14
SY5687P	Trexon™ 5/0 18" UD P-13	5/0	Undyed	1	18	12	P-13
SY5690P	Trexon™ 4/0 18" UD P-13	4/0	Undyed	1	18	12	P-13
SY632P	Trexon™ 3/0 30" VL C-23	3/0	Violet	1	30	36	C-23
SY633P	Trexon™ 4/0 30" VL C-23	4/0	Violet	1	30	36	C-23
SY691P	Trexon™ 4/0 30" UD C-13	4/0	Undyed	1	30	36	C-13
SY693P	Trexon™ 3/0 30" UD C-14	3/0	Undyed	1	30	36	C-14
SY791P	Trexon™ 4/0 30" VL C-13	4/0	Violet	1	30	36	C-13
SY822P	Trexon™ 3/0 30" UD C-13	3/0	Undyed	1	30	36	C-13
SY823P	Trexon™ 2/0 30" UD C-14	2/0	Undyed	1	30	36	C-14
SY922P	Trexon™ 3/0 30" VL C-13	3/0	Violet	1	30	36	C-13
SY923P	Trexon™ 2/0 30" VL C-14	2/0	Violet	1	30	36	C-14
UY202P	Trexon™ 5/0 30" VL CV-23	5/0	Violet	1	30	36	CV-23
UY203P	Trexon™ 4/0 30" VL CV-23	4/0	Violet	1	30	36	CV-23
UY204P	Trexon™ 3/0 30" VL CV-23	3/0	Violet	1	30	36	CV-23
UY214P	Trexon™ 4/0 30" UD CV-23	4/0	Undyed	1	30	36	CV-23

Important: Always refer to the instructions for use (IFU) supplied with products for complete instructions, indications, contraindications, warnings, and precautions.

1. Based on internal engineering report RE00582176 Rev B, Trexon™ claims report. April 2026. 2. Trexon™ monofilament synthetic absorbable suture [instructions for use]. Mansfield, MA: Medtronic; 2025.

©2026 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

All other brands are trademarks of a Medtronic company. 05/2026 - US-WC-2600011 - [WF# 21656062]

[medtronic.com/trexon](https://www.medtronic.com/trexon)

Medtronic