

Medtronic

**WOUND CLOSURE
PRODUCT CATALOGUE**



Veterinary ophthalmic and microsurgical sutures



Animal health Europe




Contents

● Velosorb™ Fast	4
● Polysorb™	6
● Maxon™	10
● Surgipro™ / Surgipro™ II	12
● Ti-Cron™ / Surgidac™	16
● Dermalon™ / Monosof™	18
● Novafil™	22
● Sofsilk™	24
● Needles chart	27


How to read a suture label



Taper point

Needles specifications	Size		Length cm	Color	Units per box	Comments	Reorder code
	USP	EP					
CV-1 3/8 Circle 9 mm 	8-0	0.4	30	Violet	12		CV-9916G

Spatula

Needles specifications	Size		Length cm	Color	Units per box	Comments	Reorder code
	USP	EP					
SE-140-8 3/8 Circle 6.3 mm Curvature 140° 	8-0	0.4	30	Violet	12		CV-501G
	7-0	0.5	30	Violet	12		CV-500G

Polysorb™

Braided absorbable suture

**Strong during the critical
wound healing period***






*Results indicative for human trials.


1. van Ramshorst GH, Nieuwenhuizen J, Hop WC, Arends P, Boom J, Jeekel J, Lange JF. Abdominal wound dehiscence in adults: development and validation of a risk model. *World J Surg.* 2010 Jan;34(1):20-7. doi: 10.1007/s00268-009-0277-y. PMID: 19898894; PMCID: PMC2795859.

Please refer to the individual package IFUs for complete instructions, indications, contraindications, warnings and precautions.








Reverse cutting

Needles specifications		Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
HE-6 3/8 Circle 7.4 mm Double-arm Curvature 103°		6-0	0.7	45	Violet	12	Surgalloy™	L-1551-K
		C-1	1	45	Undyed	36	SL-171	
C-1 3/8 Circle 12 mm		6-0	0.7	45	Undyed	36	SL-171	
		5-0	1	45	Undyed	36	SL-172	
HE-3 3/8 Circle 13 mm Double-arm Curvature 155°		5-0	1	30	Undyed	12	Surgalloy™	L-4323K




Taper point

Needles specifications		Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
MV-135-5 3/8 Circle 6 mm		8-0	0.4	13	Violet	12	Surgalloy™	L-2800

Spatula

Needles specifications	Size	Length cm	Color	Units per box	Comments	Reorder code		
							USP	EP
SE-100-8 1/4 Circle 6.3 mm Double-arm Curvature 100°		6-0	0.7	45	Violet	12	Surgalloy™	L-1750K
								
SS-29 1/4 Circle 7.4 mm Double-arm Curvature 103°		6-0	0.7	30	Violet	12	Surgalloy™	L-1756K
		6-0	0.7	45	Violet	12	Surgalloy™	L-1755K
SS-14 1/4 Circle 8.6 mm Double-arm Curvature 115°		6-0	0.7	45	Violet	12	Surgalloy™	L-1791K
		6-0	0.7	45	Undyed	12	Surgalloy™	L-1670-K
		5-0	1	45	Violet	12	Surgalloy™	L-1792K
		5-0	1	45	Undyed	12	Surgalloy™	L-1669-K
SS-24 1/4 Circle 8.6 mm Double-arm Curvature 90°		6-0	0.7	20	Violet	12	Surgalloy™	L-1740-K
		6-0	0.7	30	Violet	12	Surgalloy™	L-1742K
		6-0	0.7	30	Undyed	12	Surgalloy™	L-1745K
		6-0	0.7	45	Undyed	24	Surgalloy™	L-2752K
		5-0	1	20	Violet	12	Surgalloy™	L-1739-K
		5-0	1	30	Violet	12	Surgalloy™	L-1785K
		5-0	1	45	Undyed	24	Surgalloy™	L-2753K
SE-140-8 3/8 Circle 6.3 mm Curvature 140°		8-0	0.4	30	Violet	12	Surgalloy™	L-2750K
		7-0	0.5	30	Violet	12	Surgalloy™	L-2751K
SE-140-8 3/8 Circle 6.3 mm Double-arm Curvature 140°		8-0	0.4	20	Violet	12	Surgalloy™	L-2747K
		8-0	0.4	30	Violet	12	Surgalloy™	L-2748K
		8-0	0.4	45	Violet	12	Surgalloy™	L-5118K
		7-0	0.5	30	Violet	12	Surgalloy™	L-1796K
		7-0	0.5	45	Violet	12	Surgalloy™	L-1897K

Spatula

Needles specifications		Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
SS-28 1/2 Circle 7.4 mm Double-arm Curvature 145°		6-0	0.7	45	Violet	12	Surgalloy™	L-1562-K
SE-160-8 1/2 Circle 5.5 mm Double-arm Curvature 160°		8-0	0.4	20	Violet	12	Surgalloy™	L-2574K
		8-0	0.4	30	Violet	12	Surgalloy™	L-2749K
		7-0	0.5	45	Violet	12	Surgalloy™	L-1787K
SS-2 1/2 Circle 8.6 mm Double-arm Curvature 200°		5-0	1	30	Undyed	12	Surgalloy™	L-0223K

Maxon™

Monofilament absorbable suture



Please refer to the individual package IFUs for complete instructions, indications, contraindications, warnings and precautions.

Spatula

Needles specifications






Size		Length cm	Color	Units per box	Comments	Reorder code
USP	EP					
7-0	0.5	45	Green	36		8886662801

DGE-6




3/8 Circle 8 mm







Taper point

Needles specifications	Size	Length cm	Color	Units per box	Comments	Reorder code	
							USP
MV-70-3 3/8 Circle 4 mm 	10-0	0.2	13	Blue	12	Surgalloy™	P-2790G
MV-100-4 3/8 Circle 5 mm 	10-0	0.2	13	Blue	12	Surgalloy™	P-2794G
MV-100-4 3/8 Circle 5 mm Double-arm 	10-0	0.2	5	Blue	12	Surgalloy™	P-2795G
MV-135-5 3/8 Circle 6 mm Double-arm 	8-0	0.4	45	Blue	36	Surgalloy™	VP-900-X
	8-0	0.4	5 × 45	Blue	12	Surgalloy™, Long pack	VP-901-MX
	8-0	0.4	2 × 45	Blue	12	Surgalloy™, Long pack	VP-900-MX
	8-0	0.4	60	Blue	36	Surgalloy™	VP-902-X
	8-0	0.4	5 × 60	Blue	12	Surgalloy™, Long pack	VP-903-MX
	8-0	0.4	2 × 60	Blue	12	Surgalloy™, Long pack	VP-902-MX
	8-0	0.4	75	Blue	12	Surgalloy™, Long pack	VP-973-X
	7-0	0.5	60	Blue	36	Surgalloy™	VP-904-X
	7-0	0.5	2 × 60	Blue	12	Surgalloy™, Long pack	VP-904-MX
MVF-135-5 3/8 Circle 6 mm 	8-0	0.4	2 × 60	Blue	12	Surgalloy™, Long pack	VP-209MX

Taper point

Needles specifications	Size		Length cm	Color	Units per box	Comments	Reorder code	
	USP	EP						
MV-175-8 3/8 Circle 8 mm Double-arm		8-0	0.4	45	Blue	36	Surgalloy™	VP-744-X
		8-0	0.4	2 × 45	Blue	12	Surgalloy™, Long pack	VP-744-MX
		8-0	0.4	60	Blue	36	Surgalloy™	VP-745-X
		8-0	0.4	2 × 60	Blue	12	Surgalloy™, Long pack	VP-743-MX
		8-0	0.4	2 × 75	Blue	12	Surgalloy™, Long pack	VP-753-MX
		7-0	0.5	2 × 45	Blue	12	Long pack	VP-736-M
		7-0	0.5	45	Blue	36	Surgalloy™	VP-628
		7-0	0.5	5 × 45	Blue	12	Surgalloy™, Long pack	VP-734-MX
		7-0	0.5	60	Blue	36	Surgalloy™	VP-735-X
		7-0	0.5	2 × 60	Blue	12	Surgalloy™, Long pack	VP-737-MX
		7-0	0.5	5 × 60	Blue	12	Surgalloy™, Long pack	VP-735-MX
		7-0	0.5	75	Blue	36	Surgalloy™	VP-747-X
		7-0	0.5	2 × 75	Blue	12	Surgalloy™, Long pack	VP-747-MX
MVF-175-8 3/8 Circle 8 mm Double-arm		8-0	0.4	2 × 60	Blue	12	Surgalloy™	VPF-738-MX
		7-0	0.5	60	Blue	36	Surgalloy™	VPF-735-X
MV-175-9 3/8 Circle 10 mm Double-arm		8-0	0.4	60	Blue	36	Surgalloy™	VP-763-X
		7-0	0.5	5 × 75	Blue	12	Surgalloy™, Long pack	VP-875-MX
		7-0	0.5	75	Blue	36	Surgalloy™	VP-875-X
		7-0	0.5	5 × 60	Blue	12	Surgalloy™, Long pack	VP-759-MX

Spatula

Needles specifications	Size		Length cm	Color	Units per box	Comments	Reorder code	
	USP	EP						
SE-140-8 3/8 Circle 6.3 mm Double-arm Curvature 140°		10-0	0.2	30	Blue	12	Surgalloy™	P-2771K
								
SE-160-6 1/2 Circle 5.5 mm Double-arm Curvature 160°		10-0	0.2	30	Blue	12	Surgalloy™	P-2756K
								

Ti-Cron™ / Surgidac™

Coated and uncoated braided polyester

Ti-Cron™ / Surgidac™ is a suture brand you have trusted for over 35 years with a proven history for performance*





**A brand you have trusted
for over 35 years.**



*Results indicative for human trials.

Please refer to the individual package IFUs for complete instructions, indications, contraindications, warnings and precautions.

Spatula

Needles specifications		Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
SS-24 1/4 Circle 8.6 mm Double-arm Curvature 90°		6-0	0.7	45	White	12	Surgalloy™	CD-3113K
		5-0	1	45	White	12	Surgalloy™	CD-3123K
		5-0	1	45	White	12	Surgalloy™	D-1764K
		4-0	1.5	45	White	12	Surgalloy™	D-1763K
SS-14 1/4 Circle 8.6 mm Double-arm Curvature 115°		6-0	0.7	45	White	12	Surgalloy™	D-1722K
		5-0	1	45	White	12	Surgalloy™	D-1760K
SS-28 3/8 Circle 7.4 mm Double-arm Curvature 145°		6-0	0.7	45	White	12	Surgalloy™	D-1683K
SS-2 1/2 Circle 8.6 mm Double-arm Curvature 200°		5-0	1	45	White	12	Surgalloy™	CD-0423K
		5-0	1	45	White	12	Surgalloy™	D-1780K
		4-0	1.5	45	White	12	Surgalloy™	D-1779K

Dermalon™ / Monosof™

Monofilament nylon

Dermalon™ / Monosof™ suture is a brand you have trusted for ~ 50 years with a proven history of performance.



**High tensile strength.
Low tissue reactivity.***











*Results indicative for human trials.

Please refer to the individual package IFUs for complete instructions, indications, contraindications, warnings and precautions.








Tapercutting

Needles specifications	Size	Length cm	Color	Units per box	Comments	Reorder code	
							USP
MVK-70-3 3/8 Circle 4 mm  MVK-70-3	10-0	0.2	13	Black	12	Surgalloy™	N-2513
	9-0	0.3	13	Black	12	Surgalloy™	N-2533
MVK-100-4 3/8 Circle 5 mm  MVK-100-4	10-0	0.2	13	Black	12	Surgalloy™	N-2531
	9-0	0.3	13	Black	12	Surgalloy™	N-2533

Premium spatula

Needles specifications	Image	Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
SE-100-8 1/4 Circle 6.3 mm Double-arm Curvature 100°		8-0	0.4	30	Black	12	Surgalloy™	N-2714K
SS-24 1/4 Circle 8.6 mm Double-arm Curvature 90°		5-0	1	45	Black	12	Surgalloy™	N-2781K
SE-140-6 3/8 Circle 6.3 mm Curvature 140°		10-0	0.2	20	Black	12	Surgalloy™	N-2730-K
SE-140-6 3/8 Circle 6.3 mm Double-arm Curvature 140°		10-0	0.2	30	Black	12	Surgalloy™	N-2770-K
		9-0	0.3	30	Black	12	Surgalloy™	N-2771-K
SE-140-8 3/8 Circle 6.3 mm Curvature 140°		10-0	0.2	20	Black	12	Surgalloy™	N-2736K
SE-140-8 3/8 Circle 6.3 mm Double-arm Curvature 140°		10-0	0.2	30	Black	12	Surgalloy™	N-2718K
		9-0	0.3	30	Black	12	Surgalloy™	N-2717-K
SE-CC-6 Compound curve 4.8 mm Double-arm Curvature 130°		10-0	0.2	30	Black	12	Surgalloy™	N-2799-K
SE-160-4 1/2 Circle 5.5 mm Double-arm Curvature 160°		11-0	0.1	30	Black	12	Surgalloy™	N-2705K
		10-0	0.2	30	Black	12	Surgalloy™	N-2704K

Spatula

Needles specifications	Image	Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
SE-160-6 1/2 Circle 5.5 mm Curvature 160°	 SE-160-6	10-0	0.2	20	Black	12	Surgalloy™	N-2734K
SE-160-6 1/2 Circle 5.5 mm Double-arm Curvature 160°	 SE-160-6	10-0	0.2	15	Black	12	Surgalloy™	N-2757K
		10-0	0.2	30	Black	12	Surgalloy™	N-2756K
		9-0	0.3	30	Black	12	Surgalloy™	N-2760K
SE-160-8 1/2 Circle 5.5 mm Curvature 160°	 SE-160-8	10-0	0.2	15	Black	12	Surgalloy™	N-2766K
SE-160-8 1/2 Circle 5.5 mm Double-arm Curvature 160°	 SE-160-8	10-0	0.2	30	Black	12	Surgalloy™	N-2786K
		9-0	0.3	30	Black	12	Surgalloy™	N-2787K
SE-175-8 1/2 Circle 6.9 mm Double-arm Curvature 175°	 SE-175-8	10-0	0.2	30	Black	12	Surgalloy™	N-2719K
		8-0	0.4	30	Black	12	Surgalloy™	N-2716K
SE-175-6 1/2 Circle 7 mm Double-arm Curvature 175°	 SE-175-6	10-0	0.2	30	Black	12		N-2966K
		9-0	0.3	30	Black	12		N-2967
SS-2 1/2 Circle 8.6 mm Double-arm Curvature 200°	 SS-2	5-0	1	45	Black	12	Surgalloy™	N-2761K

Novafil™

Monofilament polybutester

Novafil™ suture is a brand you have trusted for ~ 30 years with a proven history of performance.

**Supple and flexible.
Stretches with healing.***








*Results indicative for human trials.

Please refer to the individual package IFUs for complete instructions, indications, contraindications, warnings and precautions.


Spatula

Needles Specifications	Size		Length		Color	Units per box	Comments	Reorder code	
	USP	EP	cm	Inches					
SS-1 1/6 Circle 5.6 mm Double-arm Curvature 60°		5-0	1	30	12	Blue	12	Surgalloy™	PB-7740K
									
									
SS-24 1/4 Circle 8.6 mm Double-arm Curvature 90°		6-0	0.7	30	12	Blue	12	Surgalloy™	PB-6713K
		5-0	1	30	12	Blue	12	Surgalloy™	PB-6723K
									





Reverse cutting

Needles specifications		Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
HE-6 3/8 Circle 7.4 mm Double-arm Curvature 103°		7-0	0.5	45	Black	12	Surgalloy™	S-1768K
		6-0	0.7	45	Black	12	Surgalloy™	S-1769K
HE-5 3/8 Circle 8.6 mm Double-arm Curvature 135°		6-0	0.7	30	Black	12	Surgalloy™	S-1746K
HE-1 3/8 Circle 11 mm Double-arm Curvature 135°		6-0	0.7	45	Black	12	Surgalloy™	S-1750K
HE-3 3/8 Circle 13 mm Double-arm Curvature 155°		4-0	1.5	45	Black	12	Surgalloy™	S-1789K
		4-0	1.5	45	Black	12	Surgalloy™	S-1783K
HE-7 1/2 Circle 7.4 mm Double-arm Curvature 145°		7-0	0.5	45	Black	12	Surgalloy™	S-1766K
		6-0	0.7	45	Black	12	Surgalloy™	S-1765K

Taper point

Needles specifications	Image	Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
CV-1 3/8 Circle 9 mm		6-0	0.7	45	Black	36	Surgalloy™	VS-802
CV-1 3/8 Circle 9 mm Double-arm		7-0	0.5	45	Black	36	Surgalloy™	VS-809

Spatula

Needles specifications	Image	Size		Length cm	Color	Units per box	Comments	Reorder code
		USP	EP					
SE-100-8 1/4 Circle 6.3 mm Double-arm Curvature 100°		8-0	0.4	45	Black	12	Surgalloy™	S-2176K
SS-14 1/4 Circle 8.6 mm Double-arm Curvature 115°		6-0	0.7	45	Black	12	Surgalloy™	S-1780K
SS-24 1/4 Circle 8.6 mm Double-arm Curvature 90°		6-0	0.7	45	Black	12	Surgalloy™	S-2780K
		4-0	1.5	45	White	12	Surgalloy™	S-2782K
SE-140-8 3/8 Circle 6.3 mm Double-arm Curvature 140°		8-0	0.4	30	Black	12	Surgalloy™	S-2752K
		8-0	0.4	45	Black	12	Surgalloy™	S-1734K
		7-0	0.5	45	Black	12	Surgalloy™	S-1733K
		6-0	0.7	45	Black	12	Surgalloy™	S-1732K

Needles chart

Microsurgery

Code	Length
------	--------

⊙ taper point - 3/8 circle

MV-100-3	4 mm	⌋
MV-50-3	4 mm	⌋
MV-70-3	4 mm	⌋
MV-135-3	4 mm	⌋
MV-70-4	5 mm	⌋
MV-135-4	5 mm	⌋
MV-100-4	5 mm	⌋
CV-345	6 mm	⌋
MVF-135-5	6 mm	⌋
MV-135-5	6 mm	⌋
CV-351	8 mm	⌋
MV-175-8	8 mm	⌋
MVF-175-8	8 mm	⌋
MV-175-9	10 mm	⌋

⊙ tapercutting - 3/8 circle

MVK-70-3	4 mm	⌋
MVK-100-4	5 mm	⌋

Ophthalmic

Code	Length	Curvature	Comp
------	--------	-----------	------

▼ reverse cutting - 3/8 circle

HE-10	5.6 mm	130°	G-6* 110°	⌋
HE-6	7.4 mm	110°	G-6 110°	⌋
HE-7	7.4 mm	145°	G-7 165°	⌋
HE-5	8.6 mm	135°		⌋
HE-1	11.0 mm	135°	G-1 135°	⌋
HE-3	13.0 mm	155°	G-3 155°	⌋

▼ reverse cutting - 1/2 circle

HE-2	8.6 mm	205°	G-2 180°	⌋
------	--------	------	----------	---

⚡ inverted premium spatula - 1/4 circle

SS-29	7.4 mm	110°	S-29 110°	⌋
-------	--------	------	-----------	---

⚡ inverted premium spatula - 3/8 circle

SS-28	7.4 mm	145	S-28 163°	⌋
-------	--------	-----	-----------	---

⚡ premium spatula - 1/8 circle

SS-1	5.6 mm	60°	SM-1 60°	⌋
------	--------	-----	----------	---

⚡ premium spatula - 1/4 circle

SE-100-8	6.3 mm	100°	TG-100-8 97°	⌋
SD-1	8.1 mm	90°	RD-1 90°	⌋
SS-24	8.6 mm	95°	S-24 90°	⌋
SS-14	8.6 mm	115°	S-14 112°	⌋

⚡ premium spatula - 5/16 circle

SE-110-11	8.0 mm	110°	S-29 110° (previously SLO-110-11)	⌋
-----------	--------	------	--------------------------------------	---

⚡ premium spatula - 3/8 circle

SE-140-6	6.3 mm	140°	TG-140-6/CS-140-6 140°	⌋
SE-140-8	6.3 mm	140°	TG-140-8/CS-140-8 140°	⌋
SE-140-9	6.5 mm	140°	TG-140-8 137°	⌋

⚡ premium spatula - 1/2 circle

SE-160-4	5.5 mm	160°	TG-160-4/CS-160-4 160°	⌋
SE-160-6	5.5 mm	160°	TG-160-6/CS-160-6 160°	⌋
SE-160-8	5.5 mm	160°	TG-160-8/CS-160-8 160°	⌋
SE-160-9	5.5 mm	160°	TG-160-8/CS-160-8 160°	⌋
SE-175-6	7.0 mm	175°	TG-175-6/CS-175-6 175°	⌋
SE-175-8	6.9 mm	175°	TG-175-8/CS-175-8 175°	⌋
SS-22	7.4 mm	170°	S-22 170°	⌋
SS-2	8.6 mm	200°	S-2 180°	⌋

⚡ premium spatula - compound curve

SE-CC-6	4.8 mm	130°	TG6-S	⌋
---------	--------	------	-------	---

*Suggested

Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2024 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

^{TM*} Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

EMEA-WC-2300001-ophtalmic-needles-veterinary-catalogue-en-emea-10459613

[medtronic.com/animal-health/en-gb/](https://www.medtronic.com/animal-health/en-gb/)