

Medtronic

VERSAONE™ FAMILY TROCARS:

The access without the excess

VersaOne™ Bladed Trocar
VersaOne™ Access System



Operating room
essentials: the complete
solution you can trust.

Efficacy & simplicity.
Every time.

The VersaOne™ Bladed Trocar is available in 5 mm, 11 mm, and 12 mm sizes

Smoother insertion

Unique “dolphin nose” tip is designed to effectively dilate abdominal tissue layers with minimal insertion force

Port fixation and security

Advanced fixation ribs keep cannula from moving during surgical procedure

Protected insertion

Parabolic shield is designed to retract over the blade once in the peritoneal cavity




Safety Lock

indicates when the blade is shielded and exposed




Enhanced seal system

allows for easy instrument exchange²



Uniquely designed blade

divides tissue cleanly and precisely with control upon insertion through the abdominal wall



Proprietary ribbed cannula

for better port fixation and security compared to Applied Kii™ 12 mm trocar¹

VersaOne™ bladed trocar

Sold six units per box

	5 mm Smooth	5 mm Fixation	11 mm Smooth	11 mm Fixation	12 mm Smooth	12 mm Fixation
Short Length Cannula (70 mm)	B5SHS	B5SHF	-	-	-	-
Standard Length Cannula (100 mm)	B5STS	B5STF	B11STS	B11STF	B12STS	B12STF
Long Length Cannula (150 mm)	-	-	-	-	B12LGS	B12LGF

VersaOne™ universal cannula

Sold six units per box

	5 mm Fixation	11 mm Fixation	12 mm Smooth	12 mm Fixation
Short Length Cannula (70 mm)	UNVCA5SHF	-	-	-
Standard Length Cannula (100 mm)	UNVCA5STF	UNVCA11STF	UNVCA12STS	UNVCA12STF
Long Length Cannula (150 mm)	UNVCA5LGF	-	-	UNVCA12LGF

Contact your Medtronic representative today for more information about the VersaOne™ Access System or visit [medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

1. Based on internal test report #2143-114. Applies to 12 mm trocars, when compared to Kii™* (12 mm, Z-thread cannula). March 2013.

2. Based on internal test report #2143-123. Applies to Medtronic 12 mm trocars, when compared to the Versaport™ Bladeless 12 mm trocar. March 2013.

Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2022 Medtronic. Medtronic, Medtronic logo and Engineering the extraordinary are trademarks of Medtronic.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. EMEA-AI-2200002 - 21-weu-versaone-bladed-info-sheet-5501592

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)