

Medtronic

Endo Catch™ Gold specimen retrieval pouch



Product information

Code	Description	Qty per case
173050G	Endo Catch™ Gold 10 mm 2.5" × 6" specimen pouch, volume capacity 220 ml, 29.5 cm shaft length	6

Materials of manufacture

303, 304 stainless steel, polycarbonate 20% glass, PTFE, polycarbonate GE LEXAN, rubber/silicone durometer, polypropylene red, tincl bb, RNF 150, polyurethane

Shipping information

	Shipped case dimensions (mm)
Inner box	L 152 × W 111 × H 562

Additional information

Latex free	✓
DEHP free	✓
Sterilization method	Ethylene oxide
Shelf life	3 years
CE marked	Class IIa according to directive 93/42/EEC on Medical Devices (MDD), annex II excluding (4)

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2024 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

All other brands are trademarks of a Medtronic company. emea-ai-2000020-endocatch-gold-tech-sheet-en-gb-emea-13102949

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)