

**Medtronic**

FlexCath Cross™ Transseptal Solution

**Zero**  
exchange  
workflow





## Efficiency

# Zero exchange procedure

### Zero wire and needle exchanges


Integrated needle/dilator design removes the need for needle and guidewire exchanges.

### Zero sheath exchanges


Dilators compatible with leading industry sheaths eliminate the need for separate transseptal and therapeutic sheaths.

# Efficiency Workflow

## Traditional transseptal access workflow



## FlexCath Cross workflow with FlexCath Advance



## Versatile

# Versatility for your procedures

Choose the FlexCath Cross model which fits with your existing sheath preferences for your left-heart procedures.

FlexCath Cross is the only transseptal access system on the market **cleared for both mechanical and RF crossing.**



## Compatible with industry-leading introducers



Abbott Agilis™ NxT  
Steerable Introducer



Medtronic FlexCath Cross™  
with Mini Sheath



Boston Scientific Watchman™  
TruSeal™ and FXD™  
Access Systems++



Biosense Webster CARTO  
VIZIGO® Bi-directional  
Guiding Sheath



Merit Medical HeartSpan®  
Sheath Introducer



Abbott Swartz™ Braided  
Transseptal Guiding  
Introducer, SL1

## Controlled

# Cross with control

The FlexCath Cross needle only deploys when activated, designed for ease of control during transseptal puncture.



Ordering info

Product name	Compatible sheath	Order number
FlexCath Cross – AG 61 cm	Agilis NxT – 61 cm	900300
FlexCath Cross – AG 71 cm	Agilis NxT – 71 cm	900301
FlexCath Cross – SL 63 cm	Swartz SL1 – 63 cm	900302
FlexCath Cross – SL 81 cm	Swartz SL1 – 81 cm	900309
FlexCath Cross – VZ 71 cm	CARTO VIZIGO – 71 cm	900303
FlexCath Cross – FC 65 cm	FlexCath Advance – 65 cm	900304
FlexCath Cross – MH 63 cm <sup>++</sup>	Merit Medical HeartSpan Sheath – 63 cm	900306
FlexCath Cross – MH 81 cm <sup>++</sup>	Merit Medical HeartSpan Sheath – 81 cm	900307
FlexCath Cross – WT 75 cm <sup>++</sup>	Watchman TruSeal and FXD Access Systems – 75 cm	900308
FlexCath Cross with Mini Sheath– 65 cm	65 cm, 90 degrees, 8.5 Fr ID, 11 Fr OD	012852

Competitive matrix – transseptal needles

Feature	FlexCath Cross	BRK <sup>TM*</sup>	SafeSept <sup>**</sup>	NRG <sup>**</sup>	VersaCross <sup>**</sup>
Integrated needle/dilator design reduces exchanges	✓	✗	✗	✗	✓
Spring-tensioned safety needle	✓	✗	✗	✗	✗
Indicated for use mechanically or with RF	✓	✗	✗	✗	✗
Reposition without requiring wire/needle exchange	✓	✗	✗	✗	✓
Engineered to pair with leading sheaths	✓	✗	✗	✗	✗

Product CE mark according to law.

++ FlexCath Cross configurations are not CE marked nor commercially available in Europe.

See the device manual for information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [www.medtronic.eu](http://www.medtronic.eu).

For applicable products, consult instructions for use on [manuals.medtronic.com](http://manuals.medtronic.com). Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

CE 0297

## Medtronic

### Europe

Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

### United Kingdom/Ireland

Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 (0)1923 212213  
Fax: +44 (0)1923 241004

©2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. FlexCath Cross™ is a trademark of Medtronic, manufactured by Acutus Medical, and is exclusively distributed by Medtronic. ™\*Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. UC202313327 EE-flexcath-cross-brochure-en-emea-10711888