



# Medtronic

Evolut™ FX+ TAVI System

# Product specifications



## Bioprosthesis




### Evolut™ FX+ 23 mm Valve

Model number:  
EVFXPLUS-23

Size: 23 mm

Aortic annulus  
diameter:  
17<sup>†</sup>/18–20 mm




### Evolut™ FX+ 26 mm Valve

Model number:  
EVFXPLUS-26

Size: 26 mm

Aortic annulus  
diameter:  
20–23 mm




### Evolut™ FX+ 29 mm Valve

Model number:  
EVFXPLUS-29

Size: 29 mm

Aortic annulus  
diameter:  
23–26 mm




### Evolut™ FX+ 34 mm Valve

Model number:  
EVFXPLUS-34

Size: 34 mm

Aortic annulus  
diameter:  
26–30 mm


## Loading system



### Evolut™ FX+ Loading System

Model number:  
L-EVOLUTFX-2329

23, 26, and 29 mm TAV



### Evolut™ FX+ Loading System

Model number:  
L-EVOLUTFX-34

34 mm TAV

<sup>†</sup> Measurement for TAV-in-SAV only.

All items above are single-use. Refer to Instructions for Use for a complete list of warnings, precautions, and contraindications.

# Delivery catheter system

## Evolut™ FX+


### Delivery Catheter System with InLine™ Sheath

14 Fr equivalent (6.0 mm capsule outer diameter)

Model number:

D-EVOLUTFX-2329

23, 26, and 29 mm TAV



## Evolut™ FX+


### Delivery Catheter System with InLine™ Sheath

18 Fr equivalent (7.33 mm capsule outer diameter)

Model number:

D-EVOLUTFX-34

34 mm TAV



## Medtronic

Medtronic International Trading Sarl  
Route du Molliat 31  
Case postale  
1131 Tolochenaz  
Switzerland  
Tel: +41 (0) 21 802 70 00  
Fax: +41 (0) 21 802 79 00

[medtronic.eu](http://medtronic.eu)

© 2024 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.  
13791940-en-gb-emea