


# Product specifications

## Bioprosthesis




### Evolut FX 23 mm Valve

Model number:  
EVOLUTFX-23

Size: 23 mm

Aortic annulus  
diameter:  
18-20 mm




### Evolut FX 26 mm Valve

Model number:  
EVOLUTFX-26

Size: 26 mm

Aortic annulus  
diameter:  
20-23 mm




### Evolut FX 29 mm Valve

Model number:  
EVOLUTFX-29

Size: 29 mm

Aortic annulus  
diameter:  
23-26 mm




### Evolut FX 34 mm Valve

Model number:  
EVOLUTFX-34

Size: 34 mm

Aortic annulus  
diameter:  
26-30 mm


## Loading system



### Evolut FX Loading System

Model number:  
L-EVOLUTFX-2329

23, 26, and 29 mm TAV



### Evolut FX Loading System


Model number:  
L-EVOLUTFX-34

34 mm TAV

# Delivery catheter system


## Evolut FX delivery catheter system with InLine sheath 14 Fr equivalent (6.0 mm capsule outer diameter)

Model number:  
D-EVOLUTFX-2329  
23, 26, and 29 mm TAV



## Evolut FX delivery catheter system with InLine sheath 18 Fr equivalent (7.33 mm capsule outer diameter)

Model number:  
D-EVOLUTFX-34  
34 mm TAV



See the CoreValve™ Evolut™ R, the CoreValve™ Evolut™ PRO, the Evolut™ PRO+ device and the Evolut™ FX device manuals for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [medtronic.eu](http://medtronic.eu).

For applicable products, consult instructions for use on [manuals.medtronic.com](http://manuals.medtronic.com). Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

The commercial name of the Evolut™ R device is Medtronic CoreValve™ Evolut™ R System, the commercial name of the Evolut™ PRO device is Medtronic CoreValve™ Evolut™ PRO System, the commercial name of the Evolut™ PRO+ device is Medtronic Evolut™ PRO+ System, and the commercial name of the Evolut™ FX device is Medtronic Evolut™ FX System.

CE  
0344

**Medtronic**

Europe  
Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)

UC202203007-evolut-fx-product-specs-en-we-8648857

©2023 Medtronic. All rights reserved.

Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

™\*Third-party brands are trademarks of their respective owners.

All other brands are trademarks of a Medtronic company.

[medtronic.eu](http://medtronic.eu)