

Medtronic

Engineering the extraordinary

Designed with purpose

With a long history of innovation and a reputation for quality, we design products with your critical needs in mind.


Concerto™ Detachable Coil Systems are designed to help you easily and precisely get to the embolization area and instantly detach the coil on target.¹




Proven
reliability

99.7%


Detachment
success¹



Concerto™
Detachable Coil System





Attached



Detached

The unique Concerto™ Coil design

- Ball-and-socket design allows the coil implant to move independently of the pusher in a 360° fashion, resulting in **less energy buildup and microcatheter kickout.**¹
- Bevel retainer design allows for **articulation zone flexibility.**¹
- Soft distal pusher **minimizes catheter kickback.**¹



360°

1. Test data on file at Medtronic internal reports TR13-067 rev B / TR13-081 Rev A / TR12-054 REV A / TR14-061 Rev A / TR14-062 Rev A / TR14-065 Rev A.

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.eu.

Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

UC202001818a-concerto-detachment-infosheet-en-we-8480725
© Medtronic 2023.

medtronic.eu/embolization