

McGRATH™ video laryngoscopes

Comparison guide

McGRATH™ MAC video laryngoscope



McGRATH™ MAC+ video laryngoscope

4.5x
higher screen resolution



45%
larger screen size

	McGRATH™ MAC video laryngoscope	McGRATH™ MAC+ video laryngoscope
Specification		
Weight	0.172 kg	0.180 kg
Useful life	4 years or 2000 decontamination cycles	4 years or 2000 decontamination cycles
Visualization		
Screen resolution	320x240	720x480 (4.5x increase)
Screen Size	2.4 inch	3 inch (45% increase)
Durability		
Handle	Reinforced screen hinge	Reinforced screen hinge
IP rating	IPX8	IPX8
Battery management		
Battery type	Disposable and rechargeable batteries with auto-off feature	Rechargeable battery with auto-off feature
Disposable battery life	250 minutes	N/A
Rechargeable battery life	10,000 minutes	10,000 minutes
Run time on full charge	120 min (with rechargeable battery)	60 min (with rechargeable battery)
Ordering information and SKU numbers		
Handle	301-000-000	302-000-000
3.6V disposable battery	340-000-000	N/A
Rechargeable battery	341-100-000	341-100-000
Relay	N/A	302-100-000
Connected screen	N/A	302-200-000
MAC + Features		
Streaming and file management	NO	Yes
Video recording	NO	Yes
Touch screen	NO	Yes
Still image capture	NO	Yes
File storage (video and still image)	NO	Yes
Wireless connectivity	NO	Yes

This material should not be considered the exclusive source of information, it does not replace or supersede information contained in the device manual(s).

Please note that the intended use of a product may vary depending on geographical approvals.

See the device manual(s) for detailed information regarding the intended use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events.

For a MRI compatible device(s), consult the MRI information in the device manual(s) before performing a MRI.

If a device is eligible for eIFU usage, instructions for use can be found at Medtronic's website manuals.medtronic.com.

Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Medtronic products placed on European markets bear the CE mark and the UKCA mark (if applicable).

For any further information, contact your local Medtronic representative and/or consult Medtronic's websites.

For trained personnel only. For specific indications and instructions for use, please refer to the product manual.

©2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.
emea-re-2500043-mcgrath-video-laryngoscopes-comparison-guide-en-gb-emea-20431023

Medtronic