

Medtronic

Shiley™ Adult Flexible Evac Tracheostomy Tube with TaperGuard™ Cuff Disposable Inner Cannula



Features

- Available in three configurations, Shiley™ flexible evac
- Tracheostomy tubes offer design features that help both clinicians and patients
- Transparent flange
- Fully integrated cuff inflation line and evac portal
- Fits with Shiley™ Adult Flexible disposable inner cannula

CFN	I.D, mm	O.D, mm	Cuff Ø, mm
4CN65ED	6.5	9.4	20.6
5CN70ED	7.0	10.1	23.0
6CN75ED	7.5	10.8	25.4
7CN80ED	8.0	11.4	25.4
8CN85ED	8.5	12.2	26.6
9CN90ED	9.0	12.7	27.6
10CN10ED	10.0	13.8	28.6

Additional information

Latex free	Yes
PVC	Yes
Sterilization method	Ethylene oxide (ETO)
Shelf life	60 months
Storage requirements	Store at room temperature. Do not use if the container is opened or damaged.
Regulatory	MDD class IIb
Use	Single use
Elimination	Dispose in accordance with local disposal policy for medical devices.
Package	10/10
Country of manufacture	Ireland

Medtronic

Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2025 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

All other brands are trademarks of a Medtronic company.

emea-re-2400081-shiley-adult-flexible-evac-tracheostomy-w-disposable-inner-cannula--data-sheet-en-gb-emea-13943653

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)

Medtronic Belgium NV / SA

Burgemeester Etienne
Demunterlaan 5, Brussels 1090
Belgique / Belgien / België
Tel: +32 (0)2 456 09 00
Fax: +32 (0)2 460 26 67

Medtronic AB

Hans Werthéns Gata 19
16440 Kista
Sverige
Tel: +46 (0)8 56 85 85 00
Fax: +46 (0)8 56 85 85 01

Medtronic Danmark A/S

Arne Jacobsens Alle 17,
Copenhagen S 2300
Danmark
Tel: +45 32 48 18 00
Fax: +45 32 48 18 01

Medtronic Norge AS

Technopolis Fornebu
Martin Linges vei 25
1364 Fornebu
Norge
Tel: + 47 67 10 32 00
Fax: + 47 67 10 32 10

Medtronic Finland Oy

Lentäjantie 3,
Vantaa 01530
Suomi
Tel: +358 9725 19 88
Fax: +358 9725 19 89

Medtronic Trading NL B.V.

Larixplein 4
5616 VB Eindhoven
Nederland
Tel: +31 (0)40 711 75 55
Fax: +31 (0)40 711 75 22

[medtronic.com/covidien/se](https://www.medtronic.com/covidien/se)

[medtronic.com/dk-da/index.html](https://www.medtronic.com/dk-da/index.html)

[medtronic.com/covidien/no](https://www.medtronic.com/covidien/no)

[medtronic.com/covidien/fi](https://www.medtronic.com/covidien/fi)

[medtronic.com/covidien/nl](https://www.medtronic.com/covidien/nl)