

Medtronic

Engineering the extraordinary

Midas Rex™ MR8™

High-Speed Drill System

.....
Advancing surgical
possibilities



ClearView™
Enhanced Surgical
Visibility



The **Midas Rex™ MR8™** high-speed drill system has a lower operating temperature, less chatter, improved visibility of the surgical site and even better cutting performance.¹

Navigable with StealthStation™ S8, it seamlessly integrates with other technologies in our portfolio. With our comprehensive support, education, training and service, you have a partner that's there for you from start to finish.

Midas Rex™ MR8™.

It's not just a drill.
It's advancing surgical possibilities.

ClearView™

Variety of shapes lengths, and cutting geometries for various applications

Extended, self-irrigating sleeve protects surrounding tissue from minimally exposed rotating shaft

Low-profile integrated irrigation that minimizes obstructions while delivering irrigation directly near the cutting site

Direct connect base for one-step assembly

Visibility and protected shaft




Spine solutions

- Better surgical visibility and access for minimally invasive surgeries¹
- Navigate select ClearView™ tools with the StealthStation™ S8




Spine



Cranial solutions

- Lateral Skull Base tools now available~ 5 cm and 7 cm
- Better surgical visibility and access for skull base and endoscopic surgery¹
- Navigate select transnasal and lateral skull base tools in cranial applications with the StealthStation™ S8

Cranial



Ordering information

Spine

12 cm length, distal end	12 cm length, proximal bend	14 cm length, proximal bend
MR8-SD12BA30	MR8-SP12BA30*	MR8-SP14BA30*
MR8-SD12BA40	MR8-SP12BA40*	MR8-SP14BA40*
MR8-SD12MH25	MR8-SP12MH25	MR8-SP14MH25
MR8-SD12MH25D	MR8-SP12MH25D	MR8-SP14MH25D
MR8-SD12MH25T	MR8-SP12MH25T	MR8-SP14MH25T
MR8-SD12MH30	MR8-SP12MH30	MR8-SP14MH30
MR8-SD12MH30D	MR8-SP12MH30D	MR8-SP14MH30D
MR8-SD12MH30T	MR8-SP12MH30T*	MR8-SP14MH30T*

Transnasal

13 cm length
MR8-ND13BA20D*
MR8-ND13MH30D
MR8-ND13MH30
MR8-ND13BA30DC*
MR8-ND13BA30
MR8-ND13BA40DC*
MR8-ND13BA40
MR8-ND13BA45DC

Lateral skull base

5 cm length, proximal bend	7 cm length, proximal bend
MR8-LP05BA05D	MR8-LP07BA10
MR8-LP05BA08D	MR8-LP07BA10D*
MR8-LP05BA10	MR8-LP07BA15
MR8-LP05BA10D*	MR8-LP07BA20*
MR8-LP05BA12D	MR8-LP07BA20D*
MR8-LP05BA15	MR8-LP07BA20DC
MR8-LP05BA15D	MR8-LP07BA30*
MR8-LP05BA20*	MR8-LP07BA30D*
MR8-LP05BA20D*	MR8-LP07BA30DC
MR8-LP05BA20DC	MR8-LP07BA40
MR8-LP05BA30*	MR8-LP07BA40D
MR8-LP05BA30D*	MR8-LP07BA40DC
MR8-LP05BA30DC	
MR8-LP05BA40	
MR8-LP05BA40D	
MR8-LP05BA40DC	

Irrigation

IRD870	TUBULURE IRD870 MR8™ IRIGTN CV LOW-PRFL EA
IRD875	TUBULURE IRD875 MR8™ IRIGTN CV LOW-PRFL 5PK

*Indicates Stealth-Midas MR8 compatibility

1. Medtronic data on file. Compared to previous models.

This material should not be considered the exclusive source of information, it does not replace or supersede information contained in the device manual(s). Please note that the intended use of a product may vary depending on geographical approvals. See the device manual(s) for detailed information regarding the intended use, the (implant) procedure, indications, contraindications, warnings, precautions, and potential adverse events. For a MRI compatible device(s), consult the MRI information in the device manual(s) before performing a MRI. If a device is eligible for eFU usage, instructions for use can be found at Medtronic's website manuals.medtronic.com. Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser. Medtronic products placed on European markets comply with EU and UK legislation (if applicable) on medical devices. For any further information, contact your local Medtronic representative and/or consult Medtronic's websites.

Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliou 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

United Kingdom/Ireland
Medtronic Limited
Building 9, Croxley Park
Hatters Lane, Watford
Herts WD18 8WW
www.medtronic.co.uk
Tel: +44 (0)1923 212213
Fax: +44 (0)1923 241004

15533022-en-gb-emea © Medtronic 2026.
All rights reserved.